# IKAMAG® RCT basic

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*OPERATING INSTRUCTIONS* GB 8

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D

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## CE-DECLARATION OF CONFORMITY

GB

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 60 010-1; EN 50 082; EN 55 014; EN 55 114 and EN 60 555.

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## Garantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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## Contens

|  | Page |
|--|------|
| Guarantee  | 2    |
| Safety instructions  | 8    |
| Proper use   | 8    |
| Unpacking  | 9    |
| Commissioning  | 9    |
| Starting   | 9    |
| Function Heating   | 9    |
| Controlling the Medium Temperature via ETS-D<br>or Contact Thermometer | 10   |
| Function Stirring  | 11   |
| Maintenance and Cleaning   | 11   |
| Associated standards and regulations                                   | 12   |
| Accessories  | 12   |
| Technical data   | 12   |
| List of spare parts RCT basic  | 18   |
| Spare parts diagram RCT basic  | 19   |

## Safety instructions

Caution: Never heat up any medium with ignition point higher than the fix adjusted safety temperature limit. Value for safety temperature limit see chapter „Technical data“ .

The instrument may not be operated in rooms with explosion hazard.

Furthermore, it has to be placed upon a fire-proof and/or not burnable surface.

Please make sure that the mains cable does not contact the heating plate!

**Attention - Magnetism!** Effects of the magnetic field have to be taken into account (e.g. data carriers, cardiac pacemakers ...).

When using PTFE-coated magnetic bars, the following have to be noted:

*Chemical reactions of PTFE occur in contact with molten or dissolved alkaline and alkaline - earth metals, as well as with fine-particled powders of metals of the 2. and 3. group of the periodical system at temperatures above 300-400°C.*

*Only elementary fluorine, chlorine trifluoride und alkaline metals do attack PTFE, halogen hydrocarbons have a reversibly swelling effect.*

Source: Römpps Chemie-Lexikon and „Ullmann“ Bd.19

Exercise caution when touching the housing parts and the heating plate!

### Risk of burns!

The heating plate can reach temperatures in excess of 300°C.

Please consider possible endangering due to release of toxic or flammable gas caused by heated medium.

## Proper use

The IKAMAG RCT basic is a stirring instrument with heating function. It is used in e.g. laboratories of the chemical industry, in schools as well as in pharmacies. The instrument is suitable for tempering substances which are placed on the heating plate in vessels. By means of a built-in stirring drive and with the help of a magnetic bar within the vessel, the substances can be stirred at the same time. The mixing intensity is dependent on the motor speed and on the size of the magnetic bar.


## Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers: an IKAMAG RCT basic and operating instructions.

## Starting

After the instrument has been switched on, it is in „standby“-operation. The functions „heating“ and „stirring“ are started with the help of switches on the left and right bottom of the housing (see illustration 1)

## Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged. Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

## Function Heating

The temperature of the heating plate of the instrument is kept constant by a control circuit. In addition, it is monitored by a safety circuit. The required temperature sensor - a thermo-couple - is installed in the heating plate.

The instrument disposes of a silumin heating plate with a heating power of 600 watts. The surface of the heating plate is polished in order to facilitate cleaning. This surface gets coated rather quickly with an aluminium oxide protective layer which is chemically very resistend.

The function „heating“ is switched on and/or off by its corresponding switch. When the function „heating“ is switched on, the LED „temp“ shines green. The temperature of the heating plate is set on the rotary knob „temp“ with its corresponding dial (0-300°C). The heating plate is heated up to the temperature set on the rotary knob. During the heating phase, the color of the LED „temp“ changes from green to orange (orange = supply of energy to heating plate; green = set temperature is reached).


In the event of a disturbance in the control circuit, the heating plate is **permanently cut off** by the safety circuit.

#### PLEASE NOTE!

The heating plate only heats up, if the contact plug has been plugged (the contact plug is undetachably fitted to the instrument), or if an ETS-D and/or a contact thermometer has been connected.

(see illustration 2)

illustration 2

## Controlling the Medium Temperature via ETS-D or Contact Thermometer

The control of the medium temperature via ETS-D or contact-thermometer should be preferred. Its advantages are: just a short heating-up phase after setting the rated temperature, practically no temperature drift and only a minor temperature waviness. The medium temperature is optimized with the aid of a **FUZZY Logic Controller**, that means, that in steady state any overshooting of the temperature is reduced to a minimum.

The rear of the instrument houses the diode bushing for connection the ETS-D or contact thermometer, or for plugging the contact plug (see illustration 2).

**The electronics of the instrument supply a test current that has to flow over plug pins 3 and 5 of the diode bushing in order for the heating plate to heat.**

#### Safety contact thermometers

acc. to DIN 12 878 class 2 or acc. Gerstel are connected with a three-wire cable, the test current flows through the contact thermometer.

#### Safety function:

If the test current is interrupted because of e.g. breakage of contact thermometer or falling out of the cable plug, the heating cuts off. The required plug allocation is indicated on a reference plate above the diode bushing.

#### Contact thermometer without safety circuit

acc. to DIN 12 878 class 0. The instrument only heats if the test current circuit is closed by an electrical connection of the plug pins 3 and 5.

#### 2-wire connecting cables:

Connect plug pins 3 and 5 of the instrument plug.

### 3-wire connecting cables:

Here the test current circuit can also be produced in the terminal head of the contact thermometer (connect plug pins 2 and 3).

A 3-wire cable with the required bridge is available (accessories).

The **rated temperature** is set on the **contact thermometer** by means of the rotary magnetic head, whereby the upper edge of the traveling nut is the reference point on the rated value scale.

The required **surface temperature of the heating plate** is preselected on the **temp-rotary knob** of the instrument.

If the temp-rotary knob is set to the maximally adjustable temperature, the heating-up takes place as rapidly as possible, yet the medium temperature overshoots the rated value set on the contact thermometer.

If the temp-rotary knob is set to twice the rated value (e.g. with a rated value of 60 °C the temp-rotary knob is set to 120 °C), the result is good compromise between rapid heating-up and overshooting.

If the temp-rotary knob is set precisely to the rated temperature, the medium does not reach the rated temperature owing to a drop in temperature between heating plate and medium.

The **rated temperature** is setting on the **ETS-D**.

When using the ETS-D, there is no overshooting of the medium temperature owing to an „aiming-off allowance“-technique.

## Function Stirring

The funktion stirring is started and stopped by its corresponding switch.

The motor speed is set on the rotary knob „mot“ with its dial from 1 to 10.

If the motor is switched on, it runs-up to the speed set on the rotary knob, and the LED „mot“ beside the imprint „Mot“ shines green.

The dial from 1 - 10 conforms with a speed range from 50 . 1100 1/min.

When idling, the speed is 1100 1/min at maximum setting of the rotary knob (right-hand stop).

If the viscosity increases during the process, the speed drops a little. If the viscosity decreases, e.g. due to adding thin fluid media, the speed increases a little. By a corresponding resetting of the rotary knob, both of these cases can be balanced, as long as the torque does not exceed its maximum.

Please note that mains fluctuations within the permissible tolerance also effect minor fluctuations of the speed.

## Maintenance and Cleaning

The IKAMAG RCT basic is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

**When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.**

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

|                        |  |
|------------------------|--|
| Dyes                   | isopropyl alcohol                            |
| Construction materials | water containing tenside / isopropyl alcohol |
| Cosmetics              | water containing tenside / isopropyl alcohol |
| Foodstuffs             | water containing tenside                     |
| Fuels                  | water containing tenside                     |

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

## Associated standards and regulations

### Construction in accordance with the foolwing safety standards

EN 61 010-1 UL 3101-1 CAN/CSA C22.2 (1010-1)

### Construction in accordance with the foolwing EMC standards

EN 50 082-1 EN 55 014 EN EN60 555-2: -3 EN 50 081

### Associated EU guidlines

EMC-guidlines: 89/336/EWG

Machine guidlines: 73/023/EWG

## Accessories

|              |                           |              |                       |
|--------------|---------------------------|--------------|-----------------------|
| <b>RS 1</b>  | set of stirring bars      | <b>H16V</b>  | support rod           |
| <b>RSE</b>   | PTFE stirring bar remover | <b>H 43</b>  | contact safety devive |
| <b>R 380</b> | stand support             | <b>H 15</b>  | bath top              |
| <b>R 360</b> | holding fork              | <b>ETS-D</b> | IKATRON               |
| <b>H 36</b>  | holding rod               | <b>H 44</b>  | cross sleeve          |

### Recommended IKA-stirres

**TRIKA**-stirrers PTFE, triangular  
lengths 25, 42, 55, und 80 mm

**IKAFLO**N-stirrers PTFE  
dia 6mm lengths 10 and 15 mm  
dia 7mm lengths 20, 25, 30, 40, 50 and 60 mm  
dia 10mm lengths 70 and 80 mm

## Technical data

|   |              |   |
|---|--------------|---|
| design voltage:   | <b>VAC</b>   | 230±10%   |
| design frequency:   | <b>VAC</b>   | 115±10%   |
| <b>motor:</b>   | <b>Hz</b>    | 50/60   |
|   |              | EC-motor with ball bea<br>rings, without brushes            |
| max. power input of motor   | <b>W</b>     | 17  |
| speed range:  | <b>1/min</b> | 50 bis 1100   |
|   |              | infinitely variable setting<br>silumin, GK AlSi 12 polished |
| <b>heating plate:</b>   |              |   |
| diameter of heating plate:  | <b>mm</b>    | 135   |
| heating power:  | <b>W</b>     | 600   |
| max. surface temperature:   | <b>°C</b>    | 300   |
| max. surface temperature<br>under error conditions                          | <b>°C</b>    | 370   |
| temperatur waviness<br>of heating plate                                     | <b>K</b>     | ±2,5  |
| temperature setting:  | <b>°C</b>    | infinitely variable from<br>ambient temperatur              |
| control accuracy of medium<br>in ETS-D or kontakt-<br>thermometer operation | <b>K</b>     | ±1 FUZZY-optimized  |
| contact thermometer-<br>connection  |              | DIN 12 878 Klasse 2   |
| permitted ambient<br>temperature:   | <b>°C</b>    | +5 bis +40  |
| permitted relative humidity:  | <b>%</b>     | 80  |
| permitted duration of operation:  | <b>%</b>     | 100   |
| prozection class acc. DIN 40 050:   |              | IP 42   |
| fuses:  | F1/F2        | 10AT/250V Id-Nr 2755400                                     |
|   | F3           | 1,025AT/250 Id-Nr 2756300                                   |
| contamination level:  |              | 2   |
| overvoltage category:   |              | II  |
| protection class:   |              | 1(protective earth)   |
| dimensions:   | <b>mm</b>    | 160 x 280 x 95  |
| weight:   | <b>kg</b>    | 2,4   |

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## List of spare parts RCT b


| Item                         | Designation                  | Item                         | Designation                  |
|------------------------------|------------------------------|------------------------------|------------------------------|
| 2 bushing                    | 46 plastic cover             | 2 bushing                    | 46 plastic cover             |
| 3 radiating sheet            | 47 fixing plate              | 3 radiating sheet            | 47 fixing plate              |
| 4 sealing washer             | 48 sealing washer            | 4 sealing washer             | 48 sealing washer            |
| 5 washer                     | 102 ceramic stone            | 5 washer                     | 102 ceramic stone            |
| 6 taptite cheese-head screw  | 109 receptacle (gold-plated) | 6 taptite cheese-head screw  | 109 receptacle (gold-plated) |
| 7 ceramic bushing            | 110 insulating hose          | 7 ceramic bushing            | 110 insulating hose          |
| 8 sleeve                     | 2001 heating plate           | 8 sleeve                     | 2001 heating plate           |
| 9 foil keyboard              | 2002 EC-motor                | 9 foil keyboard              | 2002 EC-motor                |
| 10 housing-upper part        | 3001 poti                    | 10 housing-upper part        | 3001 poti                    |
| 11 potentiometer shaft       | 3002 main plug               | 11 potentiometer shaft       | 3002 main plug               |
| 12 clamping bushing          | 3003 rocker switch           | 12 clamping bushing          | 3003 rocker switch           |
| 13 hood                      | 3004 bushing plug            | 13 hood                      | 3004 bushing plug            |
| 14 hexagonal nut             | 3005 rocker switch           | 14 hexagonal nut             | 3005 rocker switch           |
| 15 cover for Pt 100          | 3006 relay                   | 15 cover for Pt 100          | 3006 relay                   |
| 17 fixing clip               | 3007 fuse F1/F2              | 17 fixing clip               | 3007 fuse F1/F2              |
| 18 taptite cheese-head screw | 3008 fuse F3                 | 18 taptite cheese-head screw | 3008 fuse F3                 |
| 19 toothed washer            | 5001 ball bearing            | 19 toothed washer            | 5001 ball bearing            |
| 20 washer                    | 5002 grooved ball bearing    | 20 washer                    | 5002 grooved ball bearing    |
| 21 O-ring (Perbunan)         |                              | 21 O-ring (Perbunan)         |                              |
| 31 bottom part, coated       |                              | 31 bottom part, coated       |                              |
| 32 plastite-screw            |                              | 32 plastite-screw            |                              |
| 33 BLP analog                |                              | 33 BLP analog                |                              |
| 34 bushing                   |                              | 34 bushing                   |                              |
| 35 magnet                    |                              | 35 magnet                    |                              |
| 36 countersunk screw         |                              | 36 countersunk screw         |                              |
| 37 plastic cover             |                              | 37 plastic cover             |                              |
| 38 taptite countersunk screw |                              | 38 taptite countersunk screw |                              |
| 39 cap                       |                              | 39 cap                       |                              |
| 40 srewed plug               |                              | 40 srewed plug               |                              |
| 41 operating knob poti       |                              | 41 operating knob poti       |                              |
| 43 contact thermometer label |                              | 43 contact thermometer label |                              |
| 44 contact plug              |                              | 44 contact plug              |                              |
| 45 taptite cheese-head screw |                              | 45 taptite cheese-head screw |                              |

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**IKA® WORKS, INC.**

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

**IKA® Works, (Asia) Sdn Bhd**

LABORATORY TECHNOLOGY  
ANALYZING TECHNOLOGY  
PROCESSING EQUIPMENT

---

The Global Market of Laboratory Equipment  
<http://www.labworld-online.com>

---

Europa - Afrika

**IKA® -WERKE GMBH & CO.KG**

Janke & Kunkel-Str. 10  
D 79 219 STAUFEN  
GERMANY  
TEL: 07633/831-0  
FAX 07633/831-98  
E-mail: sales@ika.de  
Internet: <http://www.ika.net>

America

**IKA® WORKS, INC.**

2635 NORTH CHASE PKWY. SE  
WILMINGTON, NC 28405-7419  
TEL. 800/733-3037  
TEL. 910/452-7059  
FAX 910/452-7693  
E-mail: usa@ika.net

Asien - Australien

**IKA® Works (Asia) Sdn Bhd**  
(Company No. 340448-K)

Lot 2, Jalan Indah 1/2  
Taman Industri Rawang Indah  
48000 Rawang  
Selangor, Malaysia  
TEL : (603) 6093 3322  
FAX : (603) 6093 3940  
E-mail: [ika@tm.net.my](mailto:ika@tm.net.my)