

3D OCT-2000 Series*

Optical Coherence Tomography



PERFORMANCE
YOU CAN COUNT ON



The Ultimate All-in-one Product
Created to be Fast, Easy and Precise

3D OCT-2000 Series

Optical Coherence Tomography




The TOPCON 3D OCT-2000 Series is an optimal choice for all eye care professionals

The 3D OCT-2000 Series of Spectral Domain OCTs with High Resolution Fundus Cameras was designed to meet the needs of a comprehensive fundus imaging device for all eye care professionals from the single doctor practice to a large university hospital.


3D OCT-2000 *FA plus*

» OCT, Color, Red-free, FA, FAF images acquirable



3D OCT-2000 *FA*

» OCT, Color, Red-free, FA images acquirable




3D OCT-2000

» OCT, Color, Red-free images acquirable



* Display digital Red-free



Sophisticated OCT Technology with the Integrated Fundus Camera

12 mm Wide Scan & 5 Line Cross Scan!
Rich Analysis Functions Plus High Resolution Images

50,000 A-scans/sec – Greater Details in Shorter Time









The enhanced 50,000 A-scans/sec allows for faster tomography acquisition and minimizes artifacts generated by eye movement, producing clear cross-sectional retinal images. Now there are even more imaging variations by the TOPCON 3D OCT-2000 with the new 12×9mm wide scan enabling the user to capture a wider range of the image from optic disc to macula with a single shot. Additionally the 5 Line Cross Scan can be a perfect solution for a detailed screening and a quick follow-up. Moreover, Topcon's "Enhanced Choroidal Mode*" visualizes further internal structures, allowing much superior visualization of the interface between choroid and sclera. Data analysis is now selectable from 2 formats – Fine or Basic, and can be performed fast or in detail according to your purposes. Experience sophisticated examination with the new evolved TOPCON 3D OCT-2000 Series.

Stunning Retinal Images with the Integrated High Resolution Retinal Camera

Combining OCT and a Color Fundus Camera in one unit, the TOPCON 3D OCT-2000 line-up is perfected now with FA and FAF photography functions. Furthermore digital Red-free images can be displayed easily at one touch. Owing to flexibly changeable ISO sensitivity, reduced flash level with still crystal-clear fundus observation is available. This would result in reduced patient fatigue and miosis. If the OCT image is only required, simply select "Color Photography OFF".

* Choroidal Mode + Overlapping

2009A(ver. 8.00)
Print Date: 2012/01/19
Age: _____
Technician: _____
Fixation: OD(R) Wide
Scan: 3D Wide(12.0 x 8.0mm - 512 x 128)
TOPCON



4 Easy Steps for Operation Flow



1


Register / Select Patient

After Registering / Selecting a Patient, Click on Button


When registering a patient, it is possible to input eye refractive data. Based on the inputted refractive information, the software adjusts the circle diameter for the circle scan and corrects 3D papillary diameter, area and volume, while also calculating magnification compensation, which enables accurate scan performance.

Ex.) Circle Scan Diameter Correction

Axial Myopia



Axial Hyperopia



2

Select scanning pattern

Intuitive Color Touch Panel

The color touch panel allows for easy selection of a scan pattern icon. Scan icons are easily customized from the selection of over 500 different scan patterns, making it possible to create an individual combination of commonly used icons.





3

Capture

Color Fundus / OCT / FA FAF Photography

Auto functions enables anyone to take high quality images with minimal training.


» **Color Fundus / FAF Photography:**


Auto Focus / Auto Shoot

» **OCT Photography:**


Auto Focus / Auto Z & Z Lock / Auto Polarization (Once the retinal position is detected, optimize the output sensitivity of the image and the OCT image will be displayed more clearly.)

» **FA Photography:**


Touch  button simultaneously along with the dye injection. The timer starts to count. The capturing mode can be easily switched to Color Fundus or OCT during FA photography without stopping the timer.



FA photography Starts



Color Fundus



FA photography Stops



OCT



4

Automatic export to PC

View & Analysis


An instant comparison of the enface OCT projection image and the color fundus image is available.

Ex. 3D Disc




Photography Window

Click  and analysis can be implemented.



Analysis Window



Export Report



During Capture




After Capture



PC



3D OCT-2000



Auto Disc Search Function


It automatically detects the disc center referring to the IR image. Tracking the disc center position, it contributes to better accuracy and scan positioning when circle scanning.




Auto Disc Center Detection

After 3D Scanning, the disc center is automatically detected.

- » Measurable scan range
 - 6 × 6mm 512 × 128
 - 4.5 × 4.5mm 256 × 256



Before



After

Motion Correction / Compensation / Rescanning Function

Motion Correction



Corrects the Z direction movement.

Compensation Function

Tracks the ocular and then compensates the X direction movement.

Rescanning Function

Due to Y direction movement, the scanning area may be missed. In such a case, the rescanning function automatically activates.




Auto Fovea Center Detection

The 3D OCT-2000 Series can also automatically center the fovea to ensure accurate reporting and analyzing.

* The function is available in Glaucoma Analysis - Macula


3D OCT-2000 Has Been Designed to Provide Highly Accurate Data



IR Tracking

IR tracking utilizes the IR image during capture to overlap exactly the same place on the retina. If the scanned area moves, the rescanning function automatically begins referring to the first scanned image.

* The number of re-scans can be preset in the range of 0 to 4.



1
Reference Image Displayed


Follow-up Function

The new Follow-up Function defines the scanning location based on the previously captured/selected image and turns **Lock ON** before the next capture, which enables the user to scan the same position under the same conditions.

* The follow-up function is available with Line / 5 Line Cross / Radial Scans.
* It is possible to specify the scanning location manually.



2
Lock ON Before Capture





3
Immediate Comparison

GLAUCOMA & MACULA

» 12×9 mm 3D wide Scan


The capability to capture a wider image from macula to optic disc strongly contributes to better effectiveness of the examination and reduction of patient's fatigue. Macula Analysis as well as Thickness and Significance maps of NFL, GCL+IPL, NFL+GCL+IPL are useful in detecting various macular diseases and glaucoma.



GLAUCOMA

» RNFL Trend Analysis

A maximum of 8 images of both eyes can be displayed on one screen, taking the earliest capture date as a baseline. In addition, checking "register" assures analyzing the same scanned position under the same scanning conditions every time. Color Fundus/ RNFL thickness map/ OCT images/ Cup and Disc ratio can be generated and compared to the normative database.






Fully Comprehensive Analyzed Data

-Easy-to-read & easy-to-understand report templates-

» 3D Disc Report

The below image is an example of a standard glaucoma analysis. The reference plane distance can be customized as 120 μm , 60 μm or 90 μm or any other value in the range of -90~+210 μm .



» Disc Margin


Based on the edges of the RPE of each B-scan, the software automatically determines the disc margin.




» Cup Margin

The software automatically detects the edges of the RPE and using a horizontal line between the RPE edges as a reference point, creates a line 120 μm above it. The Cup Margin will be determined at the cross points of the reference plane and the ILM.

* It is possible to customize the reference plane distance as 120 μm , 60 μm or 90 μm or any other value in the range of -90~+210 μm .

» Glaucoma Analysis - Macula

7 \times 7mm Thickness Map, Significance Map over the Red-free image with a comparison of normative database, Average Data of Superior, Inferior, Total thickness with normative database, and an Asymmetry map which produces a differential value of the superior and inferior thickness can all be shown on one report. By clicking  3D Macula Analysis Report is displayed. Not only this Glaucoma Analysis (Macula) function is useful for glaucoma but for macular disease diagnosis.





Fully Comprehensive Analyzed Data

-Easy-to-read & easy-to-understand report templates-

MACULA

» 3D Macula Report



Comparison to the rich Normative database, Thickness map display, and 3D detail evaluation is available with this 3D analysis.



» Macula Drusen Analysis

Drusen counts are described on the color fundus image and report, and are color-coded according to the Drusen area in each ETDRS grid.

* Drusen Analysis is available only at 3D 6×6mm 512×128 Fine Analysis
 * Drusen defined here counts a suspicious drusen of $\phi 125 \mu\text{m}$ or larger





ANTERIOR

» Anterior Segment Analysis

Corneal thickness map, Corneal thickness distribution diagram, Curvature radius distribution diagram, Curvature radius and Peripheral corneal thickness analysis, Manual angle measurement are all available.

*In order to capture anterior segment photography, it is necessary to use the headrest attachment.




Reference : Dr. Frederique Matonti

FA • FAF


» Import Function

FA, FAF, ICG, Red-free images can be imported easily. The examiner should mark at least 3 spots on a vessel of Color fundus image and the import fundus images, click on the “registration” and the simultaneous observation of OCT and imported images become available.

Case: Central retinal vein occlusion (CRVO)



Import Image





The FA and FAF images were photographed by 3D OCT-2000 FA and 3D OCT-2000 FA Plus model.

3D OCT-2000 in a Hospital Network

IMAGeNet™ Digital Imaging System

The Topcon 3D OCT-2000 and its viewing software plays a powerful role in the management of the patient's data. The unique software enables all patient imaging and data to be collected, saved and reviewed remotely through one unified network, IMAGeNet™. OCT images can be viewed and analyzed through the network at any location; medical meetings, surgical simulation in the operating room and in a patient consultation room. Furthermore, an integrated IMAGeNet™ system allows all clinical images taken throughout the ophthalmology department to be stored in one patient file, thereby facilitating comprehensive diagnosis.





Specifications

Observation & Photography of Fundus Image	
Scan Mode	Color, FA*, FAF**(Spaide Filters), Red-free***
Observation	Near IR
Picture Angle	45° Equivalent 30° (Digital Zoom)
Diopter Scale Range ****	-13D to +12D (in fundus photography)
Operating Distance	40.7mm (in fundus photography) 63.7mm (in anterior segment photography)*****
Photographable Diameter of Pupil	45°: ϕ 4.0mm or more Small pupil diameter: ϕ 3.3mm or more

Observation & Photography of Fundus Image / Anterior Segment Tomogram	
Scanning Range	(On Fundus) [Lateral] within 3~12mm [Vertical] within 3~9mm (On Cornea) [Lateral] within 3~6mm [Vertical] within 3~6mm
Scan Patterns (Recommended)	Macula / Disc: 3D Wide 512x128 (128 horizontal scan lines comprised of 512 A-scans), 12x9mm Macula: 5 Line Cross 1,024x10 (10 scan lines comprised of 1024 A-scans), 9mm ***** Macula: 3D Scan 512x128 (128 horizontal scan lines comprised of 512 A-scans), 6x6mm Macula: Radial Scan 1024x6 or 12 (6 or 12 radial scan lines comprised of 1024 A-scans), 6mm Macula: 7 Line Raster 1024x7 (1024 A-scans per B-scanx7) , 6mm Disc: 3D Scan 512x128 (128 horizontal scan lines comprised of 512 A-scans), 6x6mm Disc: Circle Scan 1024 A-scans, ϕ 3.4mm Anterior: Radial Scan (For Cornea) 1024x12 (12 radial scan lines comprised of 1024 A-scans), 6mm Anterior: Line Scan(For Angle Chamber) 1024 (line scan line comprised of 1024 A-scans), 3mm
Scan Speed	50,000 A-scans per second
Scan Depth	2.3mm
In-depth Resolution	Below 6 μ m
Photographable Diameter of Pupil	ϕ 2.5mm or more

Observation & Photography of Fundus Image / Fundus Tomogram	
Retinal Layers Identified	Macula: ILM, IS/OS, RPE, BM Glaucoma: ILM, NFL, IPL
OCT Reference Focus	Vitreous and Choroid
Fixation	Adjustable internal matrix LCD and external fixation device (Matrix LCD : The display position can be changed and adjusted. The presenting method can be changed.)

Light Source / Power Source / Power Supply	
Light Source	Super luminescence diode(SLD) Wavelength:840nm Half Bandwidth:50nm Output on cornea \leq 0.65 mW
Power Source	Voltage: 100/110/120/220/230/240V Frequency: 50-60Hz
Power Supply	200VA (Max 400VA)

Dimensions / Weight	
Dimensions	545mm(W) \times 535mm(D) \times 600 ~ 630mm(H)
Weight	35kg(3D OCT-2000) 36kg(3D OCT-2000 FA) 37kg(3D OCT-2000 FA plus)

* Only for FA and FA plus models **** Without the diopter compensation
** Only for FA plus model ***** With anterior segment attachment
Display digital Red-free **More variable scan patterns available with a combination of different pixel and scan range.

Courtesy: Yasuyuki Suzuki MD., Department of Ophthalmology, Tokai University Hachioji Hospital
Takeya Kohno MD., Department of Ophthalmology, Osaka City University Hospital

※ Not available for sale in the United States.

*Subject to change in design and/or specifications without advanced notice.

IMPORTANT In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.



TOPCON CORPORATION
75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, Japan. Phone:3-3558-2523/2522 Fax:3-3965-6898 www.topcon.co.jp

TOPCON MEDICAL SYSTEMS, INC.
111 Bauer Drive, Oakland, NJ 07436, U.S.A.
Phone:+1-201-599-5100 Fax:+1-201-599-5250 www.topconmedical.com

TOPCON CANADA INC.
110 Provancher Avenue, Boisbriand, QC J7G 1N1 CANADA
Phone:+1-450-430-7771 Fax:+1-450-430-6457 www.topcon.ca

TOPCON EUROPE MEDICAL B.V.
Essebaan 11, 2908 LJ Capelle a/d IJssel; P.O.Box145;
2900 AC Capelle a/d IJssel; THE NETHERLANDS
Phone:+31-(0)10-4585077 Fax:+31-(0)10-4585045
E-mail: medical@topcon.nl; www.topcon.eu

ITALY OFFICE
Viale dell'Industria 6/0; 20037 Paderno Dugnano; (Milano) ITALY
Phone:+39-02-9166071 Fax:+39-02-91081091 E-mail: topconitaly@iscallit; www.topcon.it

DANMARK OFFICE
Præstemarkvej 25, 4000 Roskilde, DANMARK
Phone:+45-46-327500 Fax:+45-46-327555
E-mail: topcon@topcondanmark.dk www.topcondanmark.dk

IRELAND OFFICE
Unit 276, Blanchardstown, Corporate Park 2 Ballycoolin Dublin 15, IRELAND
Phone:+353-18975900 Fax:+353-18293915 E-mail: medical@topcon.ie; www.topcon.ie

TOPCON S.A.R.L.
BAT A1 3 route de la révolte 93206 SAINT DENIS CEDEX
Tel :+33 1 49 21 23 23 Fax :+33 1 49 21 23 24 E-mail:topcon@topcon.fr; www.topcon.fr

TOPCON DEUTSCHLAND G.m.b.H.
Hanns-Martin-Schleyer Strasse 41; D-47877 Willich, GERMANY
Phone:+49-(0)2154-8850 Fax:+49-(0)2154-885177 E-mail:med@topcon.de; www.topcon.de

TOPCON SCANDINAVIA A.B.
Neogatan 2; P.O.Box 25; 43151 Mölndal, SWEDEN
Phone:+46-(0)31-7109200 Fax:+46-(0)31-7109249 E-mail:medical@topcon.se; www.topcon.se

TOPCON ESPAÑA S.A.
HEAD OFFICE
Frederic Mompou 4 Esc. A Bajos 3, 08960 Sant Just Desvern Barcelona, SPAIN
Phone:+34-93-4734057 Fax:+34-93-4733932 E-mail: medica@topcon.es; www.topcon.es

PORTUGAL OFFICE
Rua da Forte,6-6A,L-0,22,2790-072 Camaxide,PORTUGAL
Phone:+351-210-994626 Fax:+351-210-938786 www.topcon.pt

TOPCON (GREAT BRITAIN) LTD.
Topcon House, Kennel Side, Bone Lane, Newbury, Berkshire RG14 5PX United Kingdom
Phone:+44-(0)1635-551120 Fax:+44-(0)1635-551170 E-mail:medical@topcon.co.uk; www.topcon.co.uk

TOPCON POLSKA Sp. z. o. o.
ul. Warszawska 23; 42-470 Siewierz, POLAND
Phone:+48-(0)32-6705045 Fax:+48-(0)32-6713405 www.topcon-polska.pl

TOPCON SINGAPORE MEDICAL PTE. LTD.
60 Alexandra Terrace, #08-27 The Comtech, Singapore 118502
Phone:+65-68720606 Fax:+65-67736150 E-mail:medical_sales@topcon.com.sg www.topcon.com.sg

TOPCON INSTRUMENTS (MALAYSIA) SDN.BHD.
No. D1, (Ground Floor), Jalan Excella 2, Off Jalan Ampang Putra,
Taman Ampang Hill, 55100 Kuala Lumpur, MALAYSIA
Phone:+60-(0)3-42709866 Fax:+60-(0)3-42709766

TOPCON INSTRUMENTS (THAILAND) CO.,LTD.
77/162 Sinsathorn Tower, 37th Floor, Krungthornburi Rd., Klongtonnai,
Klongsam, Bangkok 10600, THAILAND
Phone:+66(0)2-440-1152-7 Fax:+66-(0)2-440-1158

TOPCON CORPORATION BEIJING OFFICE
Block No.9, Kangding Street Beijing Economic-Technological Development Area,
Beijing, 100176, CHINA
Phone:+86-(0)10-6780-2799 Fax:+86-(0)10-6780-2790

TOPCON CORPORATION SHANGHAI OFFICE
14L Huamin Empire Plaza, No.726, Yan-an Xi Road,
Shanghai, 200050, CHINA
Phone:+86-(0)21-5238-7722 Fax:+86-(0)21-5237-0761

TOPCON CORPORATION BEIRUT OFFICE
P.O.Box 70-1002 Anhelas Beirut, LEBANON
Phone:+961-4-523525/523526 Fax:+961-4-521119

TOPCON CORPORATION DUBAI OFFICE
P.O.Box 293705 Dubai Airport Free Zone L1U, J-12, Dubai, U.A.E
Phone:+971-4-299-5900 Fax:+971-4-299-5901