# STORZ

**KARL STORZ — ENDOSKOPE**

20131520

**Kaltlicht-Fontäne XENON NOVA® 175**

20131520

**Cold light fountain XENON NOVA® 175**

20131520

**Fuente de luz fría XENON NOVA® 175**

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**Instruction Manual**

**Manual de instrucciones**

	20131520	Kaltlicht-Fontäne XENON NOVA® 175
	20131520	Cold light fountain XENON NOVA® 175
	20131520	Fuente de luz fría XENON NOVA® 175

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the XENON NOVA® 175 cold light fountain. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ GmbH & Co. KG

Utilization and reproduction by third parties require express authorization by KARL STORZ GmbH & Co. KG

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**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

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**Controls, displays,  
connectors, and  
their uses**

- ① Main power switch
- ② Lamp replacement indicator
  - glows red when operating hours exceed 450
  - flashes when operating hours exceed 500
- ③ Brightness control knob
- ④ Lamp indicator light
  - lights up blue during light output
- ⑤ Light cable adaptor, interchangeable
- ⑥ Cooling vents
- ⑦ Potential equalization connector
- ⑧ Power fuse holder
- ⑨ Power supply cord receptacle

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 Read the instructions carefully before  
operating the equipment!

 Fiber optic light output

 Brightness

 Lamp

 **Warning!** Burn hazard!

 Only use your hands. Do not use tools!

 Reset of the service hour counter

 Potential equalization connector

 Applied part of type CF

 Environmental protection use period of  
50 years (China RoHS)

 **DANGER:** Risk of explosion if  
used in the presence of flammable  
anesthetics.

 **CAUTION:** To reduce the risk of elec-  
trical shock, do not remove cover.  
Refer servicing to qualified service  
personnel.

**Keep out of reach of patients.  
Do not store liquids on or above the  
unit.**

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Photographs of the unit .....	V
Controls, displays, connectors, and their uses .....	VI
Symbols employed .....	VII

**General information**

Device description .....	2
--------------------------	---

**Safety instructions**

Warnings and cautions .....	3
Intended use .....	6
User qualification .....	7
Safety precautions at the site of installation .....	7
Safety precautions when operating the unit .....	7

**Installation and operating instructions**

Unpacking the equipment .....	8
Basic equipment .....	8
Installing and connecting up the unit .....	8
Operation .....	10
Brightness setting .....	10

**Maintenance**

Lamp replacement .....	11
Fuse replacement .....	14
Cleaning and care of the unit .....	15
Maintenance .....	15
Servicing and repair .....	15
Repair program .....	16
Limitation of liability .....	17
Manufacturer's warranty .....	17

**Technical description**

Troubleshooting .....	19
Technical data .....	20
Technical documentation .....	21
General circuit diagram .....	22

**Spare parts,  
recommended accessories**

Spare parts list .....	23
Accessories .....	23
Cleaning, disinfection and sterilization of fiberoptic/fluid light cables.....	24/26

**Appendix**

Cleaning agents and disinfectants .....	29
Electromagnetic Compatibility (EMC) Information.....	31

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### Device description

The performance specifications of the XENON NOVA® 175 cold light fountain make it suitable for virtually all endoscopic interventions, producing excellent results.

With a color temperature of 6000 K, the 175 watt xenon lamp corresponds to the color temperature of sunlight and therefore produces exceptionally brilliant illumination. Full light intensity is reached as soon as the lamp is switched on.

Brightness is regulated via a mechanically operated dimmer, whilst lamp current remains unchanged. This control method avoids instabilities of the arc and ensures maximum lamp service life.

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**Warnings and cautions**

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words, **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

**Warning:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**Caution:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

**Note:** A Note indicates special information about operating the product, or clarifies important information.

**Warning:** Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

**Note:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

**Note:** Lamps do not have an indefinite service life, so it is possible that one may fail during use. For this reason a spare lamp system or cold light source should always be kept available during therapeutic endoscopic interventions.

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**Warning:** The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

**Warning:** Keep out of reach of patients.

**Warning:** The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

**Warning:** Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

**Warning:** Grounding reliability can only be achieved when the equipment is connected to a properly installed "Hospital Only" or "Hospital Grade" receptacle (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

**Warning:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.

**Warning:** Risk of explosion if used in the presence of flammable anesthetics.

**Warning:** Always unplug the unit before doing any maintenance work on it (e. g. cleaning).

**Warning:** Risk of burns. Lamp and heat sink may be hot. Avoid bodily contact with the lamp until it has cooled sufficiently. Parts of the lamp can reach very high temperatures which cause serious burns if touched.

**Warning:** Danger of electric shock. Do not open the unit. Refer servicing only to qualified personnel. Any opening of the unit by unauthorized persons voids the guarantee. The only exception to this are activities described in this instruction manual.

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**Safety instructions****Warnings and cautions**

**Caution:** Always adjust the cold light fountain to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene, either by direct vision or when coupled to a video camera.

**Caution:** Do not store liquids on or above the unit. Avoid allowing fluids to enter the unit.

**Caution:** Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

**Caution:** Use only fuses of the correct rating.

**Caution:** Do not touch lamp ceramic or window with bare fingers; use gloves or cloth or tissue to handle lamp. Touching the lamp with bare fingers may result in starting difficulty.

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**Light cables**

**Warning:** Danger of glare! Never look into the open end of a light cable.

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**Warning:** High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and at the tip of the endoscope.

**Warning:** Never leave an illuminated light cable with or without an endoscope connected, on or near a patient; the risk of patient burns can result.

**Caution:** Do not place connected light cables on combustible materials such as textiles (drapes) or in the direct vicinity of swabs, etc. soaked with combustible fluids, such as disinfectant. These could catch fire from heat development at the open end of the light cable.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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### **Intended use**

The cold light fountain XENON NOVA® 175 serves as a standard high intensity light source for endoscopy.

Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The XENON NOVA® 175 may only be used with accessories which have been designated by KARL STORZ as suitable for the instrument or the safe use of which is proven.

**Note:** KARL STORZ light cables are recommended for use with this unit. Some manufacturer's cables may be damaged when this unit is operated at higher settings.

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Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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
### User qualification

The cold light fountain XENON NOVA® 175 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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### Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected before power is applied to the unit in compliance with current national regulations.

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### Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using the unit.

During treatment using the cold light fountain XENON NOVA® 175 the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

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### Unpacking the equipment

Carefully unpack the cold light fountain XENON NOVA® 175 and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
### Basic equipment

- 1 Cold light fountain XENON NOVA® 175
- 1 Power cord
- 1 Instruction manual
- 1 Wallet for instruction manual

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### Installing and connecting up the unit

**Note:** The cold light fountain XENON NOVA® 175 may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes.

**Warning:** This unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Set the unit on a flat surface.

**Caution:** Do not block any air vents. Danger of overheating.

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## Installation and operating instructions

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The unit is equipped with a connector ⑦ for attaching an equipotential wire (redundant grounding wire) line.

The unit's equipotential line should be installed by a qualified electrician.

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Connect the power cord.

**Warning:** Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

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Plug the light cable into socket ⑤ as far as it will go. To disconnect the light cable pull on the connector, never on the cable itself.

**Warning:** Danger of glare. Never look into the open end of a light cable.

**Warning:** Danger of burns. Contact with the open end of a light cable connected to a cold light fountain can possibly lead to burns.

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Connect the light cable to the endoscope by a quarter turn of the knurled collar.

**Note:** KARL STORZ light cables are recommended for use with this unit. Some manufacturer's cables may be damaged when this unit is operated at higher settings.

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The light outlet point ⑤ can be unscrewed if necessary and replaced by adapters for light cables from other manufacturers (see section 'Spare parts, recommended accessories', p. 28).

**Warning:** Danger of glare! Do not operate the light source without a light cable adapter.

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### Operation

Switch on main power switch ①.

**Note:** The lamp indicator light ④ lights up blue during light output.

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### Brightness setting

Brightness can be adjusted with the brightness control knob ③.

**Note:** If the cold light fountain is switched on and off frequently, this increases wear to the xenon lamp. Therefore, during short interruptions, brightness should be reduced on the unit.

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### Display: lamp change

The XENON NOVA® 175 cold light fountain has an internal service hour counter.

When the operating hours exceed 450, the lamp change light ② glows red. If the operating hours exceed 500, the lamp change light ② flashes.

The lamp must then be changed (see Lamp Replacement on page 11).

**Caution:** For safety reasons the operating time of the lamp should not exceed 500 hours.

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### Lamp replacement

Switch off the unit and remove the power cord.

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Turn the unit upside down.

Screw out the four outer screws which hold the unit cover in place.

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Turn the unit back up the right way and slide the cover off towards the back.

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**Note:** The electronics of the unit are covered with a sheet of insulating paper for safety reasons.

Fold the insulating paper to the side.

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
**Warning:** Risk of burns. Lamp and heat sink may be hot. Please wait until the lamp/heat sink has cooled down.

Undo the 4 thumbscrews on the heat sink by hand.

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**Warning:** Danger of explosion. The lamp is filled with xenon gas at very high pressure. It must therefore be handled with appropriate care and must only be replaced by skilled electrical engineers. High mechanical forces could cause the lamp to break.

**Wear safety goggles!** We urgently recommend wearing protective clothing (face protection, gloves, etc).

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
Thoroughly clean the complete heat sink (upper and lower halves) with cotton-tipped applicators and spirit or similar.

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To ensure that the lamp has direct contact with the heat sink, it is necessary to use thermal compound.

**Caution:** The application of the thermal compound must be very sparing. Too much thermal compound can cause shortened lamp life and ignition problems.


Using a cotton-tipped applicator, apply a thin film of thermal compound to all the points where the lamp comes into contact with the heat sink.

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**Caution:** Only use the genuine replacement lamp 20132026 recommended by KARL STORZ. Do not touch the new lamp with bare fingers; use gloves or a cloth or tissue to handle the lamp. Touching the lamp with bare fingers may result in operational problems.

Insert new lamp into the bottom halves of the heat sink. Fit the upper halves of the heat sink and tighten the thumbscrews by hand uniformly.

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Press the reset pin to set the lamp's operating hour counter to zero.

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Slide the cover of the unit back on and turn the unit over.

Screw the screws back into the bottom of the unit.

Reconnect power supply.

Test the instrument for proper operation.

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**Fuse replacement**

Switch off the light source and remove the power cord.

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Remove the power fuse holder ⑧ with a screwdriver or other tool.

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**! Caution:** Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

**Note:** New fuse ratings for XENON NOVA® 175 as of serial number BC0627729!

	220–240 VAC	100–125 VAC
Power fuse	2 x T 1.6 A L	2 x T 3.15 A L

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Refit the power fuse holder ⑧.  
Reconnect the power cord.  
Test the unit for proper operation.

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### Cleaning and care of the unit

**Warning:** Always pull out power plug before cleaning! Make sure the lamp has cooled down.

**Caution:** Avoid allowing fluids to enter the unit.

The exterior of the unit must be cleaned with a disposable cloth moistened with disinfectant (no solutions predominantly containing alcohol).

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
### Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year. For detailed instructions please refer to the version of the Service Manual currently valid, catalog number SV 3367 (English version).

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### Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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**Disposal**

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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**Repair program**

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str. 34  
78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

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**Important information**

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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### Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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### Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen, Germany

**Note:** The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG  
Abt. Reparaturservice  
Dr. Karl-Storz-Str. 34  
78532 Tuttlingen, Germany

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**Warranty for xenon lamps**

Independent of the general warranty conditions for our instrument, we guarantee 500 hours of life for our xenon lamps. Lamps that burn out within the first 50 hours will be replaced free of charge; for lamps that burn out after the first 50 hours, the price of the replacement lamp will be calculated according to the number of hours the original lamp was used. This is applicable only if the defective lamp, the serial number of the XENON NOVA® 175 cold light fountain, and a record of the lamp hours are presented.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall release us from any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.\*

\* The only exception to this are activities described in this instruction manual.

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**Troubleshooting**

**Warning:** Always unplug the unit before all maintenance work.

*Symptom:*  
Complete failure of the unit.

*Possible causes:*

- No power from the power line.
- Defective power fuse.

*Remedy:*

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Make sure you use the correct type of fuse.

*Symptom:*  
No light is emitted.

*Possible causes:*

- Xenon lamp is defective.
- Defect in internal circuitry.

*Remedy:*

- Have a trained electrician change the xenon lamp.
- Call service.

**Note:** Observe the warranty conditions!

*Symptom:*  
Too little light.

*Possible causes:*

- Brightness adjustment knob is set too low.
- Lamp is defective.
- Ends of light cable and/or endoscope are soiled.
- Light cable is defective.

*Remedy:*

- Turn brightness control knob to the right.
- Have xenon lamp replaced.
- Clean ends of light cables, and light entrances and exits on the endoscope.
- Replace light cable.

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- Type of protection against electric shocks:  
Protection Class I
- Degree of protection against electric shocks:  
Applied part of type CF
- Type of protection against moisture:  
drip-water protection as per IPX 1

**According to IEC 60601-1-2:**

Please read the Electromagnetic Compatibility  
Information in the appendix (p. 31-43).

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**Directive compliance** (for 20 1315 20)

According to Medical Device Directive (MDD):

This medical product belongs to Class I

This medical product bears the CE mark in  
accordance with MDD 93/42/EEC.

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### Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the instrument that have been designated as repairable by their respective manufacturers.

Supply of such technical documentation relating to the instrument shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the instrument.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

*We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.*

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
**Accessories**

**Fiberoptic light cable**

Fiberoptic light cables contain a large number of single glass fibers. The spectral transmission of a fiberoptic light cable is approximately the same for all of the wavelengths of visible light. It decreases, however, towards the blue end of the spectrum. This causes the light to take on a warmer shade – an effect which becomes more intensive the longer the light cable. Fiberoptic light cables are opaque to ultraviolet light.

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**Cleaning, disinfection and sterilization of fiberoptic light cables**

**Note:** Comprehensive information on cleaning light cables can be found in the 'Light cable' instructions.

1. **Manual (preliminary) cleaning:** The light cables may be cleaned using a cleaning agent/disinfectant which is specially approved for endoscopes. Wipe the outside of the light cables with a disinfected sponge or disposable cloth soaked in cleaning solution.\*

2. **Disinfection:** immerse in a disinfection solution (see Appendix, page 29, 30).\*

\* re 1./2. The light cable must always be rinsed until it is free of all residue after it has been cleaned or disinfected manually. Microbiologically pure/sterile water must be used for the final rinsing.

3. **Thermal decontamination by machine** (up to 93°C)

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4. **Manual final cleaning:** Clean the fiber surfaces and optical end faces with a soft disposable cloth, disinfected sponge or cotton swab, soaked with 70% alcohol. Stubborn contaminants and deposits may be removed using cleaning paste (Art. no. 27661).

5. **Sterilization:** Steam sterilization at 134°C (+3°C) STERIS SYSTEM1® processor or gas or plasma sterilization.

**Note:** The instructions for 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.

Use a cleaning agent/disinfectant approved by KARL STORZ for cleaning/disinfection (see Appendix). You will also find the current list of approved products on our Web site ([www.karlstorz.com](http://www.karlstorz.com)).

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**Fiberoptic light cable**

**Combination of light cable with endoscope**

For safety reasons, the diameters of fiberoptic/fluid\* light cables and endoscopes should be used in the combination shown in the table opposite.  
\* see page 26

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

**Fluid light cables**

In fluid light cables the light is not transmitted through glass fibers, but through a special liquid contained in the light cable. Fluid light cables are stiffer than glass fiber light cables and bending them too much will cause a severe loss of light. However, they do transmit a brighter light than fiberoptic light cables with the same diameter. For this reason, fluid light cables are ideally suited for endoscopic photography. Fluid light cables are permeable to ultraviolet light.


**Caution:** Damage to the elastomer coat will impair the light cable.

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**Cleaning, disinfection and sterilization of fluid light cables**

1. Manual (preliminary) cleaning: The fluid light cables may be cleaned using a cleaning agent/disinfectant which is specially approved for endoscopes. Remove heavy soiling with a moist disinfected sponge or lint-free disposable cloth.\*

2. Disinfection: immerse in a disinfection solution (see Appendix, page 29, 30).\*

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3. Machine decontamination

**Caution:** Carry out chemical disinfection by machine at max. 65°C.

**Caution:** When preparing light cables in a cleaning and disinfection machine, particular care is necessary in the choice of chemical agents and machine programs.

4. Manual final cleaning: Clean the fiber surfaces and optical end faces with a soft disposable cloth, disinfected sponge or cotton swab, soaked with 70% alcohol.

Stubborn contaminants and deposits may be removed using cleaning paste (Art. no. 27661).

5. Sterilization: Gas sterilization or STERIS SYSTEM1® processor.

**Caution:** Fluid light cables must not be steam or plasma sterilized under any circumstances.

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**Adaptors**

Adaptors for light cable from	Cat. no.
ACMI	487 A
Machida	487 M
Olympus	487 O
Pilling	487 P
KARL STORZ	2108191
KARL STORZ (adaptor with integrated heat sink)	<b>20 1340 80</b>
KARL STORZ, Wolf, Dyonics, Circon-ACMI, Pilling, V. Müller, Eder, Stryker, KLI, Sparta	487 U
KARL STORZ, Wolf, Dyonics, Circon-ACMI, Olympus, V. Müller, Eder, Stryker, KLI, Sparta	487 UO
KARL STORZ video endoscopes with external insufflation connector	487 VE
Wolf	487 W
Winter-IBE	487 WI
Zeiss	487 Z

More accessories are listed in our TELEPRESENCE catalogue.

For detailed information please contact:  
KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen, Germany

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**Typ**

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)
- (A) = Only for the flexible endoscopes of the new generation (year of manufacture > 2005)

**Note:** You will also find the current list of approved products on our Web site ([www.karlstorz.com](http://www.karlstorz.com)).

**Note:** When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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**Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets**  
**Manual preparation**

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Ecolab GmbH & Co. OHG	Sekusept Plus	1, 2, 3
	Sekusept Pulver Classic	1, 2, 3
	Ultrasdesmit AF	1, 2, 3
Esteer Pharma GmbH	Ultrasreptin Aktiv	1, 2, 3
	Ultrasreptin Classic	1, 2, 3
	Afid	1, 2
Fresenius AG	Afid plus	2
	Sporcid FF	1, 2
Holifa Polska Sp. zo.o	Poisept Holifa	1, 2
Indeba Industria E Comercio LTDA	Letahaldeido	2
Jose Collado S.A.	Darodor 4000 Liquido	1
	Darodor 9000	2
	Darodor Sinaldehyd 2000	2
	Aldasan 2000	1, 2
Lysoform	Almyrol	1, 2
Dr. Rosemann GmbH	Desoform	1, 2, 3
	Lysoformin 3000	1, 2
	MedDis	2
Medichem International	MediZyme	1
	Edisonite Super	1, 3
	Mucadont-IS	2, 3
Merz Hygiene GmbH	Mucadont-Zymaktiv	1, 3
	Mucocit-T neu	1, 2, 3
METASYS Medizintechnik GmbH (formerly LONZA) Orochemie	ID 50	2
Promagent AB	A 20	1, 3, A
	Wavacide	2
	Descoton Extra	1, 2
	Descoton forte	2
Schumacher, Dr. GmbH	Perfektan Endo	2, 3
	Perfektan Neu	1, 2
	Perfektan TB	1, 2, 3, A
	Gigasept	2, 3
	Gigasept AF forte	1, 2
	Gigasept FF	1, 2
Schülke & Mayr GmbH	Gigasept Instru AF	1, 2, 3, A
	Gigasept Med forte	1, 2
	Gigazyme	1, 3
	Lysetol V	2
Steris	EnzyCare 2	1, 3
	neodisher LM 2	1, 3
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
Weigert, Dr. GmbH & Co.	neodisher medizym	1, 3
	neodisher SeptoClean	1, A
	neodisher SeptoPreClean	1, 2
Whiteley Industries PTY.Ltd	Aidal Plus	2
	Matrix	1

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**\* Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set.

For flexible endoscopes, only immerse the examination sheath, not the housing/handle part.

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**Type**

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)

**Note:** You will also find the current list of approved products on our Web site ([www.karlstorz.com](http://www.karlstorz.com)).

**Note:** When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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**Machine preparation**

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles the appropriate manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schülke & Mayr GmbH	Thermosept alka clean	1
	Thermosept DK	2
	Thermosept ED	2
	Thermosept ER	1
	Thermosept RKN-zym	1
Steris	Instru-Klenz	1, 3
	neodisher DuoClean	1, A
Weigert, Dr. GmbH & Co.	neodisher FA	1
	neodisher FA forte	1, 3
	neodisher Oxivario	1
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher mediklar	5
	neodisher medizym	1
	neodisher SeptoClean	1, 2, A
	neodisher Septo DA	2
	neodisher Septo DN	2

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**Appendix****Electromagnetic Compatibility  
(EMC) Information**

**WARNING:** Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix, as well as the instruction manual, during installation and commissioning.

The XENON NOVA® 175 model 20131520 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The XENON NOVA® 175 model 20131520 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

**NOTE:** The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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**WARNING:** The XENON NOVA® 175 model 20131520 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**WARNING:** The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

**WARNING:** The use of cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the XENON NOVA® 175 model 20131520. The cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200				
Cables which have been shown to comply with EN/IEC 60601-1-2:				
Cable type	Shielded	Length [m]	Ferrite	Used for
AC Power cord	Yes	3	No	Power supply

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Table 201 Guidance and manufacturer's declaration – electromagnetic emissions		
<p>XENON NOVA® 175 model 201315 20 is intended for use in the electromagnetic environment specified below. The customer or user of the XENON NOVA® 175 model 201315 20 should ensure that it is used in such an environment.</p>		
Emissions test	Compliance	Electromagnetic environment - guidelines
RF emissions CISPR 11	Group 1	XENON NOVA® 175 model 201315 20 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The XENON NOVA® 175 model 201315 20 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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<b>Table 202</b>			
<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
<p>XENON NOVA® 175 model 20 1315 20 is intended for use in the electromagnetic environment specified below. The customer or user of the XENON NOVA® 175 model 20 1315 20 should ensure that it is used in such an environment.</p>			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T^*$ (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Complies <5% $U_T^*$ (>95% dip in $U_T$ ) for 0.5 cycle Complies 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles Complies 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles Complies <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: $U_T$ is the a.c. mains voltage prior to application of the test level.			

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Table 204 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
XENON NOVA® 175 model <b>20</b> 1315 20 is intended for use in the electromagnetic environment specified below. The customer or user of the XENON NOVA® 175 model <b>20</b> 1315 20 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V<sub>rms</sub> 150 kHz to &lt; 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V<sub>rms</sub></p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the XENON NOVA® 175 model <b>20</b> 1315 20, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:</p> <p><math>d = [3.5/3]\sqrt{P}</math> 150 kHz to &lt; 80 MHz</p> <p><math>d = [3.5/3]\sqrt{P}</math> 80 MHz to &lt; 800 MHz</p> <p><math>d = [7/3]\sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.</p> <p><b>a</b> Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.</p> <p><b>b</b> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than <b>3</b> V/m.</p>			

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Table 206 Recommended separation distances between portable and mobile RF communications equipment and the XENON NOVA® 175 model 20 1315 20			
The XENON NOVA® 175 model 20 1315 20 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance <i>d</i> [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[ \frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance <i>d</i> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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# WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

KARL STORZ GmbH & Co. KG  
Postfach 230  
78503 Tuttlingen/Germany

To be completed by the supplier/importer:  
Company stamp/signature:

Remittente/Sello de la empresa:

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Tipo del aparato: Número de serie:

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Firma/Fecha:

To be filled out by instrument owner:  
Return address/company stamp:

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KARL STORZ Endoscopy China Ltd.  
Hong Kong Representative Office  
Unit 1601, Chinachem Exchange Square  
1 Hoi Wan Street, Quarry Bay,  
Hong Kong, People's Republic of China  
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Beijing Representative Office  
Room 610, China Life Tower  
No. 6, Chaowai Street  
Beijing, 100020, People's Republic of China  
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Unit 3901-3904, Tower 1 Grand Gateway  
No. 1 Hong Qiao Road  
Shanghai, 200030, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Chengdu Representative Office  
F-5, 24/F., Chuanxing Mansion,  
No. 18 Renming Road South  
Chengdu, Sichuan, 610016, People's Rep. of China  
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.  
Guangzhou Representative Office  
Room 1119-20, Dongshan Plaza  
69 Xianlie Road Middle, Dongshan District,  
Guangzhou, Guangdong, 510095,  
People's Rep. of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd  
3791 Jalan Bukit Merah  
06-11 e-Centre @ Redhill  
Singapore 159471, Singapore  
Tel. No. +65 63761066, Fax. No. +65 63761068  
Email : infoasia@karlstorz.com.sg  
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
3791 Jalan Bukit Merah  
#06-07 e-Centre @ Redhill  
Singapore 159471, Singapore  
Tel. No. +65 65325548, Fax No. +65 65323832  
Email : infogsp@karlstorz.com.sg  
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.  
Bois Hongo Building 6Fl, 3-42-5 Hongo  
Bunkyo-ku, Tokyo 113-0033, Japan  
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988  
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ GmbH & Co. KG**

Mittelstraße 8  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.de](mailto:info@karlstorz.de)  
Web: [www.karlstorz.com](http://www.karlstorz.com)

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All rights reserved.