STORZ XENON 300 INSTRUCTION MANUAL

Quick Links

[Controls, Displays, Connectors, and Their Uses](#)

[Service Mode](#)

[Lamp Module Replacement](#)

[Reset the Lamp's Operating Hour Counter](#)

[Troubleshooting](#)

Table of Contents

- Controls, displays, connectors, and their uses
- Symbols employed
- Table of Contents
- Photographs of the unit
- General information
- Device description
- Safety instructions
- Warnings and cautions
- Intended use
- User qualification
- Safety precautions at the site of installation
- Safety precautions when operating the unit
- Installation and operating instructions
- Unpacking the equipment
- Basic equipment
- Installing and connecting up the unit
- Operation
- Brightness control
- Standby operation
- Operation of the antifog pump
- Service mode
- Maintenance
- Lamp module replacement
- Reset the lamp's operating hour counter
- Fuse replacement
- Cleaning and care of the unit
- Maintenance
- Servicing and repair
- Repair program

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**Important information
for users of KARL STORZ
instruments**

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the XENON 300 cold light fountain. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

Note: Always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene, whether by direct vision or when coupled to a video camera.

Please be aware of the heat energy generated at the distal tip of the endoscope. The higher the intensity setting of the light source, the more heat will be generated at the endoscope tip. Never leave an illuminated light cable with or without an endoscope connected on or near a patient or surgical drape: the risk of patient burns and igniting flammable material can result.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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
Controls, displays, connectors, and their uses

- ① Power switch
- ② Lamp replacement indicator
- ③ Warning light: lamp service life
- ④ Standby button, with pilot lamp
- ⑤ Bar graph of brightness
- ⑥ ± buttons for brightness control
- ⑦ Digital readout
- ⑧ M button
- ⑨ Lamp indicator light
– lights up blue during light output
- ⑩ Antifog pump ON/OFF switch
- ⑪ Antifog air pump output
- ⑫ Light exit, unscrewable
- ⑬ Cooling vents
- ⑭ SCB connectors*
- ⑮ Potential equalization connector
- ⑯ Power supply cord receptacle
- ⑰ Power fuse holder
- ⑱ Lamp module
- ⑲ Quarter-turn fasteners

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* The KARL STORZ-SCB® interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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	Read the instructions carefully before operating the equipment!
	Standby mode
	Manual brightness control mode
	Antifog air pump output
	Antifog air pump ON/OFF switch
	Light output
	Lamp wattage
	Lamp life warning light
	Replacement indicator
	Potential equalization connector
	Applied part of type CF
	Environmental protection use period of 50 years (China RoHS)

	DANGER: Risk of explosion if used in the presence of flammable anesthetics.
	CAUTION: To reduce the risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel. Keep out of reach of patients. Do not store liquids on or above the unit.

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Photographs of the unit	A 1
Controls, displays, connectors, and their uses	A 2
Symbols employed	A 3

General information

Device description	2
--------------------------	---

Safety instructions

Warnings and cautions	3
Intended use	6
User qualification	7
Safety precautions at the site of installation	7
Safety precautions when operating the unit	7

Installation and operating instructions

Unpacking the equipment	8
Basic equipment	8
Installing and connecting up the unit	8
Operation	11
Brightness control	11
Standby operation	12
Operation of the antifog pump	12
Service mode	13

Maintenance

Lamp module replacement	14
Reset the lamp's operating hour counter	15
Fuse replacement	16
Cleaning and care of the unit	17
Maintenance	17
Servicing and repair	17
Repair program	18
Limitation of liability	19
Manufacturer's warranty	19

Technical description

Troubleshooting	21
Technical data	22
Technical documentation	23
General circuit diagram	24

Spare parts, recommended accessories

Spare parts list	25
Accessories	26

Appendix

Cleaning agents and disinfectants	31
Electromagnetic Compatibility (EMC) Information	33

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Device description

The performance specifications of the XENON 300 cold light fountain make it suitable for virtually all endoscopic interventions, producing excellent results, especially for photographic and video documentation.

With a color temperature of 6000 K, the 300 W xenon lamp corresponds to the color temperature of sunlight and therefore produces exceptionally brilliant illumination. Full light intensity is reached as soon as the lamp is switched on.

The brightness is regulated via a microprocessor controlled optomechanical dimmer, while the lamp current remains unchanged. This control method avoids instabilities of the arc and ensures maximum lamp service life.

The brightness can be regulated manually or automatically via the output signal of a video camera.

In addition, an antifog air pump is available for endoscopes which have a special antifog channel to prevent the lens from misting up.

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Function of the warning and indicator lights on the front panel

Lamp life warning light ③

With a lamp service life of 450–499 hours (450 hours corresponds to 90% of the guaranteed lamp service life), the lamp life warning light ③ is lit on the front panel during operation of the cold light fountain.

Replacement indicator ②

With a xenon lamp service life of 500 hours or more, the replacement indicator ② is lit on the front panel. This indicator light only goes out when the service hour counter for the xenon lamp is set to 0 (see section "Lamp module replacement" pages 14 and 15).

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words, **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

Note: A Note indicates special information about operating the product, or clarifies important information.

Warning: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

Note: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

Note: Lamps do not have an indefinite service life, so it is possible that one may fail during use. For this reason a spare lamp system or cold light source should always be kept available during therapeutic endoscopic interventions.

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Warning: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

Warning: Keep out of reach of patients.

Warning: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

Warning: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

Warning: Grounding reliability can only be achieved when the equipment is connected to a properly installed "Hospital Only" or "Hospital Grade" receptacle (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

Warning: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

Warning: Risk of explosion if used in the presence of flammable anesthetics.

Warning: Always unplug the unit before doing any maintenance work on it (e. g. cleaning).

Warning: Avoid bodily contact with the lamp until it has cooled sufficiently. Parts of the lamp can reach very high temperatures which cause serious burns if touched.

Warning: Danger of electric shock. Do not open the unit. Refer servicing only to qualified personnel. Any opening of the unit by unauthorized persons voids the guarantee. The only exception to this are activities described in this instruction manual.

Warning: The antifog air pump is used in conjunction with special sheaths, as for example in the fields of bronchoscopy, esophagoscopy and laryngoscopy, solely for ventilation of the object lens to prevent it from fogging up. It is in no way suitable for insufflation or lung ventilation.

Warning: Turn off antifog air pump when lens is cleaned. Do not allow air to accumulate inside the patient. The antifog air pump can dilute anesthesia. This antifog air pump is intended for connection to the corresponding antifog air port on the endoscope. Do not use for insufflation.

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Caution: Always adjust the cold light fountain to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene, either by direct vision or when coupled to a video camera.

Caution: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

Caution: Before connecting the unit to the electrical supply lines, verify that the line voltage stated on the unit's identification plate agrees with that of the electrical supply lines to be used.

Caution: Use only fuses of the correct rating.

Caution: Do not touch lamp ceramic or window with bare fingers; use gloves or cloth or tissue to handle lamp. Touching the lamp with bare fingers may result in starting difficulty.

Caution: Federal law restricts this device to sale by or on the order of a physician (USA only).

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Light cables

Warning: Danger of glare. Never look into the open end of a light cable.

Warning: Danger of burns. Contact with the open end of a light cable connected to a cold light fountain can possibly lead to burns.

Warning: High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and at the tip of the endoscope.

Warning: Never leave an illuminated light cable with or without an endoscope connected, on or near a patient; the risk of patient burns can result.

Caution: Do not place connected light cables on combustible materials such as textiles (drapes) or in the direct vicinity of swabs, etc. soaked with combustible fluids, such as disinfectant. These could catch fire from heat development at the open end of the light cable.

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The cold light fountain XENON 300 model 20 1331 20-1 serves as a light source for use with rigid and flexible endoscopes in all types of endoscopic interventions. It is fitted with an antifog pump which can be used in conjunction with suitable endoscopes to prevent the lens from fogging. Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The cold light fountain XENON 300 may only be used with accessories which have been designated by KARL STORZ as suitable for the unit or the safe use of which is proven.

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Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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
User qualification

The cold light fountain XENON 300 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations. It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected before power is applied to the unit in compliance with current national regulations.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using the unit.

During treatment using the cold light fountain XENON 300 the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring sterile application conditions where required by the type of intervention.

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Unpacking the equipment

Carefully unpack the cold light fountain XENON 300 and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
Basic equipment

- 1 Cold light fountain XENON 300
- 1 Silicone tube 610 AFT, length 250 cm
- 1 SCB connecting cord 200901 70
- 1 Power cord 400 A
- 1 Instruction manual
- 1 Wallet for instruction manual

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Installing and connecting up the unit

Note: The cold light fountain XENON 300 may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes.

Warning: This unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Fitting a label strip to the front panel

Undo the left-hand front side cover with a screwdriver.

Note: The appropriate controls/indicators on the unit are designated on the label strip.

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Remove any unmarked label strips. Select a label strip and slide it between the plastic and the metal front.

Fold the strip along the marked line so that the side cover can be refitted without any obstruction.

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Set the unit on a flat surface.

Caution: Do not block any air vents. Danger of overheating.
Adequate cooling is required for proper operation of the unit.

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The unit is equipped with a connector 15 for attaching an equipotential wire (redundant grounding wire) line.
The unit's equipotential line should be installed by a qualified electrician.

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Connect power cord. Insert power cord into power cord receptacle 16 as far as it will go.

Warning: The unit should only be connected to or disconnected from the power supply outside of areas subject to explosion hazards.

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Connect the light cable to the endoscope by a quarter turn of the knurled collar.

Note: KARL STORZ light cables are recommended for use with this unit. Some manufacturer's cables may be damaged when this unit is operated at higher settings.

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KARL STORZ-SCB®

Note: To avoid mistakenly pulling out the SCB connecting cable, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets ⑭. Connect the other end of the cable to the KARL STORZ-SCB® (KARL STORZ Communication Bus) control unit or other SCB units (see KARL STORZ-SCB® System Instruction Manual).

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Operation

Switch on main power switch ①.

Immediately after switching on, the accumulated number of operating hours of the xenon lamp to date will be shown on digital display ⑦.

Note: The lamp indicator light ⑨ lights up blue during light output.

Note: The light source should never be used in surgical procedures requiring direct illumination of the eye.

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Brightness control

The unit will initialize at 5% brightness (XENON 300 as of serial number BE3660).

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The brightness can be changed upwards or downwards using the ± buttons ⑥ in 1% increments (adjustment range 5-100%).

The output brightness is additionally indicated in the bar graph ⑤.

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Standby mode

By pressing the standby button (4), the light output will be reduced to a minimum value as such, the temperature of the equipment attached to the light source will reduce to near ambient (room) temperature and will remain at this level while still in this mode of operation.

Standby mode is indicated by the pilot lamp in the standby button (4) lighting up and the bar graph (5) disappearing.

If the standby button (4) is pressed again, the light output returns to the previously set brightness.

Pressing the + button (6) also ends standby mode. The light intensity can then be newly adjusted.

Note: If the cold light fountain is switched on and off frequently, this increases wear to the xenon lamp. Therefore, during short interruptions, the unit should be left in standby mode.

Operation of the antifog pump

! Warning: The antifog air pump is intended for antifog and lens cleaning applications only. It is not intended for insufflation or ventilation. Turn off antifog air pump when lens is cleaned. Do not allow air to accumulate inside the patient. Antifog air may dilute anesthesia.

Connect endoscopes which have an integrated antifog channel to the pump output (11) via an air tubing.

Press the air pump button (10) to toggle the air pump on or off.

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Service mode

The "Service Displays" mode allows various information concerning the operational status of the cold light fountain to be interrogated.

To activate "Service Displays" mode, the Standby button ④ must be held down for 8 seconds. The number of operating hours of the xenon lamp then appears on digital display ⑦.

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If the + button ⑥ is pressed repeatedly, the following information can be called up in sequence shown below:

- Operating hours of the xenon lamp
- Lamp starts (ignitions)/10
- Operating hours of the unit as a whole/10
- Current operating temperature in °C
- Highest operating temperature reached to date in °C
- Error message codes E 1 – E 15
- ROM version number

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By pressing the – button ⑥ the displays can be called up in reverse sequence.

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To leave this operating mode, hold down the standby button ④ again for 8 seconds or switch off the main power switch ①.

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Lamp module replacement

Warning: Switch off the unit. Disconnect
the power plug.

Undo the quarter-turn fasteners ⑨ (1/4 turn
counterclockwise) and pull the lamp door out of
the unit.

Warning: Danger of burns. The lamp
and lamp module may be very hot.

Warning: Danger of explosion.
The lamp is filled with xenon gas at very
high pressure. It must therefore be han-
dled with appropriate care and must
only be replaced by skilled electrical
engineers. High mechanical forces
could cause the lamp to break.

Wear safety goggles! We urgently rec-
ommend wearing protective clothing
(face protection, gloves, etc).

Pull both card extractors on the lamp module
simultaneously outwards. Open the lamp mod-
ule lock and pull the lamp module assembly
straight out of the unit.

Insert new lamp module assembly and fasten
the card extractors.

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Insert, close and refasten the lamp door by a quarter turn of the fasteners ⑱ clockwise.

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Write down the number of the lamp hours and keep a record as to the lamp's service life (see guarantee information). Then reset lamp hours to zero. To do so, press and hold the standby-button ④ for 8 seconds. This will put the display into service mode.

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Press the standby-button ④ and the – button ⑥ simultaneously for 8 seconds. This will reset the lamp hours to zero.

Note: Counting of the hours run by the unit as a whole will not be affected by this.

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Note: The defective lamp module can be sent back to KARL STORZ for disposal or for repair. It must be shipped in the original packaging.

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Fuse replacement

Switch off the unit and remove the power cord.

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Remove the power fuse holder ⑰ with a screwdriver or other tool.

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! Caution: Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

Note: New fuse ratings for XENON 300 as of serial number BC0627916!

	220...240 VAC	100...125 VAC
Power fuse	2 x T 2.5 AL250V	2 x T 5.0 AL250V

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Refit the power fuse holder ⑰.
Reconnect the power cord.
Test the instrument for proper operation.

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Cleaning and care of the unit

Warning: Always pull out power plug before cleaning! Make sure the lamp has cooled down.

Caution: Avoid allowing fluids to enter the unit.

The exterior of the unit can be surface disinfected with a disposable cloth soaked in disinfectant (preferably alcohol-free).

Clean the antifog tubing and connection parts in a cleaning solution before immersing them in disinfectant solution, making sure that all air bubbles have been removed, and finally rinsing them off and through with sterile or completely demineralized water.

We recommend machine preparation at 93°C.

Note: Use a cleaning agent approved by KARL STORZ for cleaning. You will find the current list of approved products on our Web site (www.karlstorz.com).

Sterilization: preferably the fractionated prevacuum method at 134°C (+3°C).

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
Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year. For detailed instructions please refer to the version of the Service Manual currently valid, catalog number SV 3541 (English version).

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

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Manufacturer's warranty

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply. Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen/Germany

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen/Germany

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Warranty for xenon lamp or lamp module

Independent of the general warranty conditions for our equipment, we guarantee 500 hours of life for our xenon lamps or lamp modules. Lamps or lamp modules that burn out within the first 50 hours will be replaced free of charge; for lamps or lamp modules that burn out after the first 50 hours, the price of the replacement lamp or module will be calculated according to the number of hours the original lamp was used. This is applicable only if the defective lamp/lamp module, the serial number of the xenon cold light fountain, and a record of the lamp hours are presented.

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Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall release us from any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.*

* The only exception to this are activities described in this instruction manual.

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Troubleshooting

Warning: Always unplug the unit before all maintenance work.

Symptom:
Complete failure of the unit.

Possible causes:

- No power from the power line.
- Defective power fuse.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Make sure you use the correct type of fuse.

Symptom:
No light is emitted.

Possible causes:

- xenon lamp is defective.
- Defect in internal circuitry.

Remedy:

- Have a trained electrician change the xenon lamp.
- Call service.

Note: Observe the warranty conditions!

Symptom:
Too little light.

Possible causes:

- Brightness adjustment is set too low.
- Lamp is defective.
- Ends of light cable and/or endoscope are soiled.
- Light cable is defective.

Remedy:

- Press brightness control „+“ button ⑥.
- Have xenon lamp replaced.
- Clean ends of light cables, and light entrances and exits on the endoscope.
- Replace light cable.

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Standard compliance (for 201331 20-1)
**According to IEC 60601-1, IEC 60601-2-18,
UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:**

- Type of protection against electric shocks:
Protection Class I
- Degree of protection against electric shocks:
Applied part of type CF 
- Type of protection against moisture:
drip-water protection as per IPX 1

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility
Information in the appendix (p. 33-45).

Directive compliance (for 201331 20-1)

According to Medical Device Directive (MDD):

This medical product belongs to Class II a

This medical product bears the CE mark in
accordance with MDD 93/42/EEC.

Note: The code number after the CE mark
indicates the responsible notified body.

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the instrument that have been designated as reparable by their respective manufacturers.

Supply of such technical documentation relating to the instrument shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the instrument.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
Accessories

Fiberoptic light cable

Fiberoptic light cables contain a large number of single glass fibers. The spectral transmission of a fiberoptic light cable is approximately the same for all of the wavelengths of visible light. It decreases, however, towards the blue end of the spectrum. This causes the light to take on a warmer shade – an effect which becomes more intensive the longer the light cable. Fiberoptic light cables are opaque to ultraviolet light.

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Cleaning, disinfection and sterilization of fiberoptic light cables

Note: Comprehensive information on cleaning light cables can be found in the 'Light cable' instructions.

1. **Manual (preliminary) cleaning:** The light cables may be cleaned using a cleaning agent/disinfectant which is specially approved for endoscopes. Wipe the outside of the light cables with a sponge or cloth soaked in cleaning solution.*

2. **Disinfection:** immerse in a disinfection solution (see Appendix, page 31, 32).*

* re 1./2. The light cable must always be rinsed until it is free of all residue after it has been cleaned or disinfected manually. Microbiologically pure/sterile water must be used for the final rinsing.

3. **Thermal decontamination by machine** (up to 93°C)

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4. **Manual final cleaning:** Clean the fiber surfaces and optical end faces with a soft cloth, sponge or cotton swab, soaked with 70% alcohol. Stubborn contaminants and deposits may be removed using cleaning paste (Art. no. 27661).

5. **Sterilization:** Steam sterilization at 134°C (+3°C) STERIS SYSTEM1® processor or gas or plasma sterilization.

Note: The instructions for 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.

Use a cleaning agent/disinfectant approved by KARL STORZ for cleaning/disinfection (see Appendix). You will also find the current list of approved products on our Web site (www.karlstorz.com).

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Fiberoptic light cable

Combination of light cable with endoscope

For safety reasons, the diameters of fiberoptic/fluid* light cables and endoscopes should be used in the combination shown in the table opposite.
* see page 28

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Fluid light cables

In fluid light cables the light is not transmitted through glass fibers, but through a special liquid contained in the light cable. Fluid light cables are stiffer than glass fiber light cables and bending them too much will cause a severe loss of light. However, they do transmit a brighter light than fiberoptic light cables with the same diameter. For this reason, fluid light cables are ideally suited for endoscopic photography. Fluid light cables are permeable to ultraviolet light.


Caution: Damage to the elastomer coat will impair the light cable.

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Cleaning, disinfection and sterilization of fluid light cables

1. Manual (preliminary) cleaning: The fluid light cables may be cleaned using a cleaning agent/disinfectant which is specially approved for endoscopes. Remove heavy soiling with a moist sponge or lint-free disposable cloth.*

2. Disinfection: immerse in a disinfection solution (see Appendix, page 31, 32).*

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3. Machine decontamination

Caution: Carry out chemical disinfection by machine at max. 65°C.

Caution: When preparing light cables in a cleaning and disinfection machine, particular care is necessary in the choice of chemical agents and machine programs.

4. Manual final cleaning: Clean the fiber surfaces and optical end faces with a soft cloth, sponge or cotton swab, soaked with 70% alcohol. Stubborn contaminants and deposits may be removed using cleaning paste (Art. no. 27661).

5. Sterilization: Gas sterilization or STERIS SYSTEM1® processor.

Caution: Fluid light cables must not be steam or plasma sterilized under any circumstances.

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More accessories are listed in our TELEPRESENCE catalogue.

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Type

- 1 = Cleaner
2 = Disinfectant
3 = Suitable for ultrasound baths
A = Not suitable for flexible endoscopes (according to manufacturers' instructions)
(A)= Only for the flexible endoscopes of the new generation (year of manufacture > 2005)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets

Manual preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Ecolab GmbH & Co. OHG	Sekusept Plus Sekusept Pulver Classic Ultrasdesmit AF Ultrasseptin Aktiv Ultrasseptin Classic	1, 2, 3 1, 2, 3 1, 2, 3 1, 2, 3 1, 2, 3
Estee Pharma GmbH	Afid Afid plus Sporcid FF Polsept Holifa	1, 2 2 1, 2 1, 2
Fresenius AG	Letahdeido	2
Holifa Polska Sp. zo.o	Darodor 4000 Liquido Darodor 9000	1 2
Indeba Indústria E Comercio LTDA	Darodor Sinaldehyd 2000	2
Jose Collado S.A.	Aldasan 2000 Almyrol Desoform Lysoformin 3000	1, 2 1, 2 1, 2, 3 1, 2
Lysoform Dr. Rosemann GmbH	MediDis	2
Medichem International	MediZyme Edisonite Super Mucadont-IS Mucadont-Zymaktiv Mucocit-1 neu	1 1, 3 2, 3 1, 3 1, 2, 3
METASYS Medizintechnik GmbH (formerly LONZA)	ID 50	2
Orochemie	A 20	1, 3, A
Promagent AB	Wavacide	2
	Descoton Extra Descoton forte Perfektan Endo Perfektan Neu	1, 2 2 2, 3 1, 2
Schumacher, Dr. GmbH	Perfektan TB Gigasept Gigasept AF forte Gigasept FF	1, 2, 3, A 2, 3 1, 2 1, 2
Schülke & Mayr GmbH	Gigasept Instru AF Gigasept Med forte Gigazyme Lyselvol V	1, 2, 3, A 1, 2 1, 3 2
Steris	EnzyCare 2	1, 3
	neodisher LM 2 neodisher mediclean neodisher mediclean forte neodisher medizym	1, 3 1 1, 3 1, 3
Weigert, Dr. GmbH & Co.	neodisher SeptoClean neodisher SeptoPreClean	1, A 1, 2
Whiteley Industries PTY.Ltd	Aldal Plus Matrix	1 1

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For **flexible endoscopes**, only immerse the examination sheath, **not** the housing/handle part.

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Type

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles the appropriate manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schülke & Mayr GmbH	Thermosept alka clean	1
	Thermosept DK	2
	Thermosept ED	2
	Thermosept ER	1
	Thermosept RKN-zym	1
Steris	Instru-Klenz	1, 3
	neodisher DuoClean	1, A
Weigert, Dr. GmbH & Co.	neodisher FA	1
	neodisher FA forte	1, 3
	neodisher Oxivario	1
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher mediklar	5
	neodisher medizym	1
	neodisher SeptoClean	1, 2, A
	neodisher Septo DA	2
	neodisher Septo DN	2

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The XENON 300 model 201331 20-1 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC. These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The XENON 300 model 201331 20-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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WARNING: The XENON 300 model 20 1331 20-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the XENON 300 model 20 1331 20-1. The cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200

Cables which have been shown to comply with EN/IEC 60601-1-2:

Cable type	Shielded	Length [m]	Ferrite	Used for
AC Power cord	Yes	3	No	Power supply

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Table 201		
Guidance and manufacturer's declaration – electromagnetic emissions		
XENON 300 model 20 1331 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the XENON 300 model 20 1331 20-1 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	XENON 300 model 20 1331 20-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment. The XENON 300 model 20 1331 20-1 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202			
Guidance and manufacturer's declaration – electromagnetic immunity			
XENON 300 model 20 1331 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the XENON 300 model 20 1331 20-1 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies 40% U_T (60% dip in U_T) for 5 cycles Complies 70% U_T (30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 204 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
XENON 300 model 20 1331 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the XENON 300 model 20 1331 20-1 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 V _{rms} 150 kHz to < 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V _{rms} 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the XENON 300 model 20 1331 20-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/3]\sqrt{P}$ 150 kHz to < 80 MHz $d = [3.5/3]\sqrt{P}$ 80 MHz to < 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
XENON 300 model 201331 20-1**

The XENON 300 model 201331 20-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to < 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

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To be filled out by instrument owner:
Return address/company stamp:

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KARL STORZ Endoskope
Regional Center for Endoscopy S.A.L.
St. Charles City Center, 5th Floor
Phoenicia Street, Mina Elhoshn
2020 0908 Beirut, Lebanon
Phone: +961 1 368181, Fax +961 1 365151

KARL STORZ Endoscopy South Africa, (Pty) Ltd.
8th Floor Convention Tower
Cnr Heerengracht & Coen Steytler, Foreshore
Cape Town 8001, South Africa
P.O.Box 6061, Roggebaai Cape Town,
8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSCOPY Kazakhstan
Khodjanova 17
050060 Almaty, Kazakhstan
Phone/Fax: +7 72 72 49 43 63, +7 72 72 49 41 00
e-mail: KarlStorzKazakhstan@itte.kz

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 – 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty. Ltd.
174 Parramatta Rd.
Camperdown, N.S.W. 2050, Australia
Phone: +61 02 8594 9150, Fax: +61 02 8594 9199
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Tel. No. +65 63761066, Fax. No. +65 63761068
Email : infoasia@karlstorz.com.sg
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
10-18 e-Centre @ Redhill
Singapore 159471
Tel. No. +65 65325548, Fax No. +65 65323832
Email : infosp@karlstorz.com.sg
Email : servicesgp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6Fl, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com

Other ManualsLib Projects

 www.manualslib.com

 www.manualslib.de

 www.manualslib.es

 www.manualslib.fr

 www.manualslib.nl

 www.manualslib.mx

 www.manualslib.tech 30+ Languages