

Masimo Product Catalog

2015 - Q1



This catalog contains Masimo proprietary material and is provided for information purposes only. No part of this catalog may be reproduced or retransmitted without Masimo's prior written consent.

The purpose of the product catalog is to provide a general listing of Masimo products, and it is not intended to represent products that are offered for sale. Please see local Masimo representatives for the commercially available products in specific regions/countries with the appropriate regulatory approvals/clearances.

Caution: Federal law restricts this device to sale by or on the order of a physician.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Masimo Corporation
40 Parker
Irvine, CA 92618, USA
Tel.: 949-297-7000
Fax.: 949-297-7001
www.masimo.com




EU authorized representative for Masimo Corporation:

EC REP
MDSS GmbH
Schiffgraben 41
D-30175 Hannover, Germany

© 2015 by Masimo Corporation.

Patents: www.masimo.com/patents.htm

Masimo, , 3D Alarm, Blue, Cleanshield, DCI, E1, FastSat, iSpO₂, LNCS, LNCP, MX-1, Pronto, Pronto-7, Pulse CO-Oximeter, PVI, Rad-5, Rad-8, Rad-57, Rad-87, Radical, Radical-7, rainbow, rainbow ReSposable, RRA, SatShare, SedLine, SET, SpCO, SpHb, SpMet, TF-I, and uSpO₂ are registered trademarks of Masimo Corporation.

3D Desat Index Alarm, Iris, M-LNCS, MOC-9, Phasein, Rad-5v, Rad-57, rainbow Acoustic Monitoring, RDS, Root, ROS, RRp, S-DOS, and S-ROS are trademarks of Masimo Corporation.

Apple, iPad, iPhone, iPod and iPod touch are registered trademarks of Apple Inc. registered in the U.S. and other countries.

Lightning is a trademark of Apple Inc.

"Made for iPod," "Made for iPhone," and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, or iPad may affect wireless performance.

The use of the trademarks PATIENT SAFETYNET and PSN is under license from University HealthSystem Consortium.

All other trademarks and registered trademarks are property of their respective owners.

The devices provided in the catalog may not have been licensed in accordance with applicable Federal Laws.

CONFIDENTIAL AND PROPRIETARY

Contents

Monitoring Devices	5
Masimo SET® Pulse Oximetry Devices.....	5
Masimo rainbow® Pulse CO-Oximetry Devices.....	8
Masimo rainbow SET Parameter Options.....	11
Specialty Monitoring and Connectivity Platform	13
SedLine® Brain Function Monitoring	14
Capnography	15
Masimo Patient SafetyNet™ System.....	17
Philips® Compatible Devices	23
GE Healthcare Compatible Devices	24
iSpO2®.....	25
Sensors and Cables	29
rainbow Acoustic Monitoring™ Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRA)	29
rainbow Adhesive Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet)	30
rainbow Adhesive Sensors (for rainbow SET Devices with SpO2, SpCO, and SpMet)	32
rainbow ReSposable™ Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet)	33
rainbow Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO2, SpMet; SpHb required)	35
Pronto-7 Reusable Spot Check Sensors	38
rainbow Reusable Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet).....	39
rainbow Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet).....	40
rainbow Patient Cables (for rainbow SET Devices)	41
Masimo SET ReSposable S-DOS™ Series.....	44
Masimo SET ReSposable S-ROS™ 3U Series	46
Masimo SET ReSposable S-ROS 12U Series.....	47
ReSposable MAC S-ROS 12U Series	49
LNCS Adhesive Sensors.....	51
LNCS Adhesive Sensors with Replaceable Tapes	55
LNCS Specialty Sensors	57
LNCS Sterile Sensors	58
LNCS Reusable Sensors	60
LNCS Red Patient Cables (for use with Masimo rainbow SET Devices)	62
LNC Patient Cables (for use with Masimo SET Devices).....	63
LNCS Adapter Cables.....	67
M-LNCS™ Adhesive Sensors	69
M-LNCS Adhesive Sensors with Replaceable Tapes	72
M-LNCS Specialty Sensors.....	74
M-LNCS Cyanotic Sensors.....	75
M-LNCS Reusable Sensors	76
M-LNC Patient Cables.....	78
M-LNCS Adapter Cables.....	81

Masimo Product Catalog 2015 - Q1

LNOP Adhesive Sensors..... 84

LNOP Specialty Sensors..... 85

LNOP Cyanotic Sensor 85

LNOP Reusable Sensors 86

LNOP Red Patient Cables (for use with Masimo rainbow SET Devices)..... 88

LNOP Patient Cables (for use with Masimo SET Devices) 90

LNOP Adapter Cables..... 94

Direct Connect Sensors (for use with Masimo SET Devices)..... 95

Red Direct Connect Sensors (for Masimo rainbow SET Devices)..... 96

Mountable Patient Cable Extensions and Adapter..... 97

Accessories 99

 Training Kits..... 99

 Bedside Device Accessories 101

 Handheld Device Accessories..... 107

 Sensor Accessories 111

 Pronto-7 Spot Check Pulse CO-Oximeter Accessories..... 115




 SatShare Interface Cables 117




 SatShare Extension Cable 120



Index..... 121

Monitoring Devices




Masimo SET® Pulse Oximetry Devices



PN	Description	UM	Product Image
9190	<p>Rad-8® , Horizontal</p> <p>Comes with either an LNOP® or LNCS® Patient Cable and a power cord.</p>	EA	
9191	<p>Rad-8, Vertical</p> <p>Comes with either an LNOP or LNCS Patient Cable and a power cord.</p>	EA	
9196	<p>Rad-5® Pulse Oximeter</p> <p>Comes with a Red DCI®-DC3 Sensor and choice of patient cable, a choice of available Colors of Handheld Boots and a choice of user manual.</p> <p>Shown in optional protective boot and with DCSC Sensor.</p>	EA	




PN	Description	UM	Product Image
9198	<p>Rad-5 Pulse Oximeter with Sample Pack</p> <p>Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of user manual.</p> <p>Shown in optional protective boot and with DCSC Sensor.</p>	EA	
9214	<p>Rad-5 Pulse Oximeter (International Only)</p> <p>Comes with an LNOP DCSC Sensor, a choice of available Colors of Handheld Boots and a choice of user manual.</p> <p>Shown in optional protective boot and with DCSC Sensor.</p>	EA	
9197	<p>Rad-5v™ Pulse Oximeter</p> <p>Comes with a Red DCI-DC3 Sensor and choice of patient cable, a choice of available colors of Handheld Boots and a choice of user manual.</p> <p>Shown in optional protective boot and with DCSC Sensor.</p>	EA	

PN	Description	UM	Product Image
9199	<p>Rad-5v Pulse Oximeter with Sample Pack</p> <p>Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of user manual.</p> <p>Shown in optional protective boot and with DCSC Sensor.</p>	EA	
9215	<p>Rad-5v Pulse Oximeter (International Only)</p> <p>Comes with an LNOP DCSC Sensor, a choice of available Colors of Handheld Boots and a choice of user manual.</p> <p>Shown in optional protective boot and with DCSC Sensor.</p>	EA	


Masimo rainbow® Pulse CO-Oximetry Devices

PN	Description	UM	Product Image
9500	<p>Radical-7® Handheld</p> <p>Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNOP, Red LNCS or rainbow Patient Cable. Can be configured with optional rainbow parameters and features:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb®) • Oxygen Content (SpOC) • Methemoglobin (SpMet®) • Carboxyhemoglobin (SpCO®) • Pleth Variability Index (PVI®) • Oxygen Reserve Index (ORI®) <p>This unit ships standard with a 802.11 wireless radio.</p> <p>Can be loaded with Respiration Rate from Pleth (RRp™) and Adaptive Threshold Alarms (ATA™) in cleared countries only. Check with customer service on availability.</p> <p>Optional rainbow parameters and features are sold separately and can be added to this unit upon request. The Radical-7 is compatible with all Radical® Docking Stations.</p> <p><i>Available in countries where cleared.</i></p> <p><i>ORI® parameter requires Root Docking Station</i></p>	EA	
9021	<p>Radical® Docking Station, RDS™-1</p> <p>Comes with a choice of power cord and serial output connector, Analog output / nurse call connector and SatShare® cable connector.</p>		
9023	<p>Radical® Docking Station, RDS-3</p> <p>limited communication capabilities</p> <p>Comes with choice of power cord and serial output connector and Analog output / nurse call connector.</p>		


PN	Description	UM	Product Image
9226	<p>Rad-87[®], Horizontal</p> <p>Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters. Includes either Red LNCS, Red LNOP or rainbow RC Patient Cable.</p> <p>Can be configured with optional rainbow parameters and features:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Oxygen Content (SpOC) • Methemoglobin (SpMet) • Carboxyhemoglobin (SpCO) • Pleth Variability Index (PVI) • Acoustic Respiration (RRa™) • Adaptive Threshold Alarm (ATA): <i>Available in countries where cleared.</i> <p>Comes with serial output connector and Analog output / Nurse Call connector.</p> <p>rainbow Acoustic Monitoring capable.</p> <p>Optional rainbow parameters and features are sold separately and can be added to this unit upon request.</p> <p><i>Available in countries where cleared.</i></p>	EA	
9228	<p>Rad-87, Vertical</p> <p>Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters. Includes either a Red LNCS, Red LNOP or rainbow RC Patient Cable. Can be configured with optional rainbow parameters and features:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb) • Oxygen Content (SpOC) • Methemoglobin (SpMet) • Carboxyhemoglobin (SpCO) • Pleth Variability Index (PVI) • Acoustic Respiration (RRa) • Adaptive Threshold Alarm (ATA): <i>Available in countries where cleared.</i> <p>Comes with serial output connector and Analog output / Nurse Call connector.</p> <p>rainbow Acoustic Monitoring capable.</p> <p>Optional rainbow parameters and features are sold separately and can be added to this unit upon request.</p> <p><i>Available in countries where cleared.</i></p>	EA	

PN	Description	UM	Product Image
9221	<p>Pronto-7® Spot Check Pulse CO-Oximeter</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (PI) and SpHb parameters.</p> <p>rainbow 4D DC Sensor not included. rainbow 4D DC Reusable Sensor must be ordered separately.</p>	EA	
9167	<p>Pronto® Spot Check Pulse CO-Oximeter®</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (PI) and SpHb parameters.</p> <p>rainbow Sensor not included.</p> <p>rainbow Sensor must be ordered separately.</p>	EA	
9216-U 9216	<p>Rad-57™ Configurable Pulse CO-Oximeter</p> <p>Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters.</p> <p>The Rad-57 can ship with or be field upgraded with:</p> <ul style="list-style-type: none"> • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Total Hemoglobin (SpHb) • Oxygen Content (SpOC) • PVI <p>When shipped with rainbow parameters, rainbow Sensor must be purchased separately.</p> <p>rainbow Sensors not included.</p> <p>9216-U is for SpO₂ only Devices.</p> <p>Refer to Handheld Device accessories for available colors.</p> <p>Shown with rainbow Adult Reusable Sensor and in optional protective boot.</p>	EA	


Masimo rainbow SET Parameter Options

PN	Description	UM	Product Image
9081	<p>Field Upgrader</p> <p>Used to load optional rainbow parameters and features onto compatible monitors. Current available options include:</p> <ul style="list-style-type: none"> • Total Hemoglobin (SpHb[®]) • Oxygen Content (SpOC) • Methemoglobin (SpMet) • Carboxyhemoglobin (SpCO[®]) • Pleth Variability Index (PVI[®]) • Oxygen Reserve Index (ORI[®]) • 3D Desat Index Alarm • 3D PI Delta Alarm <p>Each optional parameter and feature sold separately. Maximum quantity of 25 of any individual parameter or feature per Field Upgrader.</p> <p>Oxygen Reserve Index (ORI[®]) Not for sale in the USA. For Export Only.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 Field Upgrader • 1 to 5 operator's manuals • Parameter upgrade indicator sticker(s) • Postage-paid envelope to return Field Upgrader 	EA	
2630	<p>Total Hemoglobin (SpHb) parameter</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field</p>	EA	No image
2296	<p>Carboxyhemoglobin (SpCO) parameter</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.</p>	EA	No image
2297	<p>Methemoglobin (SpMet) parameter</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.</p>	EA	No image
2313	<p>Pleth Variability Index (PVI)</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.</p>	EA	No image
3950	<p>Oxygen Reserve Index (ORI™)</p> <p>Can be ordered on a Radical-7 when used with Root. Not for sale in the USA. For Export Only.</p>	EA	No image
3951	<p>Oxygen Reserve Index (ORI™), Fixed Term License</p> <p>Can be ordered on a Radical-7 when used with Root. Not for sale in the USA. For Export Only.</p>	EA	No image
2940	<p>Acoustic Respiration (RRa™)</p> <p>Can be ordered on a Rad-87 (PN 9226 and 9228) or on a Field Upgrader (PN 2926) to upgrade a compatible monitor in the field.</p>	EA	No image





PN	Description	UM	Product Image
3518	Respiration Rate from Pleth (RRp) Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field. In countries where RRp is cleared.	EA	No image
3519	Respiration Rate from Pleth (RRp), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field. <i>Available in countries where cleared.</i>	EA	No image
2631	Carboxyhemoglobin (SpCO), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field	EA	No image
2632	Methemoglobin (SpMet), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.	EA	No image
2633	Total Hemoglobin (SpHb), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field	EA	No image
2636	PI Delta 3D Alarm, Fixed Term License Can be ordered on a Radical-7 on a Field Upgrader (PN 9081) to upgrade a compatible monitor in the field.	EA	No image
2637	Desat Index 3D Alarm, Fixed Term License Can be ordered on a Radical-7 on a Field Upgrader (PN 9081) to upgrade a compatible monitor in the field.	EA	No image
3402	Adaptive Threshold Alarm (ATA), Fixed Term License Can be ordered on a Rad-87 (PNs 9226 and 9228). <i>Available in countries where cleared.</i>	EA	No image
3325	rainbow Acoustic Monitoring (RRa), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.	EA	No image
2341	Service Tool, MX-1[®] Upgrades up to 100 instruments.	EA	No image
2897	Service Tool, MX-3 Upgrades up to 100 instruments.	EA	No image
3351	Service Tool, MX-3 with ATA Upgrades up to 100 instruments. <i>Available in countries where cleared.</i>	EA	No image


PN	Description	UM	Product Image
2340	<p>Demonstration Tool</p> <p>Used to exercise monitors alarm default setting and demonstrate optional parameters (SpHb, SpCO, SpMet and PVI) when optional parameters are installed on a compatible monitor.</p>	EA	 <p>A white, rectangular device with a red top handle. The front panel features the 'rainbow' logo in multi-colored letters, the text 'Demonstration Tool' in bold, and the Masimo logo at the bottom right.</p>

Specialty Monitoring and Connectivity Platform

PN	Description	UM	Product Image
9515	<p>Root™</p> <ul style="list-style-type: none"> • Screen Size: 10.1 in (25.65 cm) • Resolution: 1280 x 800 Pixels • Color: 24-bit RGB • Display Type: Backlit Active Matrix TFT LCD • Touchscreen: Multi-Touch P-Cap <p>Provides Ethernet connection, nurse call connector, USB 2.0 Ports, MOC-9™ ports, and Iris® Connectivity Ports. Comes with choice of power cords.</p> <p><i>Available in countries where cleared.</i></p>	EA	 <p>A sleek, modern medical monitor with a large, dark screen and a silver/grey bezel. It has a curved design and is shown from a three-quarter perspective.</p>




SedLine® Brain Function Monitoring

PN	Description	UM	Product Image
9513	<p>SedLine Kit For use with Root and SedLine Sensors. Includes: 1 ISA CO2 MOC-9 Module 1 ISA CO2 Bracket Available in countries where cleared.</p>	EA	
9655	<p>SedLine Demo MOC-9 Module For use with Root. Does not include SedLine Patient Cable. Not for clinical use - for product demonstration only.</p>	EA	
23967	<p>Replacement SedLine MOC-9 Module For use with Root and SedLine Sensors. Available in countries where cleared.</p>	EA	
3637	<p>Replacement Sedline Patient Cable For use with Sedline MOC-9 and and SedLine Sensors. Available in countries where cleared.</p>	EA	







PN	Description	UM	Product Image
2479	<p>SedLine Sensor Single patient use 25 / box, non-sterile Not made with natural rubber latex. <i>Available in countries where cleared.</i></p>	Box	

Capnography







ISA CO2





PN	Description	UM	Product Image
3727	<p>ISA CO₂ Kit Includes: <ul style="list-style-type: none"> • 1 ISA CO₂ MOC-9 Module • 1 ISA CO₂ Bracket <i>Available in countries where cleared.</i></p>	EA	
9685	<p>ISA OR+ Kit <ul style="list-style-type: none"> • 1 ISA OR+ MOC-9 Module • 1 ISA OR+ Bracket </p>	EA	No Image
3738	<p>Nomoline Adapter 25 / Box, Non-sterile <i>Available in countries where cleared.</i></p>	Box	
3739	<p>Nomoline Airway Adapter Set 20 / Box, Non-sterile <i>For use with intubated patients</i> <i>Available in countries where cleared.</i></p>	Box	

EMMA

PN	Description	UM	Product Image
9632	<p>EMMA Capnometer Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 EMMA Capnometer mmHg • 1 EMMA Capnometer Pouch • 1 EMMA Capnometer Lanyard 	EA	
9633	<p>EMMA Capnometer Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 EMMA Capnometer kPa • 1 EMMA Capnometer Pouch • 1 EMMA Capnometer Lanyard 	EA	
3639	<p>EMMA Capnograph Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 EMMA Capnograph mmHg • 1 EMMA Pouch • 1 EMMA Lanyard 	EA	
3678	<p>EMMA Capnograph Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 EMMA Capnograph kPa • 1 EMMA Pouch • 1 EMMA Lanyard 	EA	
17448	<p>EMMA Adult / Pediatric Airway Adapter</p> <p>25 / Box, Non-sterile</p>	Box	
17449	<p>EMMA Infant Airway Adapter</p> <p>10 / Box, Non-sterile</p>	Box	

Masimo Patient SafetyNet™ System





PN	Description	UM	Product Image
9095	SafetyNet™ System Configurable for customer needs.	System	
23468	Patient SafetyNet View Station with Software (V 3.1.0.2)	EA	
24479	Patient SafetyNet View Station with Software (V 4.X.X.X)	EA	
24252	Kit SafetyNet Appliance 1UM	EA	
24254	Kit SafetyNet Appliance 2UM	EA	
24487	Kit SafetyNet Appliance 1UM and View	EA	

PN	Description	UM	Product Image
24489	Kit SafetyNet Appliance 2UM and View	EA	
23736	Halo Appliance	EA	
3623	Upgrade to V 4.X.X.X	EA	No Image
3479	Upgrade to Patient Trend Analysis Report – First Appliance	EA	No Image
3480	Upgrade to Patient Trend Analysis Report – Additional Appliance	EA	No Image
2749	CD Patient SafetyNet View Station Software (V 3.1.0.2)	EA	No Image
2572	CD Patient SafetyNet View Station Software (V 4.X.X.X)	EA	No Image
2445	Patient SafetyNet Bed License	EA	No Image
2125	Halo Bed License	EA	No Image
2447	Notification License <i>Note: For non-proprietary notification interface</i>	EA	No Image
2448	Patient SafetyNet EMR Connectivity Site License <i>Note: One required per physical address.</i>	EA	No Image
2146	SafetyNet Compatible Radical Docking Station, RDS-1 Comes with power cord and communication capabilities. <i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA	
2147	SafetyNet Compatible Radical Docking Station, RDS-3 Comes with power cord and communication capabilities. <i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA	




PN	Description	UM	Product Image
Coming soon	<p>SafetyNet Compatible Radical-7 Touchscreen Docking Station, RDS-1</p> <p>Comes with power cord and communication capabilities.</p> <p><i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i></p>	EA	
Coming soon	<p>SafetyNet Compatible Radical-7 Touchscreen Docking Station, RDS-3</p> <p>Comes with power cord and communication capabilities.</p> <p><i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i></p>	EA	
15007	<p>Colubris Workgroup Bridge</p> <p>One required per pulse oximeter; alternative for 3Com Workgroup Bridge</p>	EA	
15275	<p>NET232-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter</p>	EA	
15274	<p>NET232-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter</p>	EA	No image
16190	<p>NET232-POE-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter</p>	EA	No image
16191	<p>NET232-POE-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter</p>	EA	No image
15719	<p>NET232-U-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter</p>	EA	No image
15720	<p>NET232-U-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter</p>	EA	No image
23416	<p>Kit HP ProCurve Serial (Male-Male)</p>	EA	No image
23415	<p>Kit HP ProCurve Serial (Male-Female)</p>	EA	No image
17100	<p>Cables to GO Patch Cable 3 ft. Blue</p>	EA	No image
15768	<p>Cables to GO Patch Cable 7ft. Blue</p>	EA	No image
15769	<p>Cables to GO Patch Cable 14 ft. Blue</p>	EA	No image


PN	Description	UM	Product Image
15770	Cables to GO Patch Cable 25 ft. Blue	EA	No image
17644	HP V1910-24G-POE Port Switch	EA	
14530	Null Modem Cable, DB-9 Male to DB-9 Male, 1" Length	EA	No Image
14857	XRackPro2 Rackmount 4U Cabinet, Noise Reduction Enclosure	EA	
16846	HP Rackmount UPS	EA	
15010	HP Rackmount UPS – Universal	EA	
14424	Paging Transmitter UHF Frequency range (450 to 470 MHz). One required per floor in hospital.	EA	
16365	Paging Transmitter VHF Frequency range (135 to 175 MHz). One required per floor in hospital.	EA	
16262	Pager 440- 444 MHz, Synthesized, POCSAG Pager tunable in the 440 to 444 MHz range.	EA	

PN	Description	UM	Product Image
16263	Pager 444 - 448 MHz, Synthesized, POCSAG Pager tunable in the 444 to 448 MHz range.	EA	
16158	Pager 448 - 452 MHz, Synthesized, POCSAG Pager tunable in the 448 to 452 MHz range.	EA	
16264	Pager 452 - 456 MHz, Synthesized, POCSAG Pager tunable in the 452 to 456 MHz range.	EA	
14253	Pager 456 - 460 MHz, Synthesized, POCSAG Pager tunable in the 456 to 460 MHz range.	EA	
14254	Pager 460 - 464 MHz, Synthesized, POCSAG Pager tunable in the 460 to 464 MHz range.	EA	
14255	Pager 464 - 468 MHz, Synthesized, POCSAG Pager tunable in the 464 to 468 MHz range.	EA	
16527	Pager 135 – 142.999 MHz, Synthesized, POCSAG VHF Band A Pager tunable in the 135 to 142.999 MHz range.	EA	

PN	Description	UM	Product Image
16528	<p>Pager 143 – 150.999 MHz, Synthesized, POCSAG VHF Band B</p> <p>Pager tunable in the 143 to 150.999 MHz range.</p>	EA	
16366	<p>Pager 151 – 158.999 MHz, Synthesized, POCSAG VHF Band C</p> <p>Pager tunable in the 151 to 158.999 MHz range.</p>	EA	
16529	<p>Pager 159 – 166.999 MHz, Synthesized, POCSAG VHF Band D</p> <p>Pager tunable in the 159 to 166.999 MHz range.</p>	EA	
16530	<p>Pager 167 – 174.999 MHz, Synthesized, POCSAG VHF Band E</p> <p>Pager tunable in the 167 to 174.999 MHz range</p>	EA	
2456	<p>Patient SafetyNet Installation Services</p>	EA	No image
2457	<p>Masimo Wireless Site Survey</p>	EA	No image



Philips® Compatible Devices

PN	Description	UM	Product Image
2267	<p>Masimo SET IntelliVue® Module <i>Patient Cable must be ordered separately.</i></p>	EA	 <p>The image shows a vertical, light-colored module with a blue base. It features a circular port on the front and a small display or indicator at the top. The text 'SpO2' is visible on the top and middle sections, and the 'Masimo SET' logo is on the right side.</p>
1881	<p>Masimo SET® VueLink® Module Comes with VueLink cable.</p>	EA	 <p>The image shows a vertical, light-colored module with a blue base. It features a circular port on the front and a small display or indicator. The text 'VueLink AUX.PLUS' is visible at the top, and 'VueLink' is on the middle section.</p>
1841	<p>VueLink® Cable For use with Masimo SET VueLink Module, 3 ft.</p>	EA	 <p>The image shows a coiled grey cable with a blue connector at one end and a metal connector at the other.</p>




PN	Description	UM	Product Image
18710	<p>USB-to-RS232 Converter Cable</p> <p>For connecting Root with IntelliBridge(R) (EC5 portion) and VueLink(R) modules.</p>	EA	



Philips®, IntelliVue® and VueLink® are registered trademarks of Philips Corporation.



GE Healthcare Compatible Devices



PN	Description	UM	Product Image
9477	<p>uSpO₂® Pulse Oximetry Cable Kit</p> <p>For use with GE ApexPro FH Telemetry Transmitters.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 LNCS Adtx-3 • 1 LNCS Pdtx-3 • 1 LNCS Inf-3 • 1 LNCS Neo-3 sensor <p><i>US Only</i></p> <p>May only be sold to End-Customers with pre-existing installed ApexPro Telemetry Systems. For more information, see GE Upgrade Profiles under OEM on RadNet.</p>	EA	
9476	<p>uSpO₂ Pulse Oximetry Cable Kit</p> <p>For use with GE ApexPro, ApexPro CH, and CARESCAPE T4 & T14 Telemetry Transmitters.</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 LNCS Adtx-3 • 1 LNCS Pdtx-3 • 1 LNCS Inf-3 • 1 LNCS Neo-3 sensor <p><i>Available in countries where cleared.</i></p> <p>May only be sold to End-Customers with pre-existing installed ApexPro Telemetry Systems. For more information, see GE Upgrade Profiles under OEM on RadNet.</p>	EA	

iSpO2®

PN	Description	UM	Product Image
3588	<p>iSpO2 Rx 30-Pin Series, Adult Compatible with:</p> <ul style="list-style-type: none"> • IOS 7 and 8 • iPhone 4 • iPhone 4S • iPad (2nd Generation) • iPad (3rd Generation) <p><i>Available in countries where cleared. Phone not included.</i></p>	EA	
3590	<p>iSpO2 Rx 30-Pin Series, Pediatric Compatible with:</p> <ul style="list-style-type: none"> • IOS 7 and 8 • iPhone 4 • iPhone 4S • iPad (2nd Generation) • iPad (3rd Generation) <p><i>Available in countries where cleared. Phone not included.</i></p>	EA	
3721	<p>iSpO2 Rx 30-Pin Series, M-LNCST™ Compatible with:</p> <ul style="list-style-type: none"> • IOS 7 and 8 • iPhone 4 • iPhone 4S • iPad (2nd Generation) • iPad (3rd Generation) <p><i>Available in countries where cleared. Phone not included. M-LNCST sensors sold separately.</i></p>	EA	





PN	Description	UM	Product Image
9575	<p>iSpO2 Rx Lightning Series, Adult Compatible with:</p> <ul style="list-style-type: none"> • OS 7 and 8 • iPhone 6 • iPhone 6 Plus • iPhone 5 • iPhone 5S • iPhone 5C • iPad (4th Generation) • iPad Mini (1st Generation) • iPad Mini (2nd Generation) • iPad Mini (3rd Generation) • iPod (5th Generation) • iPad Air • iPad Air 2 <p><i>Available in countries where cleared. Phone not included.</i></p>	EA	
9576	<p>iSpO2 Rx Lightning Series, Pediatric Compatible with:</p> <ul style="list-style-type: none"> • IOS 7 and 8 • iPhone 6 • iPhone 6 Plus • iPhone 5 • iPhone 5S • iPhone 5C • iPad (4th Generation) • iPad Mini (1st Generation) • iPad Mini (2nd Generation) • iPod (5th Generation) • iPad Air • iPad Air 2 <p><i>Available in countries where cleared. Phone not included.</i></p>	EA	

PN	Description	UM	Product Image
9569	<p>iSpO2 Rx Lightning Series, M-LNCS</p> <p>Compatible with:</p> <ul style="list-style-type: none"> • OS 7 and 8 • iPhone 6 • iPhone 6 Plus • iPhone 5 • iPhone 5S • iPhone 5C • iPad (4th Generation) • iPad Mini (1st Generation) • iPad Mini (2nd Generation) • iPod (5th Generation) • iPad Air • iPad Air 2 <p><i>Available in countries where cleared. Phone not included. M-LNCS sensors sold separately.</i></p>	EA	
9578	<p>iSpO2 Rx Micro USB Series, Adult</p> <p>Compatible with:</p> <ul style="list-style-type: none"> • Android Jelly Bean**, KitKat, and Lollipop • Samsung Galaxy Note 10 (2014 edition) • Samsung Note II • Samsung Note 3 • Samsung Note 4 • Nexus 5 • Nexus 6 • Nexus 7 • Nexus 9 • Samsung Galaxy S3 • Samsung Galaxy S4 • Samsung Galaxy S5 • Huawei Ascend Mate 7 • Droid Mini • Motorola RAZR M • Samsung S3 • Samsung S4 • Samsung Note 2 <p><i>Available in countries where cleared. Phone not included.</i></p>	EA	




PN	Description	UM	Product Image
9579	<p>iSpO2 Rx Micro USB Series, Pediatric Compatible with:</p> <ul style="list-style-type: none"> • Android Jelly Bean**, KitKat, and Lollipop • Samsung Galaxy Note 10 (2014 edition) • Samsung Note II • Samsung Note 3 • Samsung Note 4 • Nexus 5 • Nexus 6 • Nexus 7 • Nexus 9 • Samsung Galaxy S3 • Samsung Galaxy S4 • Samsung Galaxy S5 • Huawei Ascend Mate 7 • Droid Mini • Motorola RAZR M • Samsung S3 • Samsung S4 • Samsung Note 2 <p><i>Available in countries where cleared. Phone not included.</i></p>	EA	
9577	<p>iSpO2 Rx Micro USB Series, M-LNCS Compatible with:</p> <ul style="list-style-type: none"> • Android Jelly Bean**, KitKat, and Lollipop • Samsung Galaxy Note 10 (2014 edition) • Samsung Note II • Samsung Note 3 • Samsung Note 4 • Nexus 5 • Nexus 6 • Nexus 7 • Nexus 9 • Samsung Galaxy S3 • Samsung Galaxy S4 • Samsung Galaxy S5 • Huawei Ascend Mate 7 • Droid Mini • Motorola RAZR M • Samsung S3 • Samsung S4 • Samsung Note 2 <p><i>Available in countries where cleared. Phone not included. M-LNCS sensors sold separately.</i></p>	EA	


Sensors and Cables

rainbow Acoustic Monitoring™ Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRA)


PN	Description	UM	Product Image
3475	<p>RAS-125c, Acoustic Respiration Cloth Sensor Adult / Pediatric Adhesive Sensors (RRa) Single patient use 10 / box, non-sterile Weight > 10kg Not made with natural rubber latex. <i>Available in countries where cleared.</i></p>	Box	 <p>The image shows the RAS-125c sensor, which consists of a white adhesive cloth sensor connected to a white cable with a red connector. The sensor has a small white label with 'RAS-125c' and 'RRa' printed on it.</p>
3483	<p>RAS-125c, Short Term Monitoring Acoustic Respiration Cloth Sensor Adult / Pediatric Adhesive Sensors (RRa) Single patient use 10 / box, non-sterile Weight > 10kg Not made with natural rubber latex. <i>Available in countries where cleared.</i></p>	Box	 <p>The image shows the RAS-125c sensor, similar to the one above, but with a white label that says '6 Hour' attached to the cable. The sensor has a small white label with 'RAS-125c' and 'RRa' printed on it.</p>
3660	<p>RAM Dual Cable RC-10 For use with RAM RAS and rainbow or M-LNCS sensors</p>	EA	 <p>The image shows the RAM Dual Cable RC-10, which is a white cable with two red connectors at one end and a white connector at the other end.</p>
3661	<p>RAM Dual Cable LNC-10 For use with RAM RAS and SpO² LNCS sensors</p>	EA	 <p>The image shows the RAM Dual Cable LNC-10, which is a white cable with two red connectors at one end and a white connector at the other end.</p>

rainbow Adhesive Sensors (for rainbow SET Devices with SpHb, SpO₂, and SpMet)

PN	Description	UM	Product Image
3792	<p>rainbow R1 25 Adult Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight >30 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i></p>	Box	
3793	<p>rainbow R1 20 Pediatric / Slender Digit Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i></p>	Box	
2414	<p>cablerainbow R1 25L Adult / Neonatal Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight <3 kg* or >30 kg Not made with natural rubber latex. Use replacement tape 2623. <i>Requires rainbow RC patient cable</i> *SpO₂ And SpMet only</p>	Box	




PN	Description	UM	Product Image
2415	<p>rainbow R1 20L Infant Adhesive Sensors (SpHb, SpO₂, SpMet) Single patient use 10 / box, non-sterile Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement tape 2624. <i>Requires rainbow RC patient cable</i></p>	Box	 <p>The image shows a long, thin, white adhesive sensor strip. On the left end, there is a white plastic connector with a red Masimo logo. The word "rainbow" is printed in small black letters along the length of the strip. On the right end, there is a brown plastic connector with a red Masimo logo and a small black square sensor element.</p>




rainbow Adhesive Sensors (for rainbow SET Devices with SpO₂, SpCO, and SpMet)

PN	Description	UM	Product Image
2221	<p>rainbow R25 Adult Adhesive Sensors (SpO₂, SpCO, SpMet) Single patient use 10 / box, non-sterile Weight >30 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i></p>	Box	
2222	<p>rainbow R20 Pediatric Adhesive Sensors (SpO₂, SpCO, SpMet) Single patient use 10 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i></p>	Box	
2219	<p>rainbow R25-L Adult / Neonatal Adhesive Sensors (SpO₂, SpCO, SpMet) Single patient use 10 / box, non-sterile Weight <3 kg* or >30 kg Not made with natural rubber latex. Use replacement tape 2623. <i>Requires rainbow RC patient cable</i> <i>*SpO₂ And SpMet only</i></p>	Box	
2220	<p>rainbow R20-L Infant Adhesive Sensors (SpO₂, SpCO, SpMet) Single patient use 10 / box, non-sterile Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement tape 2624. <i>Requires rainbow RC patient cable</i></p>	Box	




rainbow ReSposable™ Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet)


rainbow ReSposable Sensors (for use with Masimo rainbow SET Devices with SpHb, SpO2, and SpMet parameters)



PN	Description	UM	Product Image
3457	<p>rainbow ReSposable R2-25 Sensor System</p> <p>Adult R2-25a disposable optical sensors and R2-25r reusable optical sensors (SpO₂, SpHb, SpMet)</p> <p>20 R2-25a / box (DOS)</p> <p>1 R2-25r / box (ROST™)</p> <p>Non-sterile</p> <p>Weight > 30 kg</p> <p>Not made with natural rubber latex.</p> <p><i>Requires rainbow RC patient cable</i></p>	Box	
3458	<p>rainbow® ReSposable™ R2-20 Sensor System</p> <p>Pediatric R2-20a disposable optical sensors and R2-20r reusable sensors (SpO₂, SpHb, SpMet)</p> <p>20 R2-20a / box (DOS)</p> <p>1 R2-20r / box (ROS)</p> <p>Non-sterile</p> <p>Weight 10 kg – 50 kg</p> <p>Not made with natural rubber latex.</p> <p><i>Requires rainbow RC patient cable</i></p>	Box	
3395	<p>rainbow® ReSposable™ R2-25r Sensors</p> <p>Adult R2-25r reusable optical sensors (SpO₂, SpHb, SpMet)</p> <p>5 R2-25r / box (ROS)</p> <p>Non-sterile</p> <p>Weight > 30 kg</p> <p>Not made with natural rubber latex.</p> <p><i>Requires R2-25a</i></p>	Box	

PN	Description	UM	Product Image
3411	<p>rainbow[®] ReSposable[™]R2-20r Sensors Pediatric / Slender Digit R2-20r reusable optical sensors (SpO₂, SpHb, SpMet) 5 R2-20r / box (ROS) Non-sterile Weight 10 – 50 kg Not made with natural rubber latex. <i>Requires R2-20a</i></p>	Box	
2753	<p>rainbow[®] ReSposable[™]R2-25a Sensors Adult R2-25a disposable optical sensors (SpO₂, SpHb, SpMet) 25 R2-25a / box (DOS) Non-sterile Weight > 30 kg Not made with natural rubber latex. <i>Requires R2-25r</i></p>	Box	
2755	<p>rainbow[®] ReSposable[™] R2-20a Sensors Pediatric / Slender Digit R2-20a disposable optical sensors (SpO₂, SpHb, SpMet) Single patient use 25 R2-20a / box (DOS) Non-sterile Weight 10 – 50 kg Not made with natural rubber latex. <i>Requires R2-20r</i></p>	Box	




rainbow Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO₂, SpMet; SpHb required)

PN	Description	UM	Product Image
3418	<p>rainbow DCI® SC 200 Adult reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 200 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3419	<p>rainbow DCI-P SC 200 Pediatric reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 200 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3420	<p>rainbow DCI SC 400 Adult reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 400 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	




PN	Description	UM	Product Image
3421	<p>rainbow DCI-P SC 400 Pediatric reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 400 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3647	<p>rainbow DCI SC 1000 Adult reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 1000 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3648	<p>rainbow DCI-P SC 1000 Pediatric reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 1000 Spot Checks <i>Requires rainbow RC Patient Cable</i></p>	EA	
3797	<p>rainbow DCI-mini SC 200 infant reusable sensor. (SpHb, SpO₂ and SpMet) Weight 3 – 30 kg Not made with natural rubber latex. 200 Spot Checks <i>Requires rainbow RC Patient Cable</i> <i>Not For Sale in USA For Export Only</i></p>	EA	


PN	Description	UM	Product Image
3798	<p>rainbow DCI-mini SC 400 infant reusable sensor. (SpHb, SpO₂ and SpMet) Weight 3 – 30 kg Not made with natural rubber latex. 400 Spot Checks <i>Requires rainbow RC Patient Cable</i> <i>Not For Sale in USA For Export Only</i></p>	EA	
3799	<p>rainbow DCI-mini SC 1000 infant reusable sensor. (SpHb, SpO₂ and SpMet) Weight 3 – 30 kg Not made with natural rubber latex. 1000 Spot Checks <i>Requires rainbow RC Patient Cable</i> <i>Not For Sale in USA For Export Only</i></p>	EA	

Pronto-7 Reusable Spot Check Sensors



PN	Description	UM	Product Image
9400	<p>rainbow 4D DC Sensor, Small</p> <p>Reusable Sensor</p> <p>1 / box, non-sterile</p> <p>Used with finger sizes between 8mm – 12mm diameter.</p> <p>Not made with natural rubber latex.</p> <p>Comes with 500 or 1000 Spot checks.</p> <p><i>Must add Spot checks separately (PN 3312)</i></p> <p><i>Requires a Pronto-7.</i></p>	EA	
9401	<p>rainbow 4D DC Sensor, Medium</p> <p>Reusable Sensor</p> <p>1 / box, non-sterile</p> <p>Used with finger sizes between 12mm – 16mm diameter.</p> <p>Not made with natural rubber latex.</p> <p>Comes with 500 or 1000 Spot checks.</p> <p><i>Must add Spot checks separately (PN 3312)</i></p> <p><i>Requires a Pronto-7.</i></p>	EA	
9402	<p>rainbow 4D DC Sensor, Large</p> <p>Reusable Sensor</p> <p>1 / box, non-sterile</p> <p>Used with finger sizes between 16mm – 20mm diameter.</p> <p>Not made with natural rubber latex.</p> <p>Comes with 500 or 1000 Spot checks.</p> <p><i>Must add Spot checks separately (PN 3312)</i></p> <p><i>Requires a Pronto-7.</i></p>	EA	

rainbow Reusable Sensors (for rainbow SET Devices with SpCO, SpO₂, and SpMet)





PN	Description	UM	Product Image
2696	<p>rainbow DCI Adult Reusable Sensor, 3 ft. (SpCO, SpMet, SpO₂) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i></p>	EA	
2697	<p>rainbow DCI-P Pediatric Reusable Sensor, 3 ft. (SpCO, SpMet, SpO₂) 1 / box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i></p>	EA	
3657	<p>rainbow DCI-6 Adult Reusable Sensor, 6ft (SpCO, SpMet, SpO₂) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i></p>	EA	




PN	Description	UM	Product Image
3658	<p>rainbow DCI-P-6 Pediatric Reusable Sensor, 6ft (SpCO, SpMet, SpO₂) 1 box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex <i>Requires rainbow RC Cable</i></p>	EA	


rainbow Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO₂, and SpMet)

PN	Description	UM	Product Image
2201	<p>rainbow DCI dc3 Adult Reusable Direct Connect Sensor, 3 ft. (SpO₂, SpCO, SpMet) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
2069	<p>rainbow DCI-P-dc3 Pediatric Reusable Direct Connect Sensor, 3 ft. (SpO₂, SpCO, SpMet) 1 / box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex.</p>	EA	





rainbow Patient Cables (for rainbow SET Devices)




PN	Description	UM	Product Image
2405	<p>rainbow RC-1 rainbow 20-pin Patient Cable, 1 ft For use with rainbow and M-LNCS sensors</p>	EA	
2406	<p>rainbow RC-4 rainbow 20-pin Patient Cable, 4 ft For use with rainbow and M-LNCS sensors</p>	EA	
2404	<p>rainbow RC-12 rainbow 20-pin Patient Cable, 12 ft For use with rainbow and M-LNCS sensors</p>	EA	
3659	<p>rainbow RC-EXT-8 rainbow extension cable, 8ft For use with rainbow and M-LNCS sensors</p>	EA	

PN	Description	UM	Product Image
3342	<p>rainbow RC25-12 rainbow 25-pin, Patient Cable, 12ft For use with rainbow and M-LNCS sensors</p>	EA	 <p>A 12-foot patient cable with a white 25-pin connector on one end and a red and white connector on the other. The cable is coiled.</p>
3343	<p>rainbow RC25-12RA rainbow 25-pin, right angle Patient Cable, 12ft For use with rainbow and M-LNCS sensors</p>	EA	 <p>A 12-foot patient cable with a white 25-pin connector on one end and a right-angle red and white connector on the other. The cable is coiled.</p>
3628	<p>rainbow RC25-4 rainbow 25-pin, Patient Cable, 4 ft. For use with rainbow and M-LNCS sensors</p>	EA	 <p>A 4-foot patient cable with a white 25-pin connector on one end and a red and white connector on the other. The cable is coiled.</p>



PN	Description	UM	Product Image
3375	<p>rainbow RC25-4RA rainbow 25-pin, Patient Cable, 4 ft. For use with rainbow and M-LNCS sensors</p>	EA	

Masimo SET ReSposable S-DOS™ Series


PN	Description	UM	Product Image
2434	<p>ReSposable™ S-DOS™ 25 Adult SpO₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight > 30 kg Not made with natural rubber latex. For use with S-ROS™ 3U or 12U series.</p>	Box	
2435	<p>ReSposable S-DOS 20 Pediatric SpO₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight 10 - 50 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series</p>	Box	
2469	<p>ReSposable S-DOS 25L Neonatal / Adult SpO₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight > 30 kg, < 3 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series</p>	Box	
2470	<p>ReSposable S-DOS 20L Infant SpO₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight 3-30 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series</p>	Box	



PN	Description	UM	Product Image
2471	<p>ReSposable S-DOS 30L Large Adult SpO₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight 3 - 10kg or > 40 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series</p>	Box	 <p>The image shows a long, thin, tan-colored adhesive strip. On the left end, there is a circular sensor window with a black lens. A dashed line runs horizontally across the middle of the strip. On the right end, there is a small red logo and the word 'MASIMO' in black capital letters.</p>
2476	<p>ReSposable S-DOS 25Pt-L Neonatal SpO₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight < 1 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series</p>	Box	 <p>The image shows a long, thin, white adhesive strip. On the left end, there is a circular sensor window with a black lens. A dashed red line runs horizontally across the middle of the strip. On the right end, there is a small red logo.</p>
2477	<p>ReSposable S-DOS 20L-500 Non-adhesive Neonatal SpO₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight < 1 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series</p>	Box	 <p>The image shows a long, thin, white adhesive strip. On the left end, there is a circular sensor window with a black lens. A dashed red line runs horizontally across the middle of the strip. On the right end, there is a small red logo.</p>

Masimo SET ReSposable S-ROS™ 3U Series

PN	Description	UM	Product Image
2473	<p>LNCS ReSposable S-ROS 3U SpO₂ reusable optical sensor (ROS) Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series. <i>Requires patient cable.</i></p>	EA	
2472	<p>M-LNCS ReSposable S-ROS 3U SpO₂ reusable optical sensor (ROS) Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series. <i>Requires patient cable.</i></p>	EA	



Masimo SET ReSposable S-ROS 12U Series

PN	Description	UM	Product Image
2474	<p>ReSposable Red S-ROS 12U SpO₂ reusable optical sensor (ROS), 12 ft., 20-pin Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	 <p>A black cable with a red connector at one end and a grey connector at the other. The cable is coiled and has a small circular logo on the side.</p>
2475	<p>ReSposable S-ROS 12U SpO₂ reusable optical sensor (ROS), 12 ft., 14-pin Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	 <p>A black cable with a grey connector at one end and a white connector at the other. The cable is coiled and has a small circular logo on the side.</p>
2619	<p>ReSposable SpaceLabs S-ROS 12U SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	 <p>A black cable with a grey connector at one end and a white connector at the other. The cable is coiled and has a small circular logo on the side.</p>
2620	<p>ReSposable Philips S-ROS 12U SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	 <p>A black cable with a blue connector at one end and a white connector at the other. The cable is coiled and has a small circular logo on the side.</p>



PN	Description	UM	Product Image
2621	<p>ReSposable DB9 S-ROS 12U SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	
2622	<p>ReSposable GE S-ROS 12U SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	



ReSposable MAC S-ROS 12U Series





PN	Description	UM	Product Image
2566	<p>ReSposable MAC 395 S-ROS 12U N-395 SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	
2567	<p>ReSposable MAC 180 S-ROS 12U N-180 SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	
2568	<p>ReSposable MAC GE S-ROS 12U GE SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	
2569	<p>ReSposable MAC SL S-ROS 12U SpaceLabs SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	



PN	Description	UM	Product Image
2570	<p>ReSposable MAC SL2 S-ROS 12U SpaceLabs SpO₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series</p>	EA	
2571	<p>ReSposable MAC CMS S-ROS 12U Philips SpO₂ reusable optical sensor (ROS), 12 ft., 12-pin Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.</p>	EA	

LNCS Adhesive Sensors





PN	Description	UM	Product Image
1859	<p>LNCS Amtx Adult SpO₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.</p>	Box	
2317	<p>LNCS Amtx-3 Adult SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.</p>	Box	




PN	Description	UM	Product Image
2918	<p>LNCS E1® Adult SpO₂ disposable ear sensor, 3 ft. Single patient use 10 / box, non-sterile, non-adhesive Weight > 30 kg Not made with natural rubber latex. Requires Masimo SET MSLP board, MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology</p>	Box	 <p>The image shows a white, flexible ear sensor cable. One end has a white plastic ear hook with two small sensors. The other end has a blue connector with a white circular component.</p>
3858	<p>LNCS TFA-1 Adult / Pediatric SpO₂ disposable tranreflectance forehead sensor, 3ft Single patient use 10 / box non-sterile Weight > 10 kg Not made with natural rubber latex Not for sale in USA For Export only.</p>	Box	 <p>The image shows a white, oval-shaped forehead sensor with a blue connector. The sensor has a red Masimo logo on its surface. A white cable connects the sensor to the connector.</p>

PN	Description	UM	Product Image
3857	<p>LNCS Red TFA-1 Adult / Pediatric SpO₂ disposable tranreflectance forehead sensor 3 ft Single patient use 10 / box non-sterile Weight > 10 kg Not made with natural rubber latex Requires Masimo SET MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX (Version 7.7 or higher) technology Not for sale in USA For Export only.</p>	Box	
1860	<p>LNCS Pdtx Pediatric SpO₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	Box	
2318	<p>LNCS Pdtx-3 Pediatric SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	Box	
1861	<p>LNCS Inf-L Infant SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex.</p>	Box	

PN	Description	UM	Product Image
1862	<p>LNCS Neo-L Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex.</p>	Box	
1901	<p>LNCS NeoPt-L Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex.</p>	Box	




LNCS Adhesive Sensors with Replaceable Tapes

PN	Description	UM	Product Image
2328	<p>LNCS Inf Infant SpO₂ Adhesive Sensors, 18 in. Single Patient Use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.</p>	Box	
2319	<p>LNCS Inf-3 Infant SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.</p>	Box	
2329	<p>LNCS Neo Neonatal SpO₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.</p>	Box	
2320	<p>LNCS Neo-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.</p>	Box	

PN	Description	UM	Product Image
2330	<p>LNCS NeoPt Neonatal SpO₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2321	<p>LNCS NeoPt-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2331	<p>LNCS NeoPt-500 Neonatal SpO₂ non-adhesive sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2322.</p>	Box	




LNCS Specialty Sensors




LNCS Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat® setting

PN	Description	UM	Product Image
2411	<p>LNCS Trauma Adult SpO₂ Adhesive Sensors Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2412	<p>LNCS Newborn Neonatal Neonatal SpO₂ Adhesive Sensors Single patient use 20 / box, sterilizable Weight < 3 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2413	<p>LNCS Newborn Infant / Pediatrics Infant / Pediatric SpO₂ Adhesive Sensor Single patient use 20 / box, sterilizable Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement wrap 2322.</p>	Box	




LNCS Sterile Sensors



Note: The sensors are not intended to be used for evaluations or carried as sales samples. If there is a need for samples of sterile sensors, contact your manager for approval and include the e-mail with FRM-1386.

PN	Description	UM	Product Image
2601	<p>LNCS Adtx-3 Adult SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight > 30 kg Not made with natural rubber latex.</p>	Box	
2602	<p>LNCS Pdtx-3 Pediatric SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	Box	
2603	<p>LNCS Inf-3 Infant SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.</p>	Box	





PN	Description	UM	Product Image
2604	<p>LNCS Neo-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.</p>	Box	
2605	<p>LNCS NeoPt-3 Neonatal SpO₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2606	<p>LNCS NeoPt-500 Neonatal SpO₂ non-adhesive sensor, 18 in. Single patient use 20 / box, Sterile Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2322.</p>	Box	

LNCS Reusable Sensors





PN	Description	UM	Product Image
1863	<p>LNCS DC-I Adult Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
1864	<p>LNCS DCI-P Pediatric Reusable Sensor, 3 ft. 1 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	EA	
2653	<p>LNCS DB-I Adult Soft Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
1895	<p>LNCS TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	





PN	Description	UM	Product Image
2258	<p>LNCS Y-I Multisite Reusable Sensor, 3 ft. 1 / box, non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg Not made with natural rubber latex.</p>	EA	
1896	<p>LNCS TF-I Adult Transflectance Forehead, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	





LNCS Red Patient Cables (for use with Masimo rainbow SET Devices)





PN	Description	UM	Product Image
2365	Red LNC™-01 LNCS 20-pin SpO ₂ Spot Check Patient Cable.	EA	
2055	Red LNC™-04 LNCS 20-pin SpO ₂ Patient Cable, 4 ft.	EA	
2056	Red LNC™-10 LNCS 20-pin SpO ₂ Patient Cable, 10 ft	EA	
2057	Red LNC™-14 LNCS 20-pin SpO ₂ Patient Cable, 14 ft.	EA	

LNC Patient Cables (for use with Masimo SET Devices)





PN	Description	UM	Product Image
2364	LNC-1 LNCS Series Patient Cable.	EA	 A short patient cable with a white connector on one end and a red connector on the other.
2017	LNC-4 LNCS Series Patient Cable, 4 ft.	EA	 A 4-foot patient cable with a white connector on one end and a red connector on the other.
1814	LNC-10 LNCS Series Patient Cable, 10 ft.	EA	 A 10-foot patient cable with a white connector on one end and a red connector on the other.
3630	Red 25 LNC-4 LNCS Series 25-pin SpO ₂ Patient Cable, 4 ft	EA	 A 4-foot patient cable with a white connector on one end and a red connector on the other.




PN	Description	UM	Product Image
3376	<p>Red 25 LNC-4RA Series LNCS Series 25-pin right angle SpO₂ Patient Cable, 4 ft</p>	EA	
3345	<p>Red 25 LNC-10 Series LNCS Series 25-pin SpO₂ Patient Cable, 10 ft</p>	EA	
2013	<p>LNC-14 LNCS Series Patient Cable, 14 ft.</p>	EA	
2021	<p>LNC-4-Ext LNC Extension Cable, 4 ft.</p>	EA	

PN	Description	UM	Product Image
2281	<p>LNC MP10 LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO₂ LNCS Dual Key Patient Cable, 10 ft.</p>	EA	
3755	<p>LNC MP4 LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO₂ LNCS Dual Key Patient Cable, 4 ft</p>	EA	
2430	<p>LNC-4-GE LNCS to GE Patient Cable, 4 ft.</p>	EA	
2016	<p>LNC-10-GE LNCS to GE Patient Cable, 10 ft.</p>	EA	





PN	Description	UM	Product Image
2432	<p>LNC-SL-10 LNCS to Spacelabs Patient Cable, 10 ft.</p>	EA	
2285	<p>LNC-04-DB9 LNCS to DB9 Connector Patient Cable, 4 ft.</p>	EA	
2286	<p>LNC-10-DB9 LNCS to DB9 Connector Patient Cable, 10 ft.</p>	EA	
2287	<p>LNC-14-DB9 LNCS to DB9 Connector Patient Cable, 14 ft.</p>	EA	



LNCS Adapter Cables




PN	Description	UM	Product Image
1816	<p>LNCS to PC Adapter Cable LNCS Series to LNOP PC Adapter Cable, 18 in.</p>	EA	 <p>A white cable with a red-tipped LNOP connector on one end and a white PC connector on the other.</p>
2263	<p>LNC CMS LNCS Series to Philips 12-pin Instrument Adapter Cable, 3 m.</p>	EA	 <p>A white cable with a blue Philips 12-pin connector on one end and a blue LNC CMS connector on the other.</p>
2264	<p>LNC MAC-GE LNCS Series to GE Adapter Cable, 3 m.</p>	EA	 <p>A white cable with a blue GE connector on one end and a blue LNC MAC-GE connector on the other.</p>
2266	<p>LNC MAC-SL LNCS Series to Spacelabs Adapter Cable, 3 m.</p>	EA	 <p>A white cable with a Spacelabs connector on one end and a white LNC MAC-SL connector on the other.</p>

PN	Description	UM	Product Image
2336	<p>LNC MAC-SL2 LNCS Series to Spacelabs Adapter Cable, 3 m.</p>	EA	
2268	<p>LNC MAC-180 LNCS Series to N-180 Adapter Cable, 10 ft.</p>	EA	
2269	<p>LNC MAC-395 LNCS Series to N-395 Adapter Cable, 10 ft.</p>	EA	





M-LNCS™ Adhesive Sensors




PN	Description	UM	Product Image
2508	<p>M-LNCS Adtx Adult SpO₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.</p>	Box	
2509	<p>M-LNCS Adtx-3 Adult SpO₂ Adhesive Sensor, 3 ft. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.</p>	Box	
2510	<p>M-LNCS Pdtx Pediatric SpO₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight 10 kg - 50 kg Not made with natural rubber latex.</p>	Box	
2511	<p>M-LNCS Pdtx-3 Pediatric SpO₂ Adhesive Sensor, 3 ft. Single patient use 20 / box, sterilizable Weight 10 kg - 50 kg Not made with natural rubber latex.</p>		

PN	Description	UM	Product Image
2919	<p>M-LNCS E1 Adult SpO₂ disposable ear sensor, 3 ft. Single patient use 10 / box, non-sterile, non-adhesive Weight > 30 kg Not made with natural rubber latex. Requires Masimo SET MSLP board, MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology</p>	Box	
3856	<p>M-LNCS TFA-1 Adult / Pediatric SpO₂ disposable transreflectance forehead sensor, 3 ft Single patient use 10 / box, non-sterile Weight > 10 kg Not made with natural rubber latex Not for sale in USA for Export only</p>	Box	
3855	<p>M-LNCS Red TFA-1 Adult / Pediatric SpO₂ disposable transreflectance forehead sensor, 3 ft Single patient use 10 / box, non-sterile Weight > 10 kg Not made with natural rubber latex Requires Masimo SET MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX (Version 7.7 or higher) technology Not for sale in USA for Export only</p>	Box	

PN	Description	UM	Product Image
3670	<p>M-LNCS Inf-L Infant SPO₂ adhesive sensors, 3 ft Single patient use 20 / box, sterilizable Weight 3 kg - 20 kg Not made with natural rubber latex</p>	Box	
3668	<p>M-LNCS Neo-L Neonatal SPO₂ adhesive sensors, 3 ft Single patient use 20 / box, sterilizable Weight 3 kg - 20 kg Not made with natural rubber latex</p>	Box	
3669	<p>M-LNCS NeoPt-L Neonatal SPO₂ adhesive sensors, 3 ft Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex</p>	Box	




M-LNCS Adhesive Sensors with Replaceable Tapes

PN	Description	UM	Product Image
2512	<p>M-LNCS Inf Infant SpO₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.</p>	Box	
2513	<p>M-LNCS Inf-3 Infant SpO₂ Adhesive Sensor, 3 ft. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.</p>	Box	
2514	<p>M-LNCS Neo Neonatal / Adult SpO₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.</p>	Box	
2515	<p>M-LNCS Neo-3 Neonatal / Adult SpO₂ Adhesive Sensor, 3 ft. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.</p>	Box	


PN	Description	UM	Product Image
2516	<p>M-LNCS NeoPt Neonatal SpO₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2517	<p>M-LNCS NeoPt-3 Neonatal SpO₂ Adhesive Sensor, 3 ft. Single patient use Weight < 3 kg or > 40 kg Not made with natural rubber latex. 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2518	<p>M-LNCS NeoPt-500 Neonatal SpO₂ Non-Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2322.</p>	Box	

M-LNCS Specialty Sensors





M-LNCS Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.



PN	Description	UM	Product Image
2521	<p>M-LNCS Trauma Adult SpO₂ Sensor Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2519	<p>M-LNCS Newborn Neonatal Neonatal SpO₂ Sensor Single patient use 20 / box, sterilizable Weight < 3 kg Not made with natural rubber latex. Use replacement wrap 2309.</p>	Box	
2520	<p>M-LNCS Newborn Infant / Pediatric Infant / Pediatric SpO₂ Sensor Single patient use 20 / box, sterilizable Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement wrap 2322.</p>	Box	

M-LNCS Cyanotic Sensors





PN	Description	UM	Product Image
2522	<p>M-LNCS Blue Neonatal / Infant / Pediatric SpO₂ Adhesive Sensor Single patient use 10 / box, Non-Sterile Weight 2.5 – 30kg Not made with natural rubber latex. Use replacement tape 2822.</p>	Box	

M-LNCS Reusable Sensors




PN	Description	UM	Product Image
2501	<p>M-LNCS DC-I Adult Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
2502	<p>M-LNCS DCI-P Pediatric / Slender Digit Reusable Sensor, 3 ft. 1 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	EA	
2507	<p>M-LNCS DB-I Adult Soft Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
2503	<p>M-LNCS TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	

PN	Description	UM	Product Image
2505	<p>M-LNCS Y-I Multisite Reusable Sensor, 3 ft. 1 / box, non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg Not made with natural rubber latex.</p>	EA	
2504	<p>M-LNCS TF-I Adult Transflectance Forehead, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	




M-LNC Patient Cables





PN	Description	UM	Product Image
2523	<p>M-LNC-1 M-LNCS Series Patient Cable, 1 ft.</p>	EA	
2524	<p>M-LNC-4 M-LNCS Series Patient Cable, 4 ft.</p>	EA	
2525	<p>M-LNC-10 M-LNCS Series Patient Cable, 10 ft.</p>	EA	
2526	<p>M-LNC-14 M-LNCS Series Patient Cable, 14 ft.</p>	EA	


PN	Description	UM	Product Image
2831	<p>M-LNC-MP10 M-LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO₂ M-LNC Dual Key Patient Cable, 10 ft.</p>	EA	
3754	<p>M-LNC-MP4 M-LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO₂ M-LNC Dual Key Patient Cable, 10 ft.</p>	EA	
2655	<p>M-LNC-4-GE M-LNCS to GE Patient Cable, 4 ft.</p>	EA	
2533	<p>M-LNC-10-GE M-LNCS to GE Patient Cable, 10 ft.</p>	EA	

PN	Description	UM	Product Image
2714	<p>M-LNC-SL-10 M-LNCS to Spacelabs Patient Cable, 10 ft.</p>	EA	
2534	<p>M-LNC-04-DB9 M-LNCS Series to DB9 Connector Patient Cable, 4 ft.</p>	EA	
2535	<p>M-LNC-10-DB9 M-LNCS Series to DB9 Connector Patient Cable, 10 ft.</p>	EA	






M-LNCS Adapter Cables

PN	Description	UM	Product Image
2542	<p>M-LNC MAC-SL2 M-LNC Series to Spacelabs Adapter Cable, 2 ft.</p>	EA	
2541	<p>M-LNC MAC-SL M-LNC Series to Spacelabs Adapter Cable, 3 m.</p>	EA	
2537	<p>M-LNC MAC-180 M-LNCS Series to N-180 Adapter Cable, 10 ft.</p>	EA	

PN	Description	UM	Product Image
2538	<p>M-LNC MAC-395 M-LNCS Series to N-395 Adapter Cable, 10 ft.</p>	EA	
2540	<p>M-LNC MAC-GE M-LNCS Series to GE Adapter Cable, 3 m.</p>	EA	
2539	<p>M-LNC MAC CMS M-LNCS Series to Phillips 12-pin Instruments Adapter Cable, 3 m.</p>	EA	
2654	<p>M-LNCS to LNC Adapter Cable M-LNCS Series to LNC Patient Cable, 1.5 ft.</p>	EA	




PN	Description	UM	Product Image
2750	<p>M-LNCS to LNC-MAC Adapter Cable M-LNCS Series to LNC-MAC Adapter Cable, 1.5 ft.</p>	EA	

LNOP Adhesive Sensors


PN	Description	UM	Product Image
1001	<p>LNOP Adt Adult SpO₂ Adhesive Sensors</p> <p>Single patient use 20 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	Box	
1025	<p>LNOP Pdt Pediatric SpO₂ Adhesive Sensors</p> <p>Single patient use 20 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	Box	
1800	<p>LNOP Inf-L Infant SpO₂ Adhesive Sensors</p> <p>Single patient use 20 / box, non-sterile Weight 3 kg – 20 kg Not made with natural rubber latex. Use 2307 replacement tape.</p>	Box	
1798	<p>LNOP Neo-L Neonatal SpO₂ Adhesive Sensors</p> <p>Single patient use 20 / box, non-sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use 2308 replacement tape.</p>	Box	
1651	<p>LNOP NeoPt-L Neonatal SpO₂ Adhesive Sensors</p> <p>Single patient use 20 / box, non-sterile Weight < 1 kg Not made with natural rubber latex. Use 2309 replacement tape.</p>	Box	

LNOP Specialty Sensors

LNOP Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.


PN	Description	UM	Product Image
2358	<p>LNOP Trauma Adult SpO₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight >30 kg Not made with natural rubber latex. Use 2309 replacement tape</p>	Box	
2203	<p>LNOP Newborn Neonatal Neonatal SpO₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight < 3 kg Not made with natural rubber latex. Use 2309 replacement tape</p>	Box	
2204	<p>LNOP Newborn Infant / Pediatric Infant / Pediatric SpO₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight 3 kg – 30 kg Not made with natural rubber latex. Use 2322 replacement tape</p>	Box	

LNOP Cyanotic Sensor




PN	Description	UM	Product Image
1970	<p>LNOP Blue Neonatal / Infant / Pediatric SpO₂ Adhesive Sensors Single patient use 10 / box, non-sterile Weight 2.5 kg – 30 kg Not made with natural rubber latex. Use 1801 replacement tape.</p>	Box	


LNOP Reusable Sensors

PN	Description	UM	Product Image
1269	<p>LNOP DC-I Adult SpO₂ Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
1276	<p>LNOP DCI-P Pediatric SpO₂ Reusable Sensor, 3 ft. 1 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.</p>	EA	
1794	<p>LNOP TC-I SpO₂ Reusable Tip-Clip EAR Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
1793	<p>LNOP TF-I Adult SpO₂ Transflectance Forehead Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
1560	<p>LNOP DC-195 Adult SpO₂ Reusable Sensor 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	





PN	Description	UM	Product Image
1544	<p>LNOP Y-I Multisite SpO₂ Reusable Sensor, 3 ft. 1 / box, non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg Not made with natural rubber latex.</p>	EA	




LNOP Red Patient Cables (for use with Masimo rainbow SET Devices)





PN	Description	UM	Product Image
2058	<p>Red PC-04 LNOP Series 20-pin SpO₂ Patient Cable, 4 ft.</p>	EA	
2059	<p>Red PC-08 LNOP Series 20-pin SpO₂ Patient Cable, 8 ft.</p>	EA	
2060	<p>Red PC-12 LNOP Series 20-pin SpO₂ Patient Cable, 12 ft.</p>	EA	



PN	Description	UM	Product Image
3344	<p>Red 25 PC-12 LNOP Series 25-pin, Patient Cable, 12 ft</p>	EA	

LNOP Patient Cables (for use with Masimo SET Devices)





PN	Description	UM	Product Image
1173	<p>PC04 LNOP Series Patient Cable, 4 ft.</p>	EA	
1005	<p>PC08 LNOP Series Patient Cable, 8 ft.</p>	EA	
1006	<p>PC12 LNOP Series Patient Cable, 12 ft.</p>	EA	
3629	<p>Red 25 PC-4 LNOP Series 25-pin SpO₂ Patient Cable, 4 ft</p>	EA	

PN	Description	UM	Product Image
3631	<p>Red 25 PC-4RA LNOP Series 25-pin, right angle, SpO₂ Patient Cable, 4 ft</p>	EA	
1619	<p>PC04-Ext LNOP Extension Cable, 4 ft.</p>	EA	
2282	<p>LNOP MP12 LNOP to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO₂ LNOP Dual Key Patient Cable, 12 ft.</p>	EA	
1890	<p>PC-12-GE LNOP to GE Patient Cable, 12 ft.</p>	EA	



PN	Description	UM	Product Image
3667	<p>PC-4-GE LNOP to GE Patient Cable, 4 ft.</p>	EA	
3666	<p>PC-8-GE LNOP to GE Patient Cable, 8 ft.</p>		
2288	<p>PC-04-DB9 LNOP Series to DB9 Connector Patient Cable, 4 ft.</p>	EA	
2289	<p>PC-08-DB9 LNOP Series to DB9 Connector Patient Cable, 10 ft.</p>	EA	

PN	Description	UM	Product Image
2290	<p>PC-12-DB9 LNOP Series to DB9 Connector Patient Cable, 12 ft.</p>	EA	
2433	<p>LNOP PC-SL-12 LNOP Spacelabs Patient Cable, 12 ft.</p>	EA	




LNOP Adapter Cables

PN	Description	UM	Product Image
1645	<p>MAC-1 LNOP Series to LNC Adapter Cable, 18 in.</p>	EA	
2858	<p>M-LNCS to PC Adapter Cable M-LNCS Series to LNOP PC Patient Cables, 18 in.</p>	EA	
2014	<p>Dräger Patient Cable For use with Masimo SET LNOP sensors, 2 m.</p>	EA	
1927	<p>MAC-1 LNOP SpO₂ Patient Cable, 12 ft.</p>	EA	




Direct Connect Sensors (for use with Masimo SET Devices)

PN	Description	UM	Product Image
1396	<p>DCSC Adult Reusable Direct Connect SpO₂ Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	
2651	<p>DB-I-dc3 Adult Soft Reusable Direct Connect SpO₂ Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	

Red Direct Connect Sensors (for Masimo rainbow SET Devices)




PN	Description	UM	Product Image
2053	<p>Red DCI-dc3 Adult Reusable Direct Connect™ SpO₂ Sensor, 3 ft. 1 / box, non-sterile Weight > 30kg Not made with natural rubber latex.</p>	EA	
2256	<p>Red DCIP-dc3 Pediatric Reusable Direct Connect SpO₂ Sensor, 3 ft. 1 / box, non-sterile Weight 10kg – 50kg Not made with natural rubber latex</p>	EA	
2643	<p>Red DBI-dc3 Adult Soft Reusable Direct Connect Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.</p>	EA	





Mountable Patient Cable Extensions and Adapter

PN	Description	UM	Product Image
2478	<p>14-Pin to 20-Pin Adapter Cable Connect PXC2 to Masimo SET 20-pin devices, SpO₂ only.</p>	EA	
2371	<p>PXC2-16 Mountable Extension Cable, 16 ft. For use with Masimo SET PC, LNC, M-LNC Patient Cables and DCSC, DC-8, DC-12, DBI-dc3, DBI-dc8 Direct Connect sensors or Masimo 14-pin to 20-pin Adapter Cable.</p>	EA	
2373	<p>PXC2-30 Mountable Extension Cable, 30 ft. For use with Masimo SET® PC, LNC, M-LNC Patient Cables and DCSC, DC-8, DC-12, DBI-dc3, DBI-dc8 Direct Connect sensors or Masimo 14-pin to 20-pin Adapter Cable.</p>	EA	







Accessories






Training Kits





PN	Description	UM	Product Image
3454	<p>SpO₂ ReSposable, Adult / Pediatric Sensor Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 S-DOS 25 • 1 S-DOS 20 • 1 S-DOS 25L • 1 S-DOS 20L • 1 S-DOS 30L • 1 soft wrap armband • 1 cable tether • 1 sensor application poster • Direction for use 	EA	
3455	<p>SpO₂ ReSposable, L Series Sensor Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 S-DOS 30L • 1 S-DOS 25L • 1 S-DOS 20L • 1 S-DOS 25Pt-L • 1 soft wrap armband • 1 cable tether • 1 sensor application poster • 1 neonatal foot with application card • Direction for use 	EA	
2659	<p>M-LNCS Sensor Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 M-LNCS Amtx • 1 M-LNCS Pmtx • 1 M-LNCS Inf • 1 M-LNCS Neo • 1 M-LNCS NeoPt sensor • 1 neonatal foot with application card • Direction for use 	EA	






PN	Description	UM	Product Image
3466	<p>M-LNCS Neonatal and Infant Specialty Sensors Demo Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 M-LNCS Newborn Neonatal • 1 M-LNCS Newborn Infant / Pediatric • 1 M-LNCS Neo • 1 M-LNCS NeoPt • 1 M-LNCS NeoPt-500 • 1 M-LNCS Blue • 1 Neonatal foot with application card • 1 Directions for use 	EA	
2345	<p>LNCS Sensor Training Kit</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 LNCS Actx • 1 LNCS Pctx • 1 LNCS Inf • 1 LNCS Neo • 1 LNCS NeoPt sensor • 1 neonatal foot with application card • Direction for use 	EA	
1282	<p>LNOP Sensors Application Training Packet</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 LNOP Adt • 1 LNOP Pdt • 1 LNOP Inf-L • 1 LNOP Neo-L • 1 LNOP NeoPt-L • 2 cotton swabs • 2 alcohol swabs • 1 neonatal foot with application card • Direction for use 	EA	
3456	<p>Acoustic Respiration Sample Pack</p> <p>Includes:</p> <ul style="list-style-type: none"> • 1 RAS-125 • 1 RAS-125C Sensor 	EA	






Bedside Device Accessories



PN	Description	UM	Product Image
1315	Radical Series Replacement Battery For use with Radical and Radical-7 Blue Screen Handheld.	EA	
2386	Radical-7 Color Screen Replacement Battery For use with Radical-7 Color Screen only.	EA	
23794	Radical-7 Touchscreen Replacement Battery, Rechargeable Lithium Ion	EA	
1395	Radical Series Handheld Lock To secure Radical and Radical-7 Blue Screen Handheld in Radical Docking Stations.	EA	
30534	Radical Series Handheld Lock Key To unlock Radical and Radical-7 Blue Screen Handhelds from Docking Stations.	EA	
1584	Radical Series Power Cord Lock To secure a power cord to Radical Docking Stations 5 / pack	EA	

PN	Description	UM	Product Image
2228	Rad-8 Power Cord Lock	EA	
30475	Radical Series Handheld Case Shown with Blue Screen Handheld Radical. Blue Radical not included.	EA	
3776	Masimo SET M15 Tester Tester for devices containing Masimo SET and connects to M-LNC or rainbow RC patient cables Testers SpO ₂ and PR only	EA	
2367	9-Pin Masimo SET Tester Tester for devices containing Masimo SET and connects to LNC patient cables Tests SpO ₂ and PR only.	EA	
1795	Masimo SET LNOP Tester For Oximeter & Patient Cable SpO ₂ and PR Testing. Tester for devices containing Masimo SET and connects to LNOP PC patient cables. Tests SpO ₂ and PR only.	EA	

PN	Description	UM	Product Image
2368	<p>Masimo rainbow SET Tester</p> <p>Tester for devices containing Masimo SET or Masimo rainbow SET.</p> <p>Tests SpHb, SpCO, SpMet, PVI, PI, SpO₂ and BPM.</p>	EA	
3502	<p>Masimo rainbow SET Tester</p> <p>Tester for devices containing Masimo SET or Masimo rainbow SET, 25-pin instrument connector.</p> <p>Tests SpO₂, PR, SpCO, SpMet and SpHb</p>	EA	
1593	<p>14-Pin Masimo SET Tester</p> <p>For Oximeter SpO₂ and PR Testing</p> <p>Tester for devices containing Masimo SET, 14-pin instrument connector.</p> <p>Tests SpO₂ and PR only.</p>	EA	
3776	<p>Masimo M15 Tester</p> <p>Tester for devices containing Masimo SET and connects to M-LNC or rainbow RC patient cables Testers SpO₂ and PR only</p>	EA	

PN	Description	UM	Product Image
3403	<p>Rainbow Optical Test Sensor Tester for devices containing Masimo rainbow SET.</p> <p>Requires Masimo rainbow RC patient cable tests for SpO2, PR, SpCO, SpMet and SpHb. Not for Clinical Use</p>	EA	
3379	<p>Rainbow Optical Test Sensor, 3 FT. Tester for devices containing Masimo rainbow SET, 20-pin instrument connector.</p> <p>Tests for SpO2, PR, SpCO, SpMet and SpHb.</p> <p>Not For Clinical Use</p>	EA	
3380	<p>Rainbow Optical Test Sensor, 12 FT. Tester for devices containing Masimo rainbow SET, 20-pin instrument connector.</p> <p>Tests for SpO2, PR, SpCO, SpMet and SpHb.</p> <p>Not For Clinical Use</p>	EA	
2227	<p>Rad-8 Mounting Bracket</p>	EA	
2350	<p>Radical-7 Color Handheld Lock Locks Radical-7 Color Handheld to the Docking Station. Used with PN 2351.</p>	EA	

PN	Description	UM	Product Image
2351	<p>Radical-7 Color Handheld Lock Key</p> <p>Used for installation of Radical-7 Color Handheld Lock and unlocking handheld from docking station.</p> <p>Used with PN 2350</p>	EA	
2389	<p>Mini-Demo Connector - 30 Min</p> <p>Used to demo product (best used to display trending information). Includes RAM waveform</p>	EA	
2390	<p>Mini-Demo Connector - 2 Min</p> <p>Used to demo product (best used for short demonstrations). Includes RAM waveform</p>	EA	
36517	<p>Masimo Instrument Configuration Tool (MICT) CD</p> <p>For downloading trend data from Radical-7 (touchscreen only) and MOC-9 technologies through Root. Includes PDF of Directions for Use.</p>	EA	
18313	<p>USB-Null Modem Cable</p> <p>For connecting Root to Windows-based PC for downloading trend data with MICT.</p>	EA	



PN	Description	UM	Product Image
17417	USB-to-Serial Adapter	EA	
24206	Cart, Radical, Masimo	EA	
24208	Accessory, Stack Bracket, Radical-7	EA	No Image
24209	Accessory, Bracket, Cavi Wipes	EA	No Image

Handheld Device Accessories








PN	Description	UM	Product Image
2207	<p>Water Resistant Protective Carrying Case for Battery Operated Portable Handheld, Black</p> <p>This can be used with Rad-5/5v.</p>	EA	
2208	<p>Water Resistant Protective Carrying Case for Battery Operated Portable Handheld, Red</p> <p>This can be used with Rad-5/5v.</p>	EA	
10368	<p>Radical-7 Ambulation Pouch</p> <p>Pack of 10</p>	EA	No image








PN	Description	UM	Product Image
1909	<p>Data Transfer Download Cable</p> <p>Used for connecting Rad-5 (with 14-pin connector) to Personal Computer to enable download of trending data.</p>	EA	
2063	<p>Data Transfer Download Cable, Red – 20-Pin</p> <p>Used for connecting Rad-57's and 20-pin Rad-5s to Personal Computer to enable download of trending data.</p>	EA	
1908	<p>TrendCom, Trend Download Software</p> <p>Used with appropriate Data Transfer Download cable to initiate download and de-compression of Rad-5 and Rad-57 trending data to a Personal Computer. Used with serial cable to initiate download and decompression of Radical-7, Rad-87 and Rad-8 trending data.</p>	EA	
1842	<p>Handheld Boot, Gray</p> <p>For use with Rad-5, Rad-5v and all Rad-57 Handhelds.</p>	EA	
1980	<p>Handheld Boot, Yellow</p> <p>For use with Rad-5, Rad-5v and all Rad-57 Handhelds.</p>	EA	








PN	Description	UM	Product Image
1981	<p>Handheld Boot, Red For use with Rad-5, Rad-5v and all Rad-57 Handhelds.</p>	EA	
1982	<p>Handheld Boot, Orange For use with Rad-5, Rad-5v and all Rad-57 Handhelds.</p>	EA	
2097	<p>Handheld Boot, Royal Blue For use with Rad-5, Rad-5v and all Rad-57 Handhelds.</p>	EA	
2098	<p>Handheld Boot, Light Blue For use with Rad-5, Rad-5v and all Rad-57 Handhelds.</p>	EA	
2099	<p>Handheld Boot, Pink For use with Rad-5, Rad-5v and all Rad-57 Handhelds.</p>	EA	



PN	Description	UM	Product Image
2357	<p>Light Shield For use in strobing environment and high sunlight intensity areas. 5 / Pack</p>	EA	
2815	<p>Masimo Ambient Shield Accessory For use with rainbow Sensors. 10 / pack</p>	EA	

Sensor Accessories

PN	Description	UM	Product Image
2624	<p>Replacement Tapes rainbow R1 20L and R20-L</p> <p>For use with rainbow R1 20L and R20-L</p> <p>50 / box, non-sterile</p> <p>Not made with natural rubber latex.</p>	Box	
2623	<p>Replacement Tapes For rainbow R1 25L and R25-L</p> <p>For use with rainbow R1 25L and R25-L</p> <p>50 / box, non-sterile</p> <p>Not made with natural rubber latex.</p>	Box	
1801	<p>Replacement Tapes For LNOP Blue Sensors</p> <p>For use with LNOP Blue</p> <p>100 / pack, non-sterile</p> <p>Not made with natural rubber latex.</p>	Box	
2307	<p>Replacement Tapes for Inf Series Sensors</p> <p>For use with LNCS/M-LNCS Inf/Inf-3 or LNOP Inf-L</p> <p>102 / box, non-sterile</p> <p>Not made with natural rubber latex.</p>	Box	
2308	<p>Replacement Tapes For Neo Series Sensors</p> <p>For use with LNCS / M-LNCS Neo/Neo-3 or LNOP Neo-L</p> <p>102 / box, non-sterile</p> <p>Not made with natural rubber latex.</p>	Box	
2309	<p>Replacement SofTouch Wraps for NeoPt, NeoPt-3, Trauma and Newborn Neonatal Sensors</p> <p>For use with LNCS/M-LNCS NeoPt/NeoPt-3, LNOP NeoPt-L and LNCS/M-LNCS/LNOP Trauma and Newborn Neonatal</p> <p>10 / pack, non-sterile</p> <p>Not made with natural rubber latex.</p>	Pack	
2322	<p>Replacement SofTouch Wraps for NeoPt-500 and Newborn Infant / Pediatric Sensors</p> <p>For use with LNCS/M-LNCS NeoPt-500 and LNCS/M-LNCS/LNOP Newborn Infant/Pediatric</p> <p>10 / pack, non-sterile</p> <p>Not made with natural rubber latex.</p>	Pack	







PN	Description	UM	Product Image
2822	<p>Replacement Tape for M-LNCS Blue For use with M-LNCS Blue 50 / pack, non-sterile Not made with natural rubber latex.</p>	Pack	
1596	<p>CleanShield® Multisite Wrap For use with Y-I Sensors Adult / Pediatric / Neonatal 100 / box, non-sterile Not made with natural rubber latex.</p>	Box	
1597	<p>Standard Multisite Wrap For use with Y-I Sensors Adult / Pediatric / Neonatal 100 / box, non-sterile Not made with natural rubber latex.</p>	Box	
1598	<p>Standard Petite Wrap For use with Y-I Sensors Adult slender digit / Ped / Neonatal adhesive 100 / box, non-sterile Not made with natural rubber latex.</p>	Box	
1600	<p>Accessory Adhesive Squares For use with Y-I Sensors 144 / pack, non-sterile</p>	Pack	
1903	<p>Transflectance Sensor Adhesive Pad For use with TF-I® Sensors 20 / pack, non-sterile</p>	Pack	
16179	<p>Tele Plast Philips SpO₂ Connector Plug Compatible with Philips MMS module.</p>	EA	

PN	Description	UM	Product Image
1602	<p>Foam Wraps for LNOP / LNCS® / M-LNCS Y-I Sensors</p> <p>For use with Y-I Sensors</p> <p>12 / pack, non-sterile</p> <p>Not made with natural rubber latex.</p>	Pack	
1601	<p>Accessory Clothing Clips</p> <p>For securing reusable sensor cables to patient clothing</p> <p>5 / pack, non-sterile</p>	Pack	
3444	<p>Masimo® Soft Wrap Armband</p> <p>Used to secure the sensor cable along the patients arm.</p> <p>20 / pack, non-sterile</p> <p>Not made with natural rubber latex.</p>	Pack	
3929	<p>Masimo Forehead Sensor Headbands</p> <p>For use with forehead sensors</p> <p>20 / box, non-sterile</p> <p>Not made with natural rubber latex</p>	Box	
1608	<p>Masimo Headband</p> <p>For use with TF-I Sensors</p> <p>5 / box, non-sterile</p> <p>Not made with natural rubber latex.</p>	Box	
2215	<p>Masimo Headband</p> <p>For use with TF-I Sensors</p> <p>12 / box, non-sterile</p> <p>Not made with natural rubber latex.</p>	Box	
21514	<p>LNC Replacement Cover</p> <p>For use with all LNC patient cables</p> <p>25 covers and 50 pins.</p>	Pack	

PN	Description	UM	Product Image
3679	<p>rainbow Cable Tether</p> <p>Secures rainbow ReSposable ROS to patient cable</p> <p>20 / pack, non-sterile</p> <p>Not made with natural rubber latex.</p>	Pack	
3443	<p>Masimo Cable Tether</p> <p>Secures SpO₂ ReSposable ROS or reusable sensor to patient cable.</p> <p>20 / pack, non-sterile</p> <p>Not made with natural rubber latex.</p>	Pack	








Pronto-7 Spot Check Pulse CO-Oximeter Accessories









PN	Description	UM	Product Image
3382	Boot, Pronto-7 Red	EA	
3383	Boot, Pronto-7 Yellow	EA	
3384	Boot, Pronto-7 Green	EA	
3385	Boot, Pronto-7 Light Blue	EA	
3386	Boot, Pronto-7 Pink	EA	







PN	Description	UM	Product Image
3387	Boot, Pronto-7 Dark Blue	EA	
3388	Boot, Pronto-7 Black	EA	
3389	Boot, Pronto-7 Orange	EA	
34983	Finger Gauge, Pronto-7	EA	
16782	Travel Bag, Pronto-7	EA	
16868	Bracket, Pronto-7	EA	



SatShare Interface Cables

Note: Do not interface the Radical-7 with any monitor or module not specifically listed in this chart.

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1595	AT01	Airshields <ul style="list-style-type: none"> Athena Monitor with Nellcor SpO₂ (Type 9200, Type 9267) 	EA	
1324	DO01	Datex-Ohmeda <ul style="list-style-type: none"> AS/3 Patient Monitor (M-NSAT) System 5 Compact Monitor (M-NSAT) 	EA	
1326	DO02	Datex-Ohmeda <ul style="list-style-type: none"> AS/3 Patient Monitor (M-NESTPR-00-00 / M-NESTR-00-00) CardioCap / 5 Monitor Light Monitor System 5 Compact Monitor (M-NESTPR-00-00 / M-NESTR-00-00) 	EA	
2023	DO05	Datex-Ohmeda <ul style="list-style-type: none"> System 5 Compact Monitor (M-NESTPR-00-01] 	EA	
2223	DO06	Datex-Ohmeda <ul style="list-style-type: none"> System 5 Compact Monitor (M-PRESTN..01) 	EA	
2969	DO07	Datex-Ohmeda <ul style="list-style-type: none"> System 5 Anesthesia Monitor (E-PRESTN) 	EA	
1325	DS01	Corometrics <ul style="list-style-type: none"> Maternal / Fetal Monitor Model 118 Maternal / Fetal Monitor Model 120 		


PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1325	DS01	Critikon <ul style="list-style-type: none"> Dinamap Monitor (117310) 	EA	
1325	DS01	Datascope <ul style="list-style-type: none"> Expert Monitor with Nellcor SpO₂ (HS-500) Accutor Monitor with Nellcor SpO₂ Passport Monitor with Nellcor SpO₂ 	EA	
1325	DS01	North American Dräger <ul style="list-style-type: none"> Narcomed 4 Anesthesia Machine (MP202) Vitalert 3200 Monitor (203) 	EA	
1325	DS01	Medical Data Electronics <ul style="list-style-type: none"> Escort Monitor 	EA	
1325	DS01	Witt Biomedical <ul style="list-style-type: none"> Frontend05 Patient Monitor 	EA	
1365	DS03	Dynascope <ul style="list-style-type: none"> Fukuda Denshi DS-5300W (HS-500) 	EA	
2979	DS04	Dynascope <ul style="list-style-type: none"> Fukuda Denshi DS-5400 (HS-500) 	EA	
1357	HP03	Philips / Agilent / HP <ul style="list-style-type: none"> Omnicare / Viridia Model 24 / 26 Series Monitor System (M1204A) Series 50XM Fetal Monitor 	EA	

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1321	MQ01	<p>GE / Marquette</p> <ul style="list-style-type: none"> Solar 7000 / 8000 Series Monitor (Tram 250SL, Tram 450SL, Tram 650SL, Tram 850A&SL, and SpO2 Solar Parameter Module) Dash3000 Monitor 	EA	
1789	MS01	<p>Datex-Ohmeda</p> <ul style="list-style-type: none"> S5 (E-Masimo-00) <p>Dräger</p> <ul style="list-style-type: none"> Infinity Delta (SmartPOD with Masimo SET) <p>GE / Marquette</p> <ul style="list-style-type: none"> Solar 8000 (Tram 451M, Tram 851M, Masimo SET single parameter module) <p>Philips</p> <ul style="list-style-type: none"> IntelliVue MP40 Monitor (MMS Module [Masimo SET]) <p>GE</p> <ul style="list-style-type: none"> Carescape B850 (Patient Data Module [Masimo SET]) 	EA	
1367	NK01	<p>Nihon Kohden</p> <ul style="list-style-type: none"> Model BMS-8301A/MV-831RA (AL-800PA) Model BMS-8502A (AL-800PA) 	EA	
1368	NK02	<p>Nihon Kohden</p> <ul style="list-style-type: none"> BSM-2301 (Type 9266) BSM-5105 (Type 9266) 	EA	
1646	OH02	<p>Hellige</p> <ul style="list-style-type: none"> SMU 611 Monitor (SAT) <p><i>Not for Sale in the United States</i></p>	EA	
1323	SL01	<p>SpaceLabs Medical</p> <ul style="list-style-type: none"> PC Bedside Monitor (90303) [90465] PC2 Central Monitor (90315, 90316) [90465] <p>SpaceLabs Medical</p> <ul style="list-style-type: none"> Adult/Neonatal Monitor (90621) 	EA	

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1362	SL02	SpaceLabs Medical <ul style="list-style-type: none"> • Ultraview 1700 Bedside Monitor (90387) [90496] • Ultraview UCW Monitor (90385) [90496] 	EA	
1322	SM01	Siemens Medical Systems <ul style="list-style-type: none"> • Sirecust 1281 Monitor(99731472E2529) • SC 6002XL, SC7000 	EA	

Note: Do not interface the Radical-7 with any monitor or module not specifically listed in this chart.

SatShare Extension Cable

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)*	UM	Host Connector
1605		SatShare Extension Cable	EA	

Index

A

Accessories • 99
Airshields • 1

B

Bedside Device Accessories • 101

C

Capnography • 1, 15
Corometrics • 1
Critikon • 1

D

Datascope • 1
Datex-Ohmeda • 1
Direct Connect Sensors (for use with Masimo SET Devices) • 95
Dräger • 1

E

EMMA • 1, 16

G

GE / Marquette • 1
GE Healthcare Compatible Devices • 24

H

Handheld Device Accessories • 107
Hellige • 1

I

iMagox • 1
IntelliVue • 1
ISA • 1
ISA CO2 • 15
iSpO2 • 1
iSpO2@ • 25

L

LNC Patient Cables (for use with Masimo SET Devices) • 63
LNCS Adapter Cables • 67
LNCS Adhesive Sensors • 51
LNCS Adhesive Sensors with Replaceable Tapes • 55
LNCS Red Patient Cables (for use with Masimo rainbow SET Devices) • 62

LNCS Reusable Sensors • 60
LNCS Specialty Sensors • 57
LNCS Sterile Sensors • 58
LNOP Adapter Cables • 94
LNOP Adhesive Sensors • 64
LNOP Cyanotic Sensor • 85
LNOP Patient Cables (for use with Masimo SET Devices) • 90
LNOP Red Patient Cables (for use with Masimo rainbow SET Devices) • 88
LNOP Reusable Sensors • 86
LNOP Specialty Sensors • 85

M

Masimo Patient SafetyNet™ System • 17
Masimo rainbow SET Parameter Options • 11
Masimo rainbow® Pulse CO-Oximetry Devices • 8
Masimo SET ReSposable S-DOS™ Series • 44
Masimo SET ReSposable S-ROS 12U Series • 47
Masimo SET ReSposable S-ROS™ 3U Series • 46
Masimo SET® Pulse Oximetry Devices • 5
Medical Data Electronics • 1
M-LNC Patient Cables • 78
M-LNCS Adapter Cables • 81
M-LNCS Adhesive Sensors with Replaceable Tapes • 72
M-LNCS Cyanotic Sensors • 75
M-LNCS Reusable Sensors • 76
M-LNCS Specialty Sensors • 74
M-LNCS™ Adhesive Sensors • 69
Monitoring Devices • 5
Mountable Patient Cable Extensions and Adapter • 97

N

Nellcor • 1
Nihon Kohden • 1

P

Philips • 1
Philips® Compatible Devices • 23
Pronto • 1
Pronto-7 • 1
Pronto-7 Reusable Spot Check Sensors • 38
Pronto-7 Spot Check Pulse CO-Oximeter Accessories • 115
PVI • 1

R

Rad-5 • 1
Rad-57 • 1
Rad-5v • 1
Rad-8 • 1

Masimo Product Catalog 2015 - Q1

Rad-87 • 1
Radical-7 • 1
rainbow • 1
rainbow Acoustic Monitoring™ Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRA) • 29
rainbow Adhesive Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet) • 30
rainbow Adhesive Sensors (for rainbow SET Devices with SpO2, SpCO, and SpMet) • 32
rainbow Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet) • 40
rainbow Patient Cables (for rainbow SET Devices) • 41
rainbow ReSposable™ Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet) • 33
rainbow Reusable Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet) • 39
RDS • 1
Red Direct Connect Sensors (for Masimo rainbow SET Devices) • 96
ReSposable MAC S-ROS 12U Series • 49
Root • 1
RRA • 1
RRp • 1

S

SatShare Extension Cable • 120
SatShare Interface Cables • 117
SedLine® Brain Function Monitoring • 14
Sensor Accessories • 111
Sensors and Cables • 29
Shield • 1
Siemens Medical Systems • 1
SpaceLabs • 1
SpCO • 1
Specialty Monitoring and Connectivity Platform • 13
SpHb • 1
rainbow Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO2, SpMet) • 35
SpMet • 1

T

Tape • 1
Tether • 1
Training Kits • 99

U

uSpO2 • 1

V

VueLink • 1

W

Witt Biomedical • 1
Wrap • 1



www.masimo.com

LAB-3422AA 0215