

PENTAX

Ureteroreno Cysto Color Scope



ECY-1570K

Superb Optics, Easy Handling

ECY-1570K

Ureteroreno Cysto

Color Scope




ECY-1570K, PENTAX Video Cystoscope, has 5.5mm distal end outer diameter and 5.1mm insertion tube outer diameter. Connected with Processor, EPK-700, EPK-1000, EPK-p & EPK-i.
Hygiene concept: All root braces (3 parts) are sealed up.

Specifications

Field of View	120° (in air) 80° (in water)	
Depth of Field	3 ~ 50mm	
Tip Deflection	Up	210°
	Down	130°
Rigid Distal Width	ø5.5mm	
Insertion Tube Width	ø5.1mm	
Minimum Instrument Channel Width*	ø2.0mm	
Working Length	400mm	
Total Length	660mm	



* There is no guarantee that instruments selected solely using this minimum instrument channel width will be compatible in combination.



Outstanding performance characteristics



- Large Screen, High-Resolution
- Control Unit Designed for Ease of Operation
- Large 2.0mm channel
- PVE Connector can be rotated within 180° range



PENTAX™

HOYA Corporation PENTAX Tokyo Office

2-36-9 Maeno-cho, Itabashi-ku, Tokyo 174-8639 Japan TEL: +81-3-3960-5155 FAX: +81-3-5392-6724

PENTAX Medical Company (A Division of PENTAX of America, Inc.)

102 Chestnut Ridge Road, Montvale, New Jersey 07645 U.S.A. Toll-free: +1-800-431-5880 TEL: +1-201-571-2300 FAX: +1-201-391-4189

PENTAX Europe GmbH

Julius-Vosseler-Strasse 104, 22527 Hamburg, Germany TEL: +49-40-561 920 FAX: +49-40-560 42 13

PENTAX MEDICAL SINGAPORE PTE. LTD.

60 Alexandra Terrace, #10-08 The Comtech Singapore 118502 TEL: +65-6271-1669 FAX: +65-6271-1691

URL www.medical.pentax.co.jp

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician or other appropriately licensed medical professional.

*Specification and design are subject to change without notice due to product improvement.