

Masimo Product Catalog

2015 - Q1



This catalog contains Masimo proprietary material and is provided for information purposes only. No part of this catalog may be reproduced or retransmitted without Masimo's prior written consent.

The purpose of the product catalog is to provide a general listing of Masimo products, and it is not intended to represent products that are offered for sale. Please see local Masimo representatives for the commercially available products in specific regions/countries with the appropriate regulatory approvals/clearances.


Caution: Federal law restricts this device to sale by or on the order of a physician.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Masimo Corporation
40 Parker
Irvine, CA 92618, USA
Tel.: 949-297-7000
Fax.: 949-297-7001
www.masimo.com



EU authorized representative for Masimo Corporation:



MDSS GmbH
Schiffgraben 41
D-30175 Hannover, Germany

© 2015 by Masimo Corporation.

Patents: www.masimo.com/patents.htm

Masimo, , 3D Alarm, Blue, Cleanshield, DCI, E1, FastSat, iSpO₂, LNCS, LNOP, MX-1, Pronto, Pronto-7, Pulse CO-Oximeter, PVI, Rad-5, Rad-8, Rad-57, Rad-87, Radical, Radical-7, rainbow, rainbow ReSposable, RRa, SatShare, SedLine, SET, SpCO, SpHb, SpMet, TF-I, and uSpO₂ are registered trademarks of Masimo Corporation.

3D Desat Index Alarm, Iris, M-LNCS, MOC-9, Phasein, Rad-5v, Rad-57, rainbow Acoustic Monitoring, RDS, Root, ROS, RRp, S-DOS, and S-ROS are trademarks of Masimo Corporation.

Apple, iPad, iPhone, iPod and iPod touch are registered trademarks of Apple Inc. registered in the U.S. and other countries.

Lightning is a trademark of Apple Inc.

"Made for iPod," "Made for iPhone," and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, or iPad may affect wireless performance.

The use of the trademarks PATIENT SAFETYNET and PSN is under license from University HealthSystem Consortium.

All other trademarks and registered trademarks are property of their respective owners.

The devices provided in the catalog may not have been licensed in accordance with applicable Federal Laws.

CONFIDENTIAL AND PROPRIETARY

Contents

Monitoring Devices	5
Masimo SET® Pulse Oximetry Devices.....	5
Masimo rainbow® Pulse CO-Oximetry Devices.....	8
Masimo rainbow SET Parameter Options.....	11
Specialty Monitoring and Connectivity Platform	13
SedLine® Brain Function Monitoring	14
Capnography	15
Masimo Patient SafetyNet™ System.....	17
Philips® Compatible Devices	23
GE Healthcare Compatible Devices	24
iSpO2®.....	25
Sensors and Cables	29
rainbow Acoustic Monitoring™ Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRA)	29
rainbow Adhesive Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet)	30
rainbow Adhesive Sensors (for rainbow SET Devices with SpO2, SpCO, and SpMet)	32
rainbow ReSposable™ Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet)	33
rainbow Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO2, SpMet; SpHb required)	35
Pronto-7 Reusable Spot Check Sensors	38
rainbow Reusable Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet)	39
rainbow Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet).....	40
rainbow Patient Cables (for rainbow SET Devices)	41
Masimo SET ReSposable S-DOS™ Series.....	44
Masimo SET ReSposable S-ROS™ 3U Series	46
Masimo SET ReSposable S-ROS 12U Series.....	47
ReSposable MAC S-ROS 12U Series	49
LNCS Adhesive Sensors.....	51
LNCS Adhesive Sensors with Replaceable Tapes	55
LNCS Specialty Sensors	57
LNCS Sterile Sensors	58
LNCS Reusable Sensors	60
LNCS Red Patient Cables (for use with Masimo rainbow SET Devices)	62
LNC Patient Cables (for use with Masimo SET Devices).....	63
LNCS Adapter Cables	67
M-LNCS™ Adhesive Sensors	69
M-LNCS Adhesive Sensors with Replaceable Tapes	72
M-LNCS Specialty Sensors	74
M-LNCS Cyanotic Sensors.....	75
M-LNCS Reusable Sensors	76
M-LNC Patient Cables.....	78
M-LNCS Adapter Cables.....	81

Masimo Product Catalog 2015 - Q1

LNOP Adhesive Sensors.....	84
LNOP Specialty Sensors.....	85
LNOP Cyanotic Sensor	85
LNOP Reusable Sensors	86
LNOP Red Patient Cables (for use with Masimo rainbow SET Devices).....	88
LNOP Patient Cables (for use with Masimo SET Devices)	90
LNOP Adapter Cables.....	94
Direct Connect Sensors (for use with Masimo SET Devices).....	95
Red Direct Connect Sensors (for Masimo rainbow SET Devices).....	96
Mountable Patient Cable Extensions and Adapter.....	97
Accessories	99
Training Kits.....	99
Bedside Device Accessories	101
Handheld Device Accessories.....	107
Sensor Accessories	111
Pronto-7 Spot Check Pulse CO-Oximeter Accessories.....	115
SatShare Interface Cables	117
SatShare Extension Cable	120
Index.....	121

Monitoring Devices

Masimo SET® Pulse Oximetry Devices

PN	Description	UM	Product Image
9190	Rad-8®, Horizontal Comes with either an LNOP® or LNCS® Patient Cable and a power cord.	EA	
9191	Rad-8, Vertical Comes with either an LNOP or LNCS Patient Cable and a power cord.	EA	
9196	Rad-5® Pulse Oximeter Comes with a Red DCI®-DC3 Sensor and choice of patient cable, a choice of available Colors of Handheld Boots and a choice of user manual. Shown in optional protective boot and with DCSC Sensor.	EA	

PN	Description	UM	Product Image
9198	Rad-5 Pulse Oximeter with Sample Pack Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of user manual. Shown in optional protective boot and with DCSC Sensor.	EA	
9214	Rad-5 Pulse Oximeter (International Only) Comes with an LNOP DCSC Sensor, a choice of available Colors of Handheld Boots and a choice of user manual. Shown in optional protective boot and with DCSC Sensor.	EA	
9197	Rad-5v™ Pulse Oximeter Comes with a Red DCI-DC3 Sensor and choice of patient cable, a choice of available colors of Handheld Boots and a choice of user manual. Shown in optional protective boot and with DCSC Sensor.	EA	

PN	Description	UM	Product Image
9199	Rad-5v Pulse Oximeter with Sample Pack Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of user manual. Shown in optional protective boot and with DCSC Sensor.	EA	
9215	Rad-5v Pulse Oximeter (International Only) Comes with an LNOP DCSC Sensor, a choice of available Colors of Handheld Boots and a choice of user manual. Shown in optional protective boot and with DCSC Sensor.	EA	

Masimo rainbow® Pulse CO-Oximetry Devices

PN	Description	UM	Product Image
9500	Radical-7® Handheld Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNOP, Red LNCS or rainbow Patient Cable. Can be configured with optional rainbow parameters and features: • Total Hemoglobin (SpHb®) • Oxygen Content (SpOC) • Methemoglobin (SpMet®) • Carboxyhemoglobin (SpCO®) • Pleth Variability Index (PVI®) • Oxygen Reserve Index (ORI®) This unit ships standard with a 802.11 wireless radio. Can be loaded with Respiration Rate from Pleth (RRp™) and Adaptive Threshold Alarms (ATA™) in cleared countries only. Check with customer service on availability. Optional rainbow parameters and features are sold separately and can be added to this unit upon request. The Radical-7 is compatible with all Radical® Docking Stations. <i>Available in countries where cleared.</i> <i>ORI® parameter requires Root Docking Station</i>	EA	
9021	Radical® Docking Station, RDS™-1 Comes with a choice of power cord and serial output connector, Analog output / nurse call connector and SatShare® cable connector.		
9023	Radical® Docking Station, RDS-3 limited communication capabilities Comes with choice of power cord and serial output connector and Analog output / nurse call connector.		

PN	Description	UM	Product Image
9226	Rad-87®, Horizontal Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters. Includes either Red LNCS, Red LNOP or rainbow RC Patient Cable. Can be configured with optional rainbow parameters and features: • Total Hemoglobin (SpHb) • Oxygen Content (SpOC) • Methemoglobin (SpMet) • Carboxyhemoglobin (SpCO) • Pleth Variability Index (PVI) • Acoustic Respiration (RRa™) • Adaptive Threshold Alarm (ATA): <i>Available in countries where cleared.</i> Comes with serial output connector and Analog output / Nurse Call connector. rainbow Acoustic Monitoring capable. Optional rainbow parameters and features are sold separately and can be added to this unit upon request. <i>Available in countries where cleared.</i>	EA	 A horizontal white pulse oximeter with a black display showing three large numbers: 97, 76, and 19. It has several buttons on the right side and a red cable connector on the left.
9228	Rad-87, Vertical Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters. Includes either a Red LNCS, Red LNOP or rainbow RC Patient Cable. Can be configured with optional rainbow parameters and features: • Total Hemoglobin (SpHb) • Oxygen Content (SpOC) • Methemoglobin (SpMet) • Carboxyhemoglobin (SpCO) • Pleth Variability Index (PVI) • Acoustic Respiration (RRa) • Adaptive Threshold Alarm (ATA): <i>Available in countries where cleared.</i> Comes with serial output connector and Analog output / Nurse Call connector. rainbow Acoustic Monitoring capable. Optional rainbow parameters and features are sold separately and can be added to this unit upon request. <i>Available in countries where cleared.</i>	EA	 A vertical white pulse oximeter with a black display showing three large numbers: 97, 76, and 19. It has several buttons on the front and a red cable connector at the bottom.

PN	Description	UM	Product Image
9221	Pronto-7® Spot Check Pulse CO-Oximeter Includes SpO₂, Pulse Rate, Perfusion Index (PI) and SpHb parameters. rainbow 4D DC Sensor not included. rainbow 4D DC Reusable Sensor must be ordered separately.	EA	
9167	Pronto® Spot Check Pulse CO-Oximeter® Includes SpO₂, Pulse Rate, Perfusion Index (PI) and SpHb parameters. rainbow Sensor not included. rainbow Sensor must be ordered separately.	EA	
9216-U 9216	Rad-57™ Configurable Pulse CO-Oximeter Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters. The Rad-57 can ship with or be field upgraded with: • Carboxyhemoglobin (SpCO) • Methemoglobin (SpMet) • Total Hemoglobin (SpHb) • Oxygen Content (SpOC) • PVI When shipped with rainbow parameters, rainbow Sensor must be purchased separately. rainbow Sensors not included. 9216-U is for SpO₂ only Devices. Refer to Handheld Device accessories for available colors. Shown with rainbow Adult Reusable Sensor and in optional protective boot.	EA	

Masimo rainbow SET Parameter Options

PN	Description	UM	Product Image
9081	Field Upgrader Used to load optional rainbow parameters and features onto compatible monitors. Current available options include: • Total Hemoglobin (SpHb[®]) • Oxygen Content (SpOC) • Methemoglobin (SpMet) • Carboxyhemoglobin (SpCO[®]) • Pleth Variability Index (PVI[®]) • Oxygen Reserve Index (ORI[®]) • 3D Desat Index Alarm • 3D PI Delta Alarm Each optional parameter and feature sold separately. Maximum quantity of 25 of any individual parameter or feature per Field Upgrader. Oxygen Reserve Index (ORI[®]) Not for sale in the USA. For Export Only. Includes: • 1 Field Upgrader • 1 to 5 operator's manuals • Parameter upgrade indicator sticker(s) • Postage-paid envelope to return Field Upgrader	EA	
2630	Total Hemoglobin (SpHb) parameter Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field	EA	No image
2296	Carboxyhemoglobin (SpCO) parameter Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.	EA	No image
2297	Methemoglobin (SpMet) parameter Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.	EA	No image
2313	Pleth Variability Index (PVI) Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.	EA	No image
3950	Oxygen Reserve Index (ORI™) Can be ordered on a Radical-7 when used with Root. Not for sale in the USA. For Export Only.	EA	No image
3951	Oxygen Reserve Index (ORI™), Fixed Term License Can be ordered on a Radical-7 when used with Root. Not for sale in the USA. For Export Only.	EA	No image
2940	Acoustic Respiration (RRa™) Can be ordered on a Rad-87 (PN 9226 and 9228) or on a Field Upgrader (PN 2926) to upgrade a compatible monitor in the field.	EA	No image

PN	Description	UM	Product Image
3518	Respiration Rate from Pleth (RRp) Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field. In countries where RRp is cleared.	EA	No image
3519	Respiration Rate from Pleth (RRp), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field. <i>Available in countries where cleared.</i>	EA	No image
2631	Carboxyhemoglobin (SpCO), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field	EA	No image
2632	Methemoglobin (SpMet), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.	EA	No image
2633	Total Hemoglobin (SpHb), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field	EA	No image
2636	PI Delta 3D Alarm, Fixed Term License Can be ordered on a Radical-7 on a Field Upgrader (PN 9081) to upgrade a compatible monitor in the field.	EA	No image
2637	Desat Index 3D Alarm, Fixed Term License Can be ordered on a Radical-7 on a Field Upgrader (PN 9081) to upgrade a compatible monitor in the field.	EA	No image
3402	Adaptive Threshold Alarm (ATA), Fixed Term License Can be ordered on a Rad-87 (PNs 9226 and 9228). <i>Available in countries where cleared.</i>	EA	No image
3325	rainbow Acoustic Monitoring (RRa), Fixed Term License Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.	EA	No image
2341	Service Tool, MX-1[®] Upgrades up to 100 instruments.	EA	No image
2897	Service Tool, MX-3 Upgrades up to 100 instruments.	EA	No image
3351	Service Tool, MX-3 with ATA Upgrades up to 100 instruments. <i>Available in countries where cleared.</i>	EA	No image

PN	Description	UM	Product Image
2340	Demonstration Tool Used to exercise monitors alarm default setting and demonstrate optional parameters (SpHb, SpCO, SpMet and PVI) when optional parameters are installed on a compatible monitor.	EA	 A white, rectangular device with a red top handle. The front face features the 'rainbow' logo in multi-colored text, the words 'Demonstration Tool' in black, and the Masimo logo at the bottom right.

Specialty Monitoring and Connectivity Platform

PN	Description	UM	Product Image
9515	Root™ - Screen Size: 10.1 in (25.65 cm) - Resolution: 1280 x 800 Pixels - Color: 24-bit RGB - Display Type: Backlit Active Matrix TFT LCD - Touchscreen: Multi-Touch P-Cap Provides Ethernet connection, nurse call connector, USB 2.0 Ports, MOC-9™ ports, and Iris® Connectivity Ports. Comes with choice of power cords. <i>Available in countries where cleared.</i>	EA	 A sleek, black and silver medical monitor with a large, tilted touchscreen display. The device has a modern, angular design with a silver bezel and a black body.

SedLine® Brain Function Monitoring

PN	Description	UM	Product Image
9513	SedLine Kit For use with Root and SedLine Sensors. Includes: 1 ISA CO2 MOC-9 Module 1 ISA CO2 Bracket Available in countries where cleared.	EA	
9655	SedLine Demo MOC-9 Module For use with Root. Does not include SedLine Patient Cable. Not for clinical use - for product demonstration only.	EA	
23967	Replacement SedLine MOC-9 Module For use with Root and SedLine Sensors. Available in countries where cleared.	EA	
3637	Replacement Sedline Patient Cable For use with Sedline MOC-9 and and SedLine Sensors. Available in countries where cleared.	EA	

PN	Description	UM	Product Image
2479	SedLine Sensor Single patient use 25 / box, non-sterile Not made with natural rubber latex. <i>Available in countries where cleared.</i>	Box	

Capnography

ISA CO2

PN	Description	UM	Product Image
3727	ISA CO₂ Kit Includes: 1 ISA CO₂ MOC-9 Module 1 ISA CO₂ Bracket <i>Available in countries where cleared.</i>	EA	
9685	ISA OR+ Kit 1 ISA OR+ MOC-9 Module 1 ISA OR+ Bracket	EA	No Image
3738	Nomoline Adapter 25 / Box, Non-sterile <i>Available in countries where cleared.</i>	Box	
3739	Nomoline Airway Adapter Set 20 / Box, Non-sterile <i>For use with intubated patients</i> <i>Available in countries where cleared.</i>	Box	

EMMA

PN	Description	UM	Product Image
9632	EMMA Capnometer Kit Includes: • 1 EMMA Capnometer mmHg • 1 EMMA Capnometer Pouch • 1 EMMA Capnometer Lanyard	EA	
9633	EMMA Capnometer Kit Includes: • 1 EMMA Capnometer kPa • 1 EMMA Capnometer Pouch • 1 EMMA Capnometer Lanyard	EA	
3639	EMMA Capnograph Kit Includes: • 1 EMMA Capnograph mmHg • 1 EMMA Pouch • 1 EMMA Lanyard	EA	
3678	EMMA Capnograph Kit Includes: • 1 EMMA Capnograph kPa • 1 EMMA Pouch • 1 EMMA Lanyard	EA	
17448	EMMA Adult / Pediatric Airway Adapter 25 / Box, Non-sterile	Box	
17449	EMMA Infant Airway Adapter 10 / Box, Non-sterile	Box	

Masimo Patient SafetyNet™ System

PN	Description	UM	Product Image
9095	SafetyNet™ System Configurable for customer needs.	System	
23468	Patient SafetyNet View Station with Software (V 3.1.0.2)	EA	
24479	Patient SafetyNet View Station with Software (V 4.X.X.X)	EA	
24252	Kit SafetyNet Appliance 1UM	EA	
24254	Kit SafetyNet Appliance 2UM	EA	
24487	Kit SafetyNet Appliance 1UM and View	EA	

PN	Description	UM	Product Image
24489	Kit SafetyNet Appliance 2UM and View	EA	
23736	Halo Appliance	EA	
3623	Upgrade to V 4.X.X.X	EA	No Image
3479	Upgrade to Patient Trend Analysis Report – First Appliance	EA	No Image
3480	Upgrade to Patient Trend Analysis Report – Additional Appliance	EA	No Image
2749	CD Patient SafetyNet View Station Software (V 3.1.0.2)	EA	No Image
2572	CD Patient SafetyNet View Station Software (V 4.X.X.X)	EA	No Image
2445	Patient SafetyNet Bed License	EA	No Image
2125	Halo Bed License	EA	No Image
2447	Notification License <i>Note: For non-proprietary notification interface</i>	EA	No Image
2448	Patient SafetyNet EMR Connectivity Site License <i>Note: One required per physical address.</i>	EA	No Image
2146	SafetyNet Compatible Radical Docking Station, RDS-1 Comes with power cord and communication capabilities. <i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA	
2147	SafetyNet Compatible Radical Docking Station, RDS-3 Comes with power cord and communication capabilities. <i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA	

PN	Description	UM	Product Image
Coming soon	SafetyNet Compatible Radical-7 Touchscreen Docking Station, RDS-1 Comes with power cord and communication capabilities. <i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA	
Coming soon	SafetyNet Compatible Radical-7 Touchscreen Docking Station, RDS-3 Comes with power cord and communication capabilities. <i>SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA	
15007	Colubris Workgroup Bridge One required per pulse oximeter; alternative for 3Com Workgroup Bridge	EA	
15275	NET232-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter	EA	
15274	NET232-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter	EA	No image
16190	NET232-POE-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter	EA	No image
16191	NET232-POE-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter	EA	No image
15719	NET232-U-DTE (Male) RS232 Serial to Ethernet Intelligent Cable Adapter	EA	No image
15720	NET232-U-DCE (Female) RS232 Serial to Ethernet Intelligent Cable Adapter	EA	No image
23416	Kit HP ProCurve Serial (Male-Male)	EA	No image
23415	Kit HP ProCurve Serial (Male-Female)	EA	No image
17100	Cables to GO Patch Cable 3 ft. Blue	EA	No image
15768	Cables to GO Patch Cable 7ft. Blue	EA	No image
15769	Cables to GO Patch Cable 14 ft. Blue	EA	No image

PN	Description	UM	Product Image
15770	Cables to GO Patch Cable 25 ft. Blue	EA	No image
17644	HP V1910-24G-POE Port Switch	EA	
14530	Null Modem Cable, DB-9 Male to DB-9 Male, 1" Length	EA	No Image
14857	XRackPro2 Rackmount 4U Cabinet, Noise Reduction Enclosure	EA	
16846	HP Rackmount UPS	EA	
15010	HP Rackmount UPS – Universal	EA	
14424	Paging Transmitter UHF Frequency range (450 to 470 MHz). One required per floor in hospital.	EA	
16365	Paging Transmitter VHF Frequency range (135 to 175 MHz). One required per floor in hospital.	EA	
16262	Pager 440- 444 MHz, Synthesized, POCSAG Pager tunable in the 440 to 444 MHz range.	EA	

PN	Description	UM	Product Image
16263	Pager 444 - 448 MHz, Synthesized, POCSAG Pager tunable in the 444 to 448 MHz range.	EA	
16158	Pager 448 - 452 MHz, Synthesized, POCSAG Pager tunable in the 448 to 452 MHz range.	EA	
16264	Pager 452 - 456 MHz, Synthesized, POCSAG Pager tunable in the 452 to 456 MHz range.	EA	
14253	Pager 456 - 460 MHz, Synthesized, POCSAG Pager tunable in the 456 to 460 MHz range.	EA	
14254	Pager 460 - 464 MHz, Synthesized, POCSAG Pager tunable in the 460 to 464 MHz range.	EA	
14255	Pager 464 - 468 MHz, Synthesized, POCSAG Pager tunable in the 464 to 468 MHz range.	EA	
16527	Pager 135 – 142.999 MHz, Synthesized, POCSAG VHF Band A Pager tunable in the 135 to 142.999 MHz range.	EA	

PN	Description	UM	Product Image
16528	Pager 143 – 150.999 MHz, Synthesized, POCSAG VHF Band B Pager tunable in the 143 to 150.999 MHz range.	EA	
16366	Pager 151 – 158.999 MHz, Synthesized, POCSAG VHF Band C Pager tunable in the 151 to 158.999 MHz range.	EA	
16529	Pager 159 – 166.999 MHz, Synthesized, POCSAG VHF Band D Pager tunable in the 159 to 166.999 MHz range.	EA	
16530	Pager 167 – 174.999 MHz, Synthesized, POCSAG VHF Band E Pager tunable in the 167 to 174.999 MHz range	EA	
2456	Patient SafetyNet Installation Services	EA	No image
2457	Masimo Wireless Site Survey	EA	No image

Philips® Compatible Devices

PN	Description	UM	Product Image
2267	Masimo SET IntelliVue® Module <i>Patient Cable must be ordered separately.</i>	EA	
1881	Masimo SET® VueLink® Module Comes with VueLink cable.	EA	
1841	VueLink® Cable For use with Masimo SET VueLink Module, 3 ft.	EA	

PN	Description	UM	Product Image
18710	USB-to-RS232 Converter Cable For connecting Root with IntelliBridge(R) (EC5 portion) and VueLink(R) modules.	EA	 A white USB-to-RS232 converter cable with a blue USB-A connector on one end and a white RS232 DB9 connector on the other. The cable is coiled and has a small label with the Masimo logo.

Philips®, IntelliVue® and VueLink® are registered trademarks of Philips Corporation.

GE Healthcare Compatible Devices

PN	Description	UM	Product Image
9477	uSpO₂® Pulse Oximetry Cable Kit For use with GE ApexPro FH Telemetry Transmitters. Includes: • 1 LNCS Adtx-3 • 1 LNCS Pdtx-3 • 1 LNCS Inf-3 • 1 LNCS Neo-3 sensor <i>US Only</i> May only be sold to End-Customers with pre-existing installed ApexPro Telemetry Systems. For more information, see GE Upgrade Profiles under OEM on RadNet.	EA	 A white pulse oximetry cable kit with a blue connector on one end and a red connector on the other. The cable is coiled and has a small label with the Masimo logo.
9476	uSpO₂ Pulse Oximetry Cable Kit For use with GE ApexPro, ApexPro CH, and CARESCAPE T4 & T14 Telemetry Transmitters. Includes: • 1 LNCS Adtx-3 • 1 LNCS Pdtx-3 • 1 LNCS Inf-3 • 1 LNCS Neo-3 sensor <i>Available in countries where cleared.</i> May only be sold to End-Customers with pre-existing installed ApexPro Telemetry Systems. For more information, see GE Upgrade Profiles under OEM on RadNet.	EA	 A white pulse oximetry cable kit with a blue connector on one end and a red connector on the other. The cable is coiled and has a small label with the Masimo logo.

iSpO2®

PN	Description	UM	Product Image
3588	iSpO2 Rx 30-Pin Series, Adult Compatible with: • IOS 7 and 8 • iPhone 4 • iPhone 4S • iPad (2nd Generation) • iPad (3rd Generation) <i>Available in countries where cleared. Phone not included.</i>	EA	
3590	iSpO2 Rx 30-Pin Series, Pediatric Compatible with: • IOS 7 and 8 • iPhone 4 • iPhone 4S • iPad (2nd Generation) • iPad (3rd Generation) <i>Available in countries where cleared. Phone not included.</i>	EA	
3721	iSpO2 Rx 30-Pin Series, M-LNCS™ Compatible with: • IOS 7 and 8 • iPhone 4 • iPhone 4S • iPad (2nd Generation) • iPad (3rd Generation) <i>Available in countries where cleared. Phone not included. M-LNCS sensors sold separately.</i>	EA	

PN	Description	UM	Product Image
9575	iSpO2 Rx Lightning Series, Adult Compatible with: • OS 7 and 8 • iPhone 6 • iPhone 6 Plus • iPhone 5 • iPhone 5S • iPhone 5C • iPad (4th Generation) • iPad Mini (1st Generation) • iPad Mini (2nd Generation) • iPad Mini (3rd Generation) • iPod (5th Generation) • iPad Air • iPad Air 2 <i>Available in countries where cleared. Phone not included.</i>	EA	
9576	iSpO2 Rx Lightning Series, Pediatric Compatible with: • IOS 7 and 8 • iPhone 6 • iPhone 6 Plus • iPhone 5 • iPhone 5S • iPhone 5C • iPad (4th Generation) • iPad Mini (1st Generation) • iPad Mini (2nd Generation) • iPod (5th Generation) • iPad Air • iPad Air 2 <i>Available in countries where cleared. Phone not included.</i>	EA	

PN	Description	UM	Product Image
9569	iSpO2 Rx Lightning Series, M-LNCS Compatible with: • OS 7 and 8 • iPhone 6 • iPhone 6 Plus • iPhone 5 • iPhone 5S • iPhone 5C • iPad (4th Generation) • iPad Mini (1st Generation) • iPad Mini (2nd Generation) • iPod (5th Generation) • iPad Air • iPad Air 2 <i>Available in countries where cleared. Phone not included. M-LNCS sensors sold separately.</i>	EA	
9578	iSpO2 Rx Micro USB Series, Adult Compatible with: • Android Jelly Bean**, KitKat, and Lollipop • Samsung Galaxy Note 10 (2014 edition) • Samsung Note II • Samsung Note 3 • Samsung Note 4 • Nexus 5 • Nexus 6 • Nexus 7 • Nexus 9 • Samsung Galaxy S3 • Samsung Galaxy S4 • Samsung Galaxy S5 • Huawei Ascend Mate 7 • Droid Mini • Motorola RAZR M • Samsung S3 • Samsung S4 • Samsung Note 2 <i>Available in countries where cleared. Phone not included.</i>	EA	

PN	Description	UM	Product Image
9579	iSpO2 Rx Micro USB Series, Pediatric Compatible with: • Android Jelly Bean**, KitKat, and Lollipop • Samsung Galaxy Note 10 (2014 edition) • Samsung Note II • Samsung Note 3 • Samsung Note 4 • Nexus 5 • Nexus 6 • Nexus 7 • Nexus 9 • Samsung Galaxy S3 • Samsung Galaxy S4 • Samsung Galaxy S5 • Huawei Ascend Mate 7 • Droid Mini • Motorola RAZR M • Samsung S3 • Samsung S4 • Samsung Note 2 <i>Available in countries where cleared. Phone not included.</i>	EA	
9577	iSpO2 Rx Micro USB Series, M-LNCS Compatible with: • Android Jelly Bean**, KitKat, and Lollipop • Samsung Galaxy Note 10 (2014 edition) • Samsung Note II • Samsung Note 3 • Samsung Note 4 • Nexus 5 • Nexus 6 • Nexus 7 • Nexus 9 • Samsung Galaxy S3 • Samsung Galaxy S4 • Samsung Galaxy S5 • Huawei Ascend Mate 7 • Droid Mini • Motorola RAZR M • Samsung S3 • Samsung S4 • Samsung Note 2 <i>Available in countries where cleared. Phone not included.</i> <i>M-LNCS sensors sold separately.</i>	EA	

Sensors and Cables

rainbow Acoustic Monitoring™ Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRa)

PN	Description	UM	Product Image
3475	RAS-125c, Acoustic Respiration Cloth Sensor Adult / Pediatric Adhesive Sensors (RRa) Single patient use 10 / box, non-sterile Weight > 10kg Not made with natural rubber latex. <i>Available in countries where cleared.</i>	Box	 The image shows the RAS-125c sensor, which consists of a white adhesive cloth sensor connected to a white cable. The cable has a red connector at one end and a white connector at the other. The sensor is labeled with 'RAS-125c' and 'RRa'.
3483	RAS-125c, Short Term Monitoring Acoustic Respiration Cloth Sensor Adult / Pediatric Adhesive Sensors (RRa) Single patient use 10 / box, non-sterile Weight > 10kg Not made with natural rubber latex. <i>Available in countries where cleared.</i>	Box	 The image shows the RAS-125c sensor, which consists of a white adhesive cloth sensor connected to a white cable. The cable has a red connector at one end and a white connector at the other. The sensor is labeled with 'RAS-125c' and 'RRa'. A small white label with '6 Hour' is attached to the cable.
3660	RAM Dual Cable RC-10 For use with RAM RAS and rainbow or M-LNCS sensors	EA	 The image shows the RAM Dual Cable RC-10, which is a white cable with two red connectors at one end and two white connectors at the other.
3661	RAM Dual Cable LNC-10 For use with RAM RAS and SpO ₂ LNCS sensors	EA	 The image shows the RAM Dual Cable LNC-10, which is a white cable with two red connectors at one end and two white connectors at the other.

rainbow Adhesive Sensors (for rainbow SET Devices with SpHb, SpO₂, and SpMet)

PN	Description	UM	Product Image
3792	rainbow R1 25 Adult Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10 / box, non-sterile Weight >30 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i>	Box	
3793	rainbow R1 20 Pediatric / Slender Digit Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i>	Box	
2414	cablerainbow R1 25L Adult / Neonatal Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10 / box, non-sterile Weight <3 kg* or >30 kg Not made with natural rubber latex. Use replacement tape 2623. <i>Requires rainbow RC patient cable</i> *SpO ₂ And SpMet only	Box	

PN	Description	UM	Product Image
2415	rainbow R1 20L Infant Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10 / box, non-sterile Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement tape 2624. <i>Requires rainbow RC patient cable</i>	Box	 The image shows a single Rainbow R1 20L infant adhesive sensor. It consists of a white, rectangular adhesive pad with a small red and white logo on the left side. A thin, white cable extends from the center of the pad, with the word "rainbow" printed in small letters along its length. The cable ends in a brown, L-shaped connector with a small red and white logo on the right side.


rainbow Adhesive Sensors (for rainbow SET Devices with SpO₂, SpCO, and SpMet)

PN	Description	UM	Product Image
2221	rainbow R25 Adult Adhesive Sensors (SpO ₂ , SpCO, SpMet) Single patient use 10 / box, non-sterile Weight >30 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i>	Box	
2222	rainbow R20 Pediatric Adhesive Sensors (SpO ₂ , SpCO, SpMet) Single patient use 10 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i>	Box	
2219	rainbow R25-L Adult / Neonatal Adhesive Sensors (SpO ₂ , SpCO, SpMet) Single patient use 10 / box, non-sterile Weight <3 kg* or >30 kg Not made with natural rubber latex. Use replacement tape 2623. <i>Requires rainbow RC patient cable</i> <i>*SpO₂ And SpMet only</i>	Box	
2220	rainbow R20-L Infant Adhesive Sensors (SpO ₂ , SpCO, SpMet) Single patient use 10 / box, non-sterile Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement tape 2624. <i>Requires rainbow RC patient cable</i>	Box	


rainbow ReSposable™ Sensors (for rainbow SET Devices with SpHb, SpO₂, and SpMet)



rainbow ReSposable Sensors (for use with Masimo rainbow SET Devices with SpHb, SpO₂, and SpMet parameters)

PN	Description	UM	Product Image
3457	rainbow ReSposable R2-25 Sensor System Adult R2-25a disposable optical sensors and R2-25r reusable optical sensors (SpO₂, SpHb, SpMet) 20 R2-25a / box (DOS) 1 R2-25r / box (ROS™) Non-sterile Weight > 30 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i>	Box	
3458	rainbow® ReSposable™ R2-20 Sensor System Pediatric R2-20a disposable optical sensors and R2-20r reusable sensors (SpO₂, SpHb, SpMet) 20 R2-20a / box (DOS) 1 R2-20r / box (ROS) Non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex. <i>Requires rainbow RC patient cable</i>	Box	
3395	rainbow® ReSposable™ R2-25r Sensors Adult R2-25r reusable optical sensors (SpO₂, SpHb, SpMet) 5 R2-25r / box (ROS) Non-sterile Weight > 30 kg Not made with natural rubber latex. <i>Requires R2-25a</i>	Box	

PN	Description	UM	Product Image
3411	rainbow® ReSposable™ R2-20r Sensors Pediatric / Slender Digit R2-20r reusable optical sensors (SpO ₂ , SpHb, SpMet) 5 R2-20r / box (ROS) Non-sterile Weight 10 – 50 kg Not made with natural rubber latex. <i>Requires R2-20a</i>	Box	 A long, thin, white reusable optical sensor with a pink connector at one end and a black connector at the other. The word "rainbow" is printed on the white body.
2753	rainbow® ReSposable™ R2-25a Sensors Adult R2-25a disposable optical sensors (SpO ₂ , SpHb, SpMet) 25 R2-25a / box (DOS) Non-sterile Weight > 30 kg Not made with natural rubber latex. <i>Requires R2-25r</i>	Box	 A white, butterfly-shaped disposable optical sensor with two adhesive pads. The central sensor unit is visible. The word "rainbow" and "R2 25" are printed on the right pad.
2755	rainbow® ReSposable™ R2-20a Sensors Pediatric / Slender Digit R2-20a disposable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 25 R2-20a / box (DOS) Non-sterile Weight 10 – 50 kg Not made with natural rubber latex. <i>Requires R2-20r</i>	Box	 A white, butterfly-shaped disposable optical sensor with two adhesive pads. The central sensor unit is visible. The word "rainbow" and "R2 20" are printed on the right pad.

rainbow Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO₂, SpMet; SpHb required)

PN	Description	UM	Product Image
3418	rainbow DCI® SC 200 Adult reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 200 Spot Checks <i>Requires rainbow RC Patient Cable</i>	EA	
3419	rainbow DCI-P SC 200 Pediatric reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 200 Spot Checks <i>Requires rainbow RC Patient Cable</i>	EA	
3420	rainbow DCI SC 400 Adult reusable sensor, 3 ft. (SpHb, SpO₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 400 Spot Checks <i>Requires rainbow RC Patient Cable</i>	EA	



PN	Description	UM	Product Image
3421	rainbow DCI-P SC 400 Pediatric reusable sensor, 3 ft. (SpHb, SpO ₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 400 Spot Checks <i>Requires rainbow RC Patient Cable</i>	EA	
3647	rainbow DCI SC 1000 Adult reusable sensor, 3 ft. (SpHb, SpO ₂ and SpMet) Weight > 30 kg Not made with natural rubber latex. 1000 Spot Checks <i>Requires rainbow RC Patient Cable</i>	EA	
3648	rainbow DCI-P SC 1000 Pediatric reusable sensor, 3 ft. (SpHb, SpO ₂ and SpMet) Weight 10 – 50 kg Not made with natural rubber latex. 1000 Spot Checks <i>Requires rainbow RC Patient Cable</i>	EA	
3797	rainbow DCI-mini SC 200 infant reusable sensor. (SpHb, SpO ₂ and SpMet) Weight 3 – 30 kg Not made with natural rubber latex. 200 Spot Checks <i>Requires rainbow RC Patient Cable</i> <i>Not For Sale in USA For Export Only</i>	EA	

PN	Description	UM	Product Image
3798	rainbow DCI-mini SC 400 infant reusable sensor. (SpHb, SpO ₂ and SpMet) Weight 3 – 30 kg Not made with natural rubber latex. 400 Spot Checks <i>Requires rainbow RC Patient Cable</i> <i>Not For Sale in USA For Export Only</i>	EA	
3799	rainbow DCI-mini SC 1000 infant reusable sensor. (SpHb, SpO ₂ and SpMet) Weight 3 – 30 kg Not made with natural rubber latex. 1000 Spot Checks <i>Requires rainbow RC Patient Cable</i> <i>Not For Sale in USA For Export Only</i>	EA	

Pronto-7 Reusable Spot Check Sensors


PN	Description	UM	Product Image
9400	rainbow 4D DC Sensor, Small Reusable Sensor 1 / box, non-sterile Used with finger sizes between 8mm – 12mm diameter. Not made with natural rubber latex. Comes with 500 or 1000 Spot checks. <i>Must add Spot checks separately (PN 3312)</i> <i>Requires a Pronto-7.</i>	EA	
9401	rainbow 4D DC Sensor, Medium Reusable Sensor 1 / box, non-sterile Used with finger sizes between 12mm – 16mm diameter. Not made with natural rubber latex. Comes with 500 or 1000 Spot checks. <i>Must add Spot checks separately (PN 3312)</i> <i>Requires a Pronto-7.</i>	EA	
9402	rainbow 4D DC Sensor, Large Reusable Sensor 1 / box, non-sterile Used with finger sizes between 16mm – 20mm diameter. Not made with natural rubber latex. Comes with 500 or 1000 Spot checks. <i>Must add Spot checks separately (PN 3312)</i> <i>Requires a Pronto-7.</i>	EA	

rainbow Reusable Sensors (for rainbow SET Devices with SpCO, SpO₂, and SpMet)

PN	Description	UM	Product Image
2696	rainbow DCI Adult Reusable Sensor, 3 ft. (SpCO, SpMet, SpO ₂) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i>	EA	
2697	rainbow DCI-P Pediatric Reusable Sensor, 3 ft. (SpCO, SpMet, SpO ₂) 1 / box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i>	EA	
3657	rainbow DCI-6 Adult Reusable Sensor, 6ft (SpCO, SpMet, SpO ₂) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex <i>Requires rainbow RC Patient Cable</i>	EA	

PN	Description	UM	Product Image
3658	rainbow DCI-P-6 Pediatric Reusable Sensor, 6ft (SpCO, SpMet, SpO ₂) 1 box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex <i>Requires rainbow RC Cable</i>	EA	

rainbow Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO₂, and SpMet)

PN	Description	UM	Product Image
2201	rainbow DCI dc3 Adult Reusable Direct Connect Sensor, 3 ft. (SpO ₂ , SpCO, SpMet) 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
2069	rainbow DCI-P-dc3 Pediatric Reusable Direct Connect Sensor, 3 ft. (SpO ₂ , SpCO, SpMet) 1 / box, non-sterile Weight 10 kg - 50 kg Not made with natural rubber latex.	EA	



rainbow Patient Cables (for rainbow SET Devices)

PN	Description	UM	Product Image
2405	rainbow RC-1 rainbow 20-pin Patient Cable, 1 ft For use with rainbow and M-LNCS sensors	EA	
2406	rainbow RC-4 rainbow 20-pin Patient Cable, 4 ft For use with rainbow and M-LNCS sensors	EA	
2404	rainbow RC-12 rainbow 20-pin Patient Cable, 12 ft For use with rainbow and M-LNCS sensors	EA	
3659	rainbow RC-EXT-8 rainbow extension cable, 8ft For use with rainbow and M-LNCS sensors	EA	

PN	Description	UM	Product Image
3342	rainbow RC25-12 rainbow 25-pin, Patient Cable, 12ft For use with rainbow and M-LNCS sensors	EA	
3343	rainbow RC25-12RA rainbow 25-pin, right angle Patient Cable, 12ft For use with rainbow and M-LNCS sensors	EA	
3628	rainbow RC25-4 rainbow 25-pin, Patient Cable, 4 ft. For use with rainbow and M-LNCS sensors	EA	

PN	Description	UM	Product Image
3375	rainbow RC25-4RA rainbow 25-pin, Patient Cable, 4 ft. For use with rainbow and M-LNCS sensors	EA	 A coiled white patient cable with two red connectors. Each connector has a white plastic housing with a red latch and a small rainbow logo. The cable is neatly coiled in the center of the image.

Masimo SET ReSposable S-DOS™ Series

PN	Description	UM	Product Image
2434	ReSposable™ S-DOS™ 25 Adult SpO ₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight > 30 kg Not made with natural rubber latex. For use with S-ROS™ 3U or 12U series.	Box	 A white, butterfly-shaped adhesive sensor with two gold-colored electrical contacts in the center. The text 'S-DOS 25 SpO2' is printed on the left wing, and the Masimo logo is on the right wing.
2435	ReSposable S-DOS 20 Pediatric SpO ₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight 10 - 50 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series	Box	 A white, butterfly-shaped adhesive sensor similar to the 25 model, but with 'S-DOS 20 SpO2' printed on the left wing and a cross symbol on the right wing.
2469	ReSposable S-DOS 25L Neonatal / Adult SpO ₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight > 30 kg, < 3 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series	Box	 A long, thin, tan-colored adhesive sensor with a circular gold-colored contact on the left and a dashed line indicating the placement of the second contact on the right. The Masimo logo is on the right end.
2470	ReSposable S-DOS 20L Infant SpO ₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight 3-30 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series	Box	 A long, thin, tan-colored adhesive sensor similar to the 25L model, with a circular gold-colored contact on the left and a dashed line on the right. The Masimo logo is on the right end.

PN	Description	UM	Product Image
2471	ReSposable S-DOS 30L Large Adult SpO ₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight 3 - 10kg or > 40 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series	Box	 A long, thin, tan-colored adhesive sensor strip with a black circular sensor window in the center and a small white tab on the left. The Masimo logo is visible on the right end.
2476	ReSposable S-DOS 25Pt-L Neonatal SpO ₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight < 1 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series	Box	 A long, thin, white adhesive sensor strip with a black circular sensor window in the center and a small white tab on the left. A red dashed line runs along the center, and a red Masimo logo is on the right end.
2477	ReSposable S-DOS 20L-500 Non-adhesive Neonatal SpO ₂ disposable optical sensors (DOS) Single patient use 20 / box, non-sterile Weight < 1 kg Not made with natural rubber latex. For use with S-ROS 3U or 12U series	Box	 A long, thin, white adhesive sensor strip with a black circular sensor window in the center and a small white tab on the left. A red dashed line runs along the center, and a red Masimo logo is on the right end.

Masimo SET ReSposable S-ROS™ 3U Series

PN	Description	UM	Product Image
2473	LNCS ReSposable S-ROS 3U SpO₂ reusable optical sensor (ROS) Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series. <i>Requires patient cable.</i>	EA	
2472	M-LNCS ReSposable S-ROS 3U SpO₂ reusable optical sensor (ROS) Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series. <i>Requires patient cable.</i>	EA	

Masimo SET ReSposable S-ROS 12U Series

PN	Description	UM	Product Image
2474	ReSposable Red S-ROS 12U SpO ₂ reusable optical sensor (ROS), 12 ft., 20-pin Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	 A black cable with a red connector at one end and a 20-pin connector at the other. It features a coiled section and a small circular label with the Masimo logo.
2475	ReSposable S-ROS 12U SpO ₂ reusable optical sensor (ROS), 12 ft., 14-pin Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	 A black cable with a white connector at one end and a 14-pin connector at the other. It features a coiled section and a small circular label with the Masimo logo.
2619	ReSposable SpaceLabs S-ROS 12U SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	 A black cable with a white connector at one end and a 14-pin connector at the other. It features a coiled section and a small circular label with the Masimo logo.
2620	ReSposable Philips S-ROS 12U SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	 A black cable with a blue connector at one end and a 14-pin connector at the other. It features a coiled section and a small circular label with the Masimo logo.

PN	Description	UM	Product Image
2621	ReSposable DB9 S-ROS 12U SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	
2622	ReSposable GE S-ROS 12U SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	

ReSposable MAC S-ROS 12U Series

PN	Description	UM	Product Image
2566	ReSposable MAC 395 S-ROS 12U N-395 SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	
2567	ReSposable MAC 180 S-ROS 12U N-180 SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	
2568	ReSposable MAC GE S-ROS 12U GE SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	
2569	ReSposable MAC SL S-ROS 12U SpaceLabs SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	

PN	Description	UM	Product Image
2570	ReSposable MAC SL2 S-ROS 12U SpaceLabs SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series	EA	 A black cable with a blue and white connector at one end and a black connector with a circular logo at the other. The cable is coiled.
2571	ReSposable MAC CMS S-ROS 12U Philips SpO ₂ reusable optical sensor (ROS), 12 ft., 12-pin Reusable 1 / box, non-sterile Not made with natural rubber latex. For use with S-DOS series.	EA	 A black cable with a blue and white connector at one end and a black connector with a circular logo at the other. The cable is coiled.

LNCS Adhesive Sensors

PN	Description	UM	Product Image
1859	LNCS Adtx Adult SpO ₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.	Box	
2317	LNCS Adtx-3 Adult SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.	Box	

PN	Description	UM	Product Image
2918	LNCS E1® Adult SpO₂ disposable ear sensor, 3 ft. Single patient use 10 / box, non-sterile, non-adhesive Weight > 30 kg Not made with natural rubber latex. Requires Masimo SET MSLP board, MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology	Box	 The image shows a white, coiled ear sensor cable. At one end is a white ear clip with a small red Masimo logo. At the other end is a blue connector with a white plug.
3858	LNCS TFA-1 Adult / Pediatric SpO₂ disposable tranreflectance forehead sensor, 3ft Single patient use 10 / box non-sterile Weight > 10 kg Not made with natural rubber latex Not for sale in USA For Export only.	Box	 The image shows a white, oval-shaped forehead sensor with a red Masimo logo. It is connected to a white cable that ends in a blue connector with a white plug.

PN	Description	UM	Product Image
3857	LNCS Red TFA-1 Adult / Pediatric SpO ₂ disposable tranflectance forehead sensor 3 ft Single patient use 10 / box non-sterile Weight > 10 kg Not made with natural rubber latex Requires Masimo SET MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX (Version 7.7 or higher) technology Not for sale in USA For Export only.	Box	
1860	LNCS Pdtx Pediatric SpO ₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight 10 kg – 50 kg Not made with natural rubber latex.	Box	
2318	LNCS Pdtx-3 Pediatric SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 10 kg – 50 kg Not made with natural rubber latex.	Box	
1861	LNCS Inf-L Infant SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex.	Box	

PN	Description	UM	Product Image
1862	LNCS Neo-L Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex.	Box	
1901	LNCS NeoPt-L Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex.	Box	

LNCS Adhesive Sensors with Replaceable Tapes

PN	Description	UM	Product Image
2328	LNCS Inf Infant SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.	Box	
2319	LNCS Inf-3 Infant SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.	Box	
2329	LNCS Neo Neonatal SpO ₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.	Box	
2320	LNCS Neo-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.	Box	

PN	Description	UM	Product Image
2330	LNCS NeoPt Neonatal SpO ₂ Adhesive Sensors, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	
2321	LNCS NeoPt-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	
2331	LNCS NeoPt-500 Neonatal SpO ₂ non-adhesive sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2322.	Box	

LNCS Specialty Sensors

LNCS Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat® setting

PN	Description	UM	Product Image
2411	LNCS Trauma Adult SpO ₂ Adhesive Sensors Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	
2412	LNCS Newborn Neonatal Neonatal SpO ₂ Adhesive Sensors Single patient use 20 / box, sterilizable Weight < 3 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	
2413	LNCS Newborn Infant / Pediatrics Infant / Pediatric SpO ₂ Adhesive Sensor Single patient use 20 / box, sterilizable Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement wrap 2322.	Box	

LNCS Sterile Sensors

Note: The sensors are not intended to be used for evaluations or carried as sales samples. If there is a need for samples of sterile sensors, contact your manager for approval and include the e-mail with FRM-1386.

PN	Description	UM	Product Image
2601	LNCS Adtx-3 Adult SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight > 30 kg Not made with natural rubber latex.	Box	
2602	LNCS Pdtx-3 Pediatric SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight 10 kg – 50 kg Not made with natural rubber latex.	Box	
2603	LNCS Inf-3 Infant SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.	Box	

PN	Description	UM	Product Image
2604	LNCS Neo-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.	Box	 A photograph of the LNCS Neo-3 sensor, which consists of a small, rectangular, tan-colored adhesive pad with a black connector at one end, attached to a coiled white cable.
2605	LNCS NeoPt-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single patient use 20 / box, Sterile Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	 A photograph of the LNCS NeoPt-3 sensor, which consists of a small, rectangular, white adhesive pad with a black connector at one end, attached to a coiled white cable.
2606	LNCS NeoPt-500 Neonatal SpO ₂ non-adhesive sensor, 18 in. Single patient use 20 / box, Sterile Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2322.	Box	 A photograph of the LNCS NeoPt-500 sensor, which consists of a small, rectangular, white non-adhesive sensor with a black connector at one end, attached to a coiled white cable.

LNCS Reusable Sensors

PN	Description	UM	Product Image
1863	LNCS DC-I Adult Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
1864	LNCS DCI-P Pediatric Reusable Sensor, 3 ft. 1 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.	EA	
2653	LNCS DB-I Adult Soft Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
1895	LNCS TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	

PN	Description	UM	Product Image
2258	LNCS Y-I Multisite Reusable Sensor, 3 ft. 1 / box, non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg Not made with natural rubber latex.	EA	
1896	LNCS TF-I Adult Transflectance Forehead, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	

LNCS Red Patient Cables (for use with Masimo rainbow SET Devices)

PN	Description	UM	Product Image
2365	Red LNC™-01 LNCS 20-pin SpO ₂ Spot Check Patient Cable.	EA	 A short, curved patient cable with a clear 20-pin connector on one end and a red 20-pin connector on the other.
2055	Red LNC™-04 LNCS 20-pin SpO ₂ Patient Cable, 4 ft.	EA	 A coiled patient cable with a clear 20-pin connector on one end and a red 20-pin connector on the other.
2056	Red LNC™-10 LNCS 20-pin SpO ₂ Patient Cable, 10 ft	EA	 A coiled patient cable with a clear 20-pin connector on one end and a red 20-pin connector on the other.
2057	Red LNC™-14 LNCS 20-pin SpO ₂ Patient Cable, 14 ft.	EA	 A coiled patient cable with a clear 20-pin connector on one end and a red 20-pin connector on the other.

LNC Patient Cables (for use with Masimo SET Devices)

PN	Description	UM	Product Image
2364	LNC-1 LNCS Series Patient Cable.	EA	
2017	LNC-4 LNCS Series Patient Cable, 4 ft.	EA	
1814	LNC-10 LNCS Series Patient Cable, 10 ft.	EA	
3630	Red 25 LNC-4 LNCS Series 25-pin SpO ₂ Patient Cable, 4 ft	EA	

PN	Description	UM	Product Image
3376	Red 25 LNC-4RA Series LNCS Series 25-pin right angle SpO ₂ Patient Cable, 4 ft	EA	
3345	Red 25 LNC-10 Series LNCS Series 25-pin SpO ₂ Patient Cable, 10 ft	EA	
2013	LNC-14 LNCS Series Patient Cable, 14 ft.	EA	
2021	LNC-4-Ext LNC Extension Cable, 4 ft.	EA	

PN	Description	UM	Product Image
2281	LNC MP10 LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO ₂ LNCS Dual Key Patient Cable, 10 ft.	EA	
3755	LNC MP4 LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO ₂ LNCS Dual Key Patient Cable, 4 ft	EA	
2430	LNC-4-GE LNCS to GE Patient Cable, 4 ft.	EA	
2016	LNC-10-GE LNCS to GE Patient Cable, 10 ft.	EA	


PN	Description	UM	Product Image
2432	LNC-SL-10 LNCS to Spacelabs Patient Cable, 10 ft.	EA	
2285	LNC-04-DB9 LNCS to DB9 Connector Patient Cable, 4 ft.	EA	
2286	LNC-10-DB9 LNCS to DB9 Connector Patient Cable, 10 ft.	EA	
2287	LNC-14-DB9 LNCS to DB9 Connector Patient Cable, 14 ft.	EA	

LNCS Adapter Cables

PN	Description	UM	Product Image
1816	LNCS to PC Adapter Cable LNCS Series to LNOP PC Adapter Cable, 18 in.	EA	
2263	LNC CMS LNCS Series to Philips 12-pin Instrument Adapter Cable, 3 m.	EA	
2264	LNC MAC-GE LNCS Series to GE Adapter Cable, 3 m.	EA	
2266	LNC MAC-SL LNCS Series to Spacelabs Adapter Cable, 3 m.	EA	

PN	Description	UM	Product Image
2336	LNC MAC-SL2 LNCS Series to Spacelabs Adapter Cable, 3 m.	EA	 A coiled grey cable with a blue Spacelabs connector on one end and a blue LNC connector on the other.
2268	LNC MAC-180 LNCS Series to N-180 Adapter Cable, 10 ft.	EA	 A coiled grey cable with a grey N-180 connector on one end and a grey LNC connector on the other.
2269	LNC MAC-395 LNCS Series to N-395 Adapter Cable, 10 ft.	EA	 A coiled grey cable with a white N-395 connector on one end and a grey LNC connector on the other.

M-LNCS™ Adhesive Sensors

PN	Description	UM	Product Image
2508	M-LNCS Adtx Adult SpO ₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.	Box	 A white, butterfly-shaped adhesive sensor with a black cable and a black connector. The sensor has a central circular window with a star symbol.
2509	M-LNCS Adtx-3 Adult SpO ₂ Adhesive Sensor, 3 ft. Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex.	Box	 A white, butterfly-shaped adhesive sensor with a black cable and a black connector. The sensor has a central circular window with a star symbol.
2510	M-LNCS Pdtx Pediatric SpO ₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight 10 kg - 50 kg Not made with natural rubber latex.	Box	 A white, butterfly-shaped adhesive sensor with a black cable and a black connector. The sensor has a central circular window with a star symbol.
2511	M-LNCS Pdtx-3 Pediatric SpO ₂ Adhesive Sensor, 3 ft. Single patient use 20 / box, sterilizable Weight 10 kg - 50 kg Not made with natural rubber latex.		

PN	Description	UM	Product Image
2919	M-LNCS E1 Adult SpO ₂ disposable ear sensor, 3 ft. Single patient use 10 / box, non-sterile, non-adhesive Weight > 30 kg Not made with natural rubber latex. Requires Masimo SET MSLP board, MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology	Box	
3856	M-LNCS TFA-1 Adult / Pediatric SpO ₂ disposable transreflectance forehead sensor, 3 ft Single patient use 10 / box, non-sterile Weight > 10 kg Not made with natural rubber latex Not for sale in USA for Export only	Box	
3855	M-LNCS Red TFA-1 Adult / Pediatric SpO ₂ disposable transreflectance forehead sensor, 3 ft Single patient use 10 / box, non-sterile Weight > 10 kg Not made with natural rubber latex Requires Masimo SET MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX (Version 7.7 or higher) technology Not for sale in USA for Export only	Box	

PN	Description	UM	Product Image
3670	M-LNCS Inf-L Infant SPO ₂ adhesive sensors, 3 ft Single patient use 20 / box, sterilizable Weight 3 kg - 20 kg Not made with natural rubber latex	Box	
3668	M-LNCS Neo-L Neonatal SPO ₂ adhesive sensors, 3 ft Single patient use 20 / box, sterilizable Weight 3 kg - 20 kg Not made with natural rubber latex	Box	
3669	M-LNCS NeoPt-L Neonatal SPO ₂ adhesive sensors, 3 ft Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex	Box	

M-LNCS Adhesive Sensors with Replaceable Tapes

PN	Description	UM	Product Image
2512	M-LNCS Inf Infant SpO ₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.	Box	
2513	M-LNCS Inf-3 Infant SpO ₂ Adhesive Sensor, 3 ft. Single patient use 20 / box, sterilizable Weight 3 kg – 20 kg Not made with natural rubber latex. Use replacement tape 2307.	Box	
2514	M-LNCS Neo Neonatal / Adult SpO ₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.	Box	
2515	M-LNCS Neo-3 Neonatal / Adult SpO ₂ Adhesive Sensor, 3 ft. Single patient use 20 / box, sterilizable Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use replacement tape 2308.	Box	

PN	Description	UM	Product Image
2516	M-LNCS NeoPt Neonatal SpO ₂ Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	
2517	M-LNCS NeoPt-3 Neonatal SpO ₂ Adhesive Sensor, 3 ft. Single patient use Weight < 3 kg or > 40 kg Not made with natural rubber latex. 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	
2518	M-LNCS NeoPt-500 Neonatal SpO ₂ Non-Adhesive Sensor, 18 in. Single patient use 20 / box, sterilizable Weight < 1 kg Not made with natural rubber latex. Use replacement wrap 2322.	Box	

M-LNCS Specialty Sensors

M-LNCS Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.

PN	Description	UM	Product Image
2521	M-LNCS Trauma Adult SpO ₂ Sensor Single patient use 20 / box, sterilizable Weight > 30 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	 A white, rectangular, adhesive SpO2 sensor with a yellow connector and a coiled white cable. The sensor has a red crosshair marking on its surface.
2519	M-LNCS Newborn Neonatal Neonatal SpO ₂ Sensor Single patient use 20 / box, sterilizable Weight < 3 kg Not made with natural rubber latex. Use replacement wrap 2309.	Box	 A white, rectangular, adhesive SpO2 sensor with a yellow connector and a coiled white cable. The sensor has a red crosshair marking on its surface.
2520	M-LNCS Newborn Infant / Pediatric Infant / Pediatric SpO ₂ Sensor Single patient use 20 / box, sterilizable Weight 3 kg – 30 kg Not made with natural rubber latex. Use replacement wrap 2322.	Box	 A white, rectangular, adhesive SpO2 sensor with a yellow connector and a coiled white cable. The sensor has a red crosshair marking on its surface.

M-LNCS Cyanotic Sensors

PN	Description	UM	Product Image
2522	M-LNCS Blue Neonatal / Infant / Pediatric SpO ₂ Adhesive Sensor Single patient use 10 / box, Non-Sterile Weight 2.5 – 30kg Not made with natural rubber latex. Use replacement tape 2822.	Box	

M-LNCS Reusable Sensors

PN	Description	UM	Product Image
2501	M-LNCS DC-I Adult Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
2502	M-LNCS DCI-P Pediatric / Slender Digit Reusable Sensor, 3 ft. 1 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.	EA	
2507	M-LNCS DB-I Adult Soft Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
2503	M-LNCS TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	

PN	Description	UM	Product Image
2505	M-LNCS Y-I Multisite Reusable Sensor, 3 ft. 1 / box, non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg Not made with natural rubber latex.	EA	 A white, coiled cable with a white connector at one end and two white sensors at the other end. The sensors have small red and green LEDs.
2504	M-LNCS TF-I Adult Transflectance Forehead, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	 A white, coiled cable with a white connector at one end and a single black sensor at the other end. The sensor has a small red LED.

M-LNC Patient Cables

PN	Description	UM	Product Image
2523	M-LNC-1 M-LNCS Series Patient Cable, 1 ft.	EA	
2524	M-LNC-4 M-LNCS Series Patient Cable, 4 ft.	EA	
2525	M-LNC-10 M-LNCS Series Patient Cable, 10 ft.	EA	
2526	M-LNC-14 M-LNCS Series Patient Cable, 14 ft.	EA	

PN	Description	UM	Product Image
2831	M-LNC-MP10 M-LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO ₂ M-LNC Dual Key Patient Cable, 10 ft.	EA	 A coiled white patient cable with a blue connector on one end and a white connector with a black key on the other.
3754	M-LNC-MP4 M-LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO ₂ M-LNC Dual Key Patient Cable, 10 ft.	EA	 A coiled white patient cable with a blue connector on one end and a white connector with a black key on the other.
2655	M-LNC-4-GE M-LNCS to GE Patient Cable, 4 ft.	EA	 A coiled white patient cable with a blue connector on one end and a white connector with a red key on the other.
2533	M-LNC-10-GE M-LNCS to GE Patient Cable, 10 ft.	EA	 A coiled white patient cable with a blue connector on one end and a white connector with a red key on the other.

PN	Description	UM	Product Image
2714	M-LNC-SL-10 M-LNCS to Spacelabs Patient Cable, 10 ft.	EA	 A coiled white patient cable with a white connector on one end and a red connector on the other.
2534	M-LNC-04-DB9 M-LNCS Series to DB9 Connector Patient Cable, 4 ft.	EA	 A coiled white patient cable with a red connector on one end and a DB9 connector on the other.
2535	M-LNC-10-DB9 M-LNCS Series to DB9 Connector Patient Cable, 10 ft.	EA	 A coiled white patient cable with a red connector on one end and a DB9 connector on the other.

M-LNCS Adapter Cables

PN	Description	UM	Product Image
2542	M-LNC MAC-SL2 M-LNC Series to Spacelabs Adapter Cable, 2 ft.	EA	 A white cable with a grey connector on one end and a black connector on the other, coiled in a loose loop.
2541	M-LNC MAC-SL M-LNC Series to Spacelabs Adapter Cable, 3 m.	EA	 A white cable with a black connector on one end and a grey connector on the other, coiled in a loose loop.
2537	M-LNC MAC-180 M-LNCS Series to N-180 Adapter Cable, 10 ft.	EA	 A white cable with a grey connector on one end and a black connector on the other, coiled in a loose loop.

PN	Description	UM	Product Image
2538	M-LNC MAC-395 M-LNCS Series to N-395 Adapter Cable, 10 ft.	EA	 A white cable with a large white connector at one end and a smaller white connector at the other, coiled in a loose loop.
2540	M-LNC MAC-GE M-LNCS Series to GE Adapter Cable, 3 m.	EA	 A white cable with a white connector at one end and a blue and white connector at the other, coiled in a loose loop.
2539	M-LNC MAC CMS M-LNCS Series to Phillips 12-pin Instruments Adapter Cable, 3 m.	EA	 A white cable with a white connector at one end and a blue and white connector at the other, coiled in a loose loop.
2654	M-LNCS to LNC Adapter Cable M-LNCS Series to LNC Patient Cable, 1.5 ft.	EA	 A white cable with a black connector at one end and a red and white connector at the other, coiled in a loose loop.



PN	Description	UM	Product Image
2750	M-LNCS to LNC-MAC Adapter Cable M-LNCS Series to LNC-MAC Adapter Cable, 1.5 ft.	EA	

LNOP Adhesive Sensors

PN	Description	UM	Product Image
1001	LNOP Adt Adult SpO ₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	Box	
1025	LNOP Pdt Pediatric SpO ₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.	Box	
1800	LNOP Inf-L Infant SpO ₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight 3 kg – 20 kg Not made with natural rubber latex. Use 2307 replacement tape.	Box	
1798	LNOP Neo-L Neonatal SpO ₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex. Use 2308 replacement tape.	Box	
1651	LNOP NeoPt-L Neonatal SpO ₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight < 1 kg Not made with natural rubber latex. Use 2309 replacement tape.	Box	

LNOP Specialty Sensors

LNOP Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.

PN	Description	UM	Product Image
2358	LNOP Trauma Adult SpO ₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight >30 kg Not made with natural rubber latex. Use 2309 replacement tape	Box	
2203	LNOP Newborn Neonatal Neonatal SpO ₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight < 3 kg Not made with natural rubber latex. Use 2309 replacement tape	Box	
2204	LNOP Newborn Infant / Pediatric Infant / Pediatric SpO ₂ Adhesive Sensors Single patient use 20 / box, non-sterile Weight 3 kg – 30 kg Not made with natural rubber latex. Use 2322 replacement tape	Box	

LNOP Cyanotic Sensor

PN	Description	UM	Product Image
1970	LNOP Blue Neonatal / Infant / Pediatric SpO ₂ Adhesive Sensors Single patient use 10 / box, non-sterile Weight 2.5 kg – 30 kg Not made with natural rubber latex. Use 1801 replacement tape.	Box	

LNOP Reusable Sensors

PN	Description	UM	Product Image
1269	LNOP DC-I Adult SpO ₂ Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
1276	LNOP DCI-P Pediatric SpO ₂ Reusable Sensor, 3 ft. 1 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex.	EA	
1794	LNOP TC-I SpO ₂ Reusable Tip-Clip EAR Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
1793	LNOP TF-I Adult SpO ₂ Transflectance Forehead Reusable Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
1560	LNOP DC-195 Adult SpO ₂ Reusable Sensor 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	

PN	Description	UM	Product Image
1544	LNOP Y-I Multisite SpO ₂ Reusable Sensor, 3 ft. 1 / box, non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg Not made with natural rubber latex.	EA	

LNOP Red Patient Cables (for use with Masimo rainbow SET Devices)

PN	Description	UM	Product Image
2058	Red PC-04 LNOP Series 20-pin SpO ₂ Patient Cable, 4 ft.	EA	
2059	Red PC-08 LNOP Series 20-pin SpO ₂ Patient Cable, 8 ft.	EA	
2060	Red PC-12 LNOP Series 20-pin SpO ₂ Patient Cable, 12 ft.	EA	

PN	Description	UM	Product Image
3344	Red 25 PC-12 LNOP Series 25-pin, Patient Cable, 12 ft	EA	

LNOP Patient Cables (for use with Masimo SET Devices)

PN	Description	UM	Product Image
1173	PC04 LNOP Series Patient Cable, 4 ft.	EA	
1005	PC08 LNOP Series Patient Cable, 8 ft.	EA	
1006	PC12 LNOP Series Patient Cable, 12 ft.	EA	
3629	Red 25 PC-4 LNOP Series 25-pin SpO ₂ Patient Cable, 4 ft	EA	

PN	Description	UM	Product Image
3631	Red 25 PC-4RA LNOP Series 25-pin, right angle, SpO ₂ Patient Cable, 4 ft	EA	
1619	PC04-Ext LNOP Extension Cable, 4 ft.	EA	
2282	LNOP MP12 LNOP to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO ₂ LNOP Dual Key Patient Cable, 12 ft.	EA	
1890	PC-12-GE LNOP to GE Patient Cable, 12 ft.	EA	

PN	Description	UM	Product Image
3667	PC-4-GE LNOP to GE Patient Cable, 4 ft.	EA	
3666	PC-8-GE LNOP to GE Patient Cable, 8 ft.		
2288	PC-04-DB9 LNOP Series to DB9 Connector Patient Cable, 4 ft.	EA	
2289	PC-08-DB9 LNOP Series to DB9 Connector Patient Cable, 10 ft.	EA	

PN	Description	UM	Product Image
2290	PC-12-DB9 LNOP Series to DB9 Connector Patient Cable, 12 ft.	EA	
2433	LNOP PC-SL-12 LNOP Spacelabs Patient Cable, 12 ft.	EA	

LNOP Adapter Cables

PN	Description	UM	Product Image
1645	MAC-1 LNOP Series to LNC Adapter Cable, 18 in.	EA	
2858	M-LNCS to PC Adapter Cable M-LNCS Series to LNOP PC Patient Cables, 18 in.	EA	
2014	Dräger Patient Cable For use with Masimo SET LNOP sensors, 2 m.	EA	
1927	MAC-1 LNOP SpO ₂ Patient Cable, 12 ft.	EA	

Direct Connect Sensors (for use with Masimo SET Devices)

PN	Description	UM	Product Image
1396	DCSC Adult Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	
2651	DB-I-dc3 Adult Soft Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	

Red Direct Connect Sensors (for Masimo rainbow SET Devices)

PN	Description	UM	Product Image
2053	Red DCI-dc3 Adult Reusable Direct Connect™ SpO ₂ Sensor, 3 ft. 1 / box, non-sterile Weight > 30kg Not made with natural rubber latex.	EA	
2256	Red DCIP-dc3 Pediatric Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1 / box, non-sterile Weight 10kg – 50kg Not made with natural rubber latex	EA	
2643	Red DBI-dc3 Adult Soft Reusable Direct Connect Sensor, 3 ft. 1 / box, non-sterile Weight > 30 kg Not made with natural rubber latex.	EA	

Mountable Patient Cable Extensions and Adapter

PN	Description	UM	Product Image
2478	14-Pin to 20-Pin Adapter Cable Connect PXC2 to Masimo SET 20-pin devices, SpO ₂ only.	EA	
2371	PXC2-16 Mountable Extension Cable, 16 ft. For use with Masimo SET PC, LNC, M-LNC Patient Cables and DCSC, DC-8, DC-12, DBI-dc3, DBI-dc8 Direct Connect sensors or Masimo 14-pin to 20-pin Adapter Cable.	EA	
2373	PXC2-30 Mountable Extension Cable, 30 ft. For use with Masimo SET® PC, LNC, M-LNC Patient Cables and DCSC, DC-8, DC-12, DBI-dc3, DBI-dc8 Direct Connect sensors or Masimo 14-pin to 20-pin Adapter Cable.	EA	


Accessories

Training Kits

PN	Description	UM	Product Image
3454	SpO₂ ReSposable, Adult / Pediatric Sensor Training Kit Includes: • 1 S-DOS 25 • 1 S-DOS 20 • 1 S-DOS 25L • 1 S-DOS 20L • 1 S-DOS 30L • 1 soft wrap armband • 1 cable tether • 1 sensor application poster • Direction for use	EA	
3455	SpO₂ ReSposable, L Series Sensor Training Kit Includes: • 1 S-DOS 30L • 1 S-DOS 25L • 1 S-DOS 20L • 1 S-DOS 25Pt-L • 1 soft wrap armband • 1 cable tether • 1 sensor application poster • 1 neonatal foot with application card • Direction for use	EA	
2659	M-LNCS Sensor Training Kit Includes: • 1 M-LNCS Adtx • 1 M-LNCS Pdtx • 1 M-LNCS Inf • 1 M-LNCS Neo • 1 M-LNCS NeoPt sensor • 1 neonatal foot with application card • Direction for use	EA	

PN	Description	UM	Product Image
3466	M-LNCS Neonatal and Infant Specialty Sensors Demo Kit Includes: • 1 M-LNCS Newborn Neonatal • 1 M-LNCS Newborn Infant / Pediatric • 1 M-LNCS Neo • 1 M-LNCS NeoPt • 1 M-LNCS NeoPt-500 • 1 M-LNCS Blue • 1 Neonatal foot with application card • 1 Directions for use	EA	
2345	LNCS Sensor Training Kit Includes: • 1 LNCS Actx • 1 LNCS Pctx • 1 LNCS Inf • 1 LNCS Neo • 1 LNCS NeoPt sensor • 1 neonatal foot with application card • Direction for use	EA	
1282	LNOP Sensors Application Training Packet Includes: • 1 LNOP Adt • 1 LNOP Pdt • 1 LNOP Inf-L • 1 LNOP Neo-L • 1 LNOP NeoPt-L • 2 cotton swabs • 2 alcohol swabs • 1 neonatal foot with application card • Direction for use	EA	
3456	Acoustic Respiration Sample Pack Includes: • 1 RAS-125 • 1 RAS-125C Sensor	EA	

Bedside Device Accessories

PN	Description	UM	Product Image
1315	Radical Series Replacement Battery For use with Radical and Radical-7 Blue Screen Handheld.	EA	
2386	Radical-7 Color Screen Replacement Battery For use with Radical-7 Color Screen only.	EA	
23794	Radical-7 Touchscreen Replacement Battery, Rechargeable Lithium Ion	EA	
1395	Radical Series Handheld Lock To secure Radical and Radical-7 Blue Screen Handheld in Radical Docking Stations.	EA	
30534	Radical Series Handheld Lock Key To unlock Radical and Radical-7 Blue Screen Handhelds from Docking Stations.	EA	
1584	Radical Series Power Cord Lock To secure a power cord to Radical Docking Stations 5 / pack	EA	

PN	Description	UM	Product Image
2228	Rad-8 Power Cord Lock	EA	
30475	Radical Series Handheld Case Shown with Blue Screen Handheld Radical. Blue Radical not included.	EA	
3776	Masimo SET M15 Tester Tester for devices containing Masimo SET and connects to M-LNC or rainbow RC patient cables Testers SpO2 and PR only	EA	
2367	9-Pin Masimo SET Tester Tester for devices containing Masimo SET and connects to LNC patient cables Tests SpO ₂ and PR only.	EA	
1795	Masimo SET LNOP Tester For Oximeter & Patient Cable SpO ₂ and PR Testing. Tester for devices containing Masimo SET and connects to LNOP PC patient cables. Tests SpO ₂ and PR only.	EA	

PN	Description	UM	Product Image
2368	Masimo rainbow SET Tester Tester for devices containing Masimo SET or Masimo rainbow SET. Tests SpHb, SpCO, SpMet, PVI, PI, SpO ₂ and BPM.	EA	
3502	Masimo rainbow SET Tester Tester for devices containing Masimo SET or Masimo rainbow SET, 25-pin instrument connector. Tests SpO ₂ , PR, SpCO, SpMet and SpHb	EA	
1593	14-Pin Masimo SET Tester For Oximeter SpO ₂ and PR Testing Tester for devices containing Masimo SET, 14-pin instrument connector. Tests SpO ₂ and PR only.	EA	
3776	Masimo M15 Tester Tester for devices containing Masimo SET and connects to M-LNC or rainbow RC patient cables Testers SpO ₂ and PR only	EA	

PN	Description	UM	Product Image
3403	Rainbow Optical Test Sensor Tester for devices containing Masimo rainbow SET. Requires Masimo rainbow RC patient cable tests for SpO2, PR, SpCO, SpMet and SpHb. Not for Clinical Use	EA	
3379	Rainbow Optical Test Sensor, 3 FT. Tester for devices containing Masimo rainbow SET, 20-pin instrument connector. Tests for SpO2, PR, SpCO, SpMet and SpHb. Not For Clinical Use	EA	
3380	Rainbow Optical Test Sensor, 12 FT. Tester for devices containing Masimo rainbow SET, 20-pin instrument connector. Tests for SpO2, PR, SpCO, SpMet and SpHb. Not For Clinical Use	EA	
2227	Rad-8 Mounting Bracket	EA	
2350	Radical-7 Color Handheld Lock Locks Radical-7 Color Handheld to the Docking Station. Used with PN 2351.	EA	

PN	Description	UM	Product Image
2351	Radical-7 Color Handheld Lock Key Used for installation of Radical-7 Color Handheld Lock and unlocking handheld from docking station. Used with PN 2350	EA	
2389	Mini-Demo Connector - 30 Min Used to demo product (best used to display trending information). Includes RAM waveform	EA	
2390	Mini-Demo Connector - 2 Min Used to demo product (best used for short demonstrations). Includes RAM waveform	EA	
36517	Masimo Instrument Configuration Tool (MICT) CD For downloading trend data from Radical-7 (touchscreen only) and MOC-9 technologies through Root. Includes PDF of Directions for Use.	EA	
18313	USB-Null Modem Cable For connecting Root to Windows-based PC for downloading trend data with MICT.	EA	

PN	Description	UM	Product Image
17417	USB-to-Serial Adapter	EA	
24206	Cart, Radical, Masimo	EA	
24208	Accessory, Stack Bracket, Radical-7	EA	No Image
24209	Accessory, Bracket, Cavi Wipes	EA	No Image

Handheld Device Accessories



PN	Description	UM	Product Image
2207	Water Resistant Protective Carrying Case for Battery Operated Portable Handheld, Black This can be used with Rad-5/5v.	EA	
2208	Water Resistant Protective Carrying Case for Battery Operated Portable Handheld, Red This can be used with Rad-5/5v.	EA	
10368	Radical-7 Ambulation Pouch Pack of 10	EA	No image


PN	Description	UM	Product Image
1909	Data Transfer Download Cable Used for connecting Rad-5 (with 14-pin connector) to Personal Computer to enable download of trending data.	EA	
2063	Data Transfer Download Cable, Red – 20-Pin Used for connecting Rad-57's and 20-pin Rad-5s to Personal Computer to enable download of trending data.	EA	
1908	TrendCom, Trend Download Software Used with appropriate Data Transfer Download cable to initiate download and de-compression of Rad-5 and Rad-57 trending data to a Personal Computer. Used with serial cable to initiate download and decompression of Radical-7, Rad-87 and Rad-8 trending data.	EA	
1842	Handheld Boot, Gray For use with Rad-5, Rad-5v and all Rad-57 Handhelds.	EA	
1980	Handheld Boot, Yellow For use with Rad-5, Rad-5v and all Rad-57 Handhelds.	EA	

PN	Description	UM	Product Image
1981	Handheld Boot, Red For use with Rad-5, Rad-5v and all Rad-57 Handhelds.	EA	
1982	Handheld Boot, Orange For use with Rad-5, Rad-5v and all Rad-57 Handhelds.	EA	
2097	Handheld Boot, Royal Blue For use with Rad-5, Rad-5v and all Rad-57 Handhelds.	EA	
2098	Handheld Boot, Light Blue For use with Rad-5, Rad-5v and all Rad-57 Handhelds.	EA	
2099	Handheld Boot, Pink For use with Rad-5, Rad-5v and all Rad-57 Handhelds.	EA	

PN	Description	UM	Product Image
2357	Light Shield For use in strobing environment and high sunlight intensity areas. 5 / Pack	EA	
2815	Masimo Ambient Shield Accessory For use with rainbow Sensors. 10 / pack	EA	

Sensor Accessories

PN	Description	UM	Product Image
2624	Replacement Tapes rainbow R1 20L and R20-L For use with rainbow R1 20L and R20-L 50 / box, non-sterile Not made with natural rubber latex.	Box	
2623	Replacement Tapes For rainbow R1 25L and R25-L For use with rainbow R1 25L and R25-L 50 / box, non-sterile Not made with natural rubber latex.	Box	
1801	Replacement Tapes For LNOP Blue Sensors For use with LNOP Blue 100 / pack, non-sterile Not made with natural rubber latex.	Box	
2307	Replacement Tapes for Inf Series Sensors For use with LNCS/M-LNCS Inf/Inf-3 or LNOP Inf-L 102 / box, non-sterile Not made with natural rubber latex.	Box	
2308	Replacement Tapes For Neo Series Sensors For use with LNCS / M-LNCS Neo/Neo-3 or LNOP Neo-L 102 / box, non-sterile Not made with natural rubber latex.	Box	
2309	Replacement SoftTouch Wraps for NeoPt, NeoPt-3, Trauma and Newborn Neonatal Sensors For use with LNCS/M-LNCS NeoPt/NeoPt-3, LNOP NeoPt-L and LNCS/M-LNCS/LNOP Trauma and Newborn Neonatal 10 / pack, non-sterile Not made with natural rubber latex.	Pack	
2322	Replacement SoftTouch Wraps for NeoPt-500 and Newborn Infant / Pediatric Sensors For use with LNCS/M-LNCS NeoPt-500 and LNCS/M-LNCS/LNOP Newborn Infant/Pediatric 10 / pack, non-sterile Not made with natural rubber latex.	Pack	


PN	Description	UM	Product Image
2822	Replacement Tape for M-LNCS Blue For use with M-LNCS Blue 50 / pack, non-sterile Not made with natural rubber latex.	Pack	
1596	CleanShield® Multisite Wrap For use with Y-I Sensors Adult / Pediatric / Neonatal 100 / box, non-sterile Not made with natural rubber latex.	Box	
1597	Standard Multisite Wrap For use with Y-I Sensors Adult / Pediatric / Neonatal 100 / box, non-sterile Not made with natural rubber latex.	Box	
1598	Standard Petite Wrap For use with Y-I Sensors Adult slender digit / Ped / Neonatal adhesive 100 / box, non-sterile Not made with natural rubber latex.	Box	
1600	Accessory Adhesive Squares For use with Y-I Sensors 144 / pack, non-sterile	Pack	
1903	Transflectance Sensor Adhesive Pad For use with TF-I® Sensors 20 / pack, non-sterile	Pack	
16179	Tele Plast Philips SpO ₂ Connector Plug Compatible with Philips MMS module.	EA	

PN	Description	UM	Product Image
1602	Foam Wraps for LNOP / LNCS® / M-LNCS Y-I Sensors For use with Y-I Sensors 12 / pack, non-sterile Not made with natural rubber latex.	Pack	
1601	Accessory Clothing Clips For securing reusable sensor cables to patient clothing 5 / pack, non-sterile	Pack	
3444	Masimo® Soft Wrap Armband Used to secure the sensor cable along the patients arm. 20 / pack, non-sterile Not made with natural rubber latex.	Pack	
3929	Masimo Forehead Sensor Headbands For use with forehead sensors 20 / box, non-sterile Not made with natural rubber latex	Box	
1608	Masimo Headband For use with TF-I Sensors 5 / box, non-sterile Not made with natural rubber latex.	Box	
2215	Masimo Headband For use with TF-I Sensors 12 / box, non-sterile Not made with natural rubber latex.	Box	
21514	LNC Replacement Cover For use with all LNC patient cables 25 covers and 50 pins.	Pack	

PN	Description	UM	Product Image
3679	rainbow Cable Tether Secures rainbow ReSposable ROS to patient cable 20 / pack, non-sterile Not made with natural rubber latex.	Pack	
3443	Masimo Cable Tether Secures SpO ₂ ReSposable ROS or reusable sensor to patient cable. 20 / pack, non-sterile Not made with natural rubber latex.	Pack	







Pronto-7 Spot Check Pulse CO-Oximeter Accessories








PN	Description	UM	Product Image
3382	Boot, Pronto-7 Red	EA	
3383	Boot, Pronto-7 Yellow	EA	
3384	Boot, Pronto-7 Green	EA	
3385	Boot, Pronto-7 Light Blue	EA	
3386	Boot, Pronto-7 Pink	EA	






PN	Description	UM	Product Image
3387	Boot, Pronto-7 Dark Blue	EA	
3388	Boot, Pronto-7 Black	EA	
3389	Boot, Pronto-7 Orange	EA	
34983	Finger Gauge, Pronto-7	EA	
16782	Travel Bag, Pronto-7	EA	
16868	Bracket, Pronto-7	EA	



SatShare Interface Cables

Note: Do not interface the Radical-7 with any monitor or module not specifically listed in this chart.

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1595	AT01	Airshields Athena Monitor with Nellcor SpO₂ (Type 9200, Type 9267)	EA	
1324	DO01	Datex-Ohmeda - AS/3 Patient Monitor (M-NSAT) - System 5 Compact Monitor (M-NSAT)	EA	
1326	DO02	Datex-Ohmeda - AS/3 Patient Monitor (M-NESTPR-00-00 / M-NESTR-00-00) - CardioCap / 5 Monitor - Light Monitor - System 5 Compact Monitor (M-NESTPR-00-00 / M-NESTR-00-00)	EA	
2023	DO05	Datex-Ohmeda System 5 Compact Monitor (M-NESTPR-00-01)	EA	
2223	DO06	Datex-Ohmeda System 5 Compact Monitor (M-PRESTN..01)	EA	
2969	DO07	Datex-Ohmeda System 5 Anesthesia Monitor (E-PRESTN)	EA	
1325	DS01	Corometrics - Maternal / Fetal Monitor Model 118 - Maternal / Fetal Monitor Model 120		

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1325	DS01	Critikon Dinamap Monitor (117310)	EA	
1325	DS01	Datascope - Expert Monitor with Nellcor SpO₂ (HS-500) - Accutor Monitor with Nellcor SpO₂ - Passport Monitor with Nellcor SpO₂	EA	
1325	DS01	North American Dräger - Narcomed 4 Anesthesia Machine (MP202) - Vitalert 3200 Monitor (203)	EA	
1325	DS01	Medical Data Electronics Escort Monitor	EA	
1325	DS01	Witt Biomedical Frontend05 Patient Monitor	EA	
1365	DS03	Dynascope Fukuda Denshi DS-5300W (HS-500)	EA	
2979	DS04	Dynascope Fukuda Denshi DS-5400 (HS-500)	EA	
1357	HP03	Philips / Agilent / HP - Omnicare / Viridia Model 24 / 26 Series Monitor System (M1204A) - Series 50XM Fetal Monitor	EA	

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1321	MQ01	GE / Marquette - Solar 7000 / 8000 Series Monitor (Tram 250SL, Tram 450SL, Tram 650SL, Tram 850A&SL, and SpO2 Solar Parameter Module) - Dash3000 Monitor	EA	
1789	MS01	Datex-Ohmeda - S5 (E-Masimo-00) Dräger - Infinity Delta (SmartPOD with Masimo SET) GE / Marquette - Solar 8000 (Tram 451M, Tram 851M, Masimo SET single parameter module) Philips - IntelliVue MP40 Monitor (MMS Module [Masimo SET]) GE - Carescape B850 (Patient Data Module [Masimo SET])	EA	
1367	NK01	Nihon Kohden - Model BMS-8301A/MV-831RA (AL-800PA) - Model BMS-8502A (AL-800PA)	EA	
1368	NK02	Nihon Kohden - BSM-2301 (Type 9266) - BSM-5105 (Type 9266)	EA	
1646	OH02	Hellige SMU 611 Monitor (SAT) <i>Not for Sale in the United States</i>	EA	
1323	SL01	SpaceLabs Medical - PC Bedside Monitor (90303) [90465] - PC2 Central Monitor (90315, 90316) [90465] SpaceLabs Medical - Adult/Neonatal Monitor (90621)	EA	

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)	UM	Host Connector
1362	SL02	SpaceLabs Medical - Ultraview 1700 Bedside Monitor (90387) [90496] - Ultraview UCW Monitor (90385) [90496]	EA	
1322	SM01	Siemens Medical Systems - Sirecust 1281 Monitor(99731472E2529) - SC 6002XL, SC7000	EA	

Note: Do not interface the Radical-7 with any monitor or module not specifically listed in this chart.

SatShare Extension Cable

PN	Model	Validated Monitor Interfaces (Plug-In Parameter Modules)*	UM	Host Connector
1605		SatShare Extension Cable	EA	

Index

A

Accessories • 99
Airshields • 1

B

Bedside Device Accessories • 101

C

Capnography • 1, 15
Corometrics • 1
Critikon • 1

D

Datascope • 1
Datex-Ohmeda • 1
Direct Connect Sensors (for use with Masimo SET Devices) • 95
Dräger • 1

E

EMMA • 1, 16

G

GE / Marquette • 1
GE Healthcare Compatible Devices • 24

H

Handheld Device Accessories • 107
Hellige • 1

I

iMagox • 1
IntelliVue • 1
ISA • 1
ISA CO2 • 15
iSpO2 • 1
iSpO2® • 25

L

LNC Patient Cables (for use with Masimo SET Devices) • 63
LNCS Adapter Cables • 67
LNCS Adhesive Sensors • 51
LNCS Adhesive Sensors with Replaceable Tapes • 55
LNCS Red Patient Cables (for use with Masimo rainbow SET Devices) • 62

LNCS Reusable Sensors • 60
LNCS Specialty Sensors • 57
LNCS Sterile Sensors • 58
LNOP Adapter Cables • 94
LNOP Adhesive Sensors • 84
LNOP Cyanotic Sensor • 85
LNOP Patient Cables (for use with Masimo SET Devices) • 90
LNOP Red Patient Cables (for use with Masimo rainbow SET Devices) • 88
LNOP Reusable Sensors • 86
LNOP Specialty Sensors • 85

M

Masimo Patient SafetyNet™ System • 17
Masimo rainbow SET Parameter Options • 11
Masimo rainbow® Pulse CO-Oximetry Devices • 8
Masimo SET ReSposable S-DOS™ Series • 44
Masimo SET ReSposable S-ROS 12U Series • 47
Masimo SET ReSposable S-ROS™ 3U Series • 46
Masimo SET® Pulse Oximetry Devices • 5
Medical Data Electronics • 1
M-LNC Patient Cables • 78
M-LNCS Adapter Cables • 81
M-LNCS Adhesive Sensors with Replaceable Tapes • 72
M-LNCS Cyanotic Sensors • 75
M-LNCS Reusable Sensors • 76
M-LNCS Specialty Sensors • 74
M-LNCS™ Adhesive Sensors • 69
Monitoring Devices • 5
Mountable Patient Cable Extensions and Adapter • 97

N

Nellcor • 1
Nihon Kohden • 1

P

Philips • 1
Philips® Compatible Devices • 23
Pronto • 1
Pronto-7 • 1
Pronto-7 Reusable Spot Check Sensors • 38
Pronto-7 Spot Check Pulse CO-Oximeter Accessories • 115
PVI • 1

R

Rad-5 • 1
Rad-57 • 1
Rad-5v • 1
Rad-8 • 1

Masimo Product Catalog 2015 - Q1

Rad-87 • 1

Radical-7 • 1

rainbow • 1

rainbow Acoustic Monitoring™ Sensor and Cables (for rainbow Acoustic Monitoring Devices with RRa) • 29

rainbow Adhesive Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet) • 30

rainbow Adhesive Sensors (for rainbow SET Devices with SpO2, SpCO, and SpMet) • 32

rainbow Direct Connect Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet) • 40

rainbow Patient Cables (for rainbow SET Devices) • 41

rainbow ReSposable™ Sensors (for rainbow SET Devices with SpHb, SpO2, and SpMet) • 33

rainbow Reusable Sensors (for rainbow SET Devices with SpCO, SpO2, and SpMet) • 39

RDS • 1

Red Direct Connect Sensors (for Masimo rainbow SET Devices) • 96

ReSposable MAC S-ROS 12U Series • 49

Root • 1

RRa • 1

RRp • 1

S

SatShare Extension Cable • 120

SatShare Interface Cables • 117

SedLine® Brain Function Monitoring • 14

Sensor Accessories • 111

Sensors and Cables • 29

Shield • 1

Siemens Medical Systems • 1

SpaceLabs • 1

SpCO • 1

Specialty Monitoring and Connectivity Platform • 13

SpHb • 1

rainbow Reusable Spot Check Sensors (for rainbow SET Devices with SpHb, SpO2, SpMet) • 35

SpMet • 1

T

Tape • 1

Tether • 1

Training Kits • 99

U

uSpO2 • 1

V

VueLink • 1

W

Witt Biomedical • 1

Wrap • 1



www.masimo.com

LAB-3422AA 0215