

JUPITER VAPORISER

Featuring..

- » For “out of circuit” use in continuous flow techniques of inhalation anaesthesia.
- » Fitted with interlock block for greater safety and ability to mount in series.
- » Calibrated specifically for the specified anaesthetic drug using a laser refractometer.
- » Temperature, flow and duration of use compensation.
- » Safety lock to help prevent accidental turn-on of dial.
- » Unique colour coded parts relevant to drug bottle colour to prevent errors when administering the anaesthetic agent.
- » Dimensions: 108mm x 200-215mm (filler dependent) x 212.45mm.
- » Manufactured in the UK.



More information...

- » The new Jupiter Vaporiser features a larger sump, interlock block and a non-spill mechanism. It also comes with a 5 year service period and it is considerably cheaper than other versions on the market. Available in all drug types in quick fill and key fill.



FOR MORE INFORMATION OR PRICING, PLEASE CONTACT OUR SALES TEAM



+44 (0)1756 700687



sales@norvap.com



www.norvap.com



Units 1 & 2 Union Business Park, Snaygill
Industrial Estate, Skipton, BD23 2QR

