

VS-800
Vital Signs Monitor

Technical Specifications

Table with 2 columns: Specification Category and Details. Categories include Safety, Dimension and Weight, Operation Environment, Storage Environment, Performance Specifications, and NIBP.

Table with 2 columns: Specification Category and Details. Categories include Accuracy of blood pressure measurement, Over-pressure protection, MINDRAY SPO2, MASIMO SPO2, NELLCOR SPO2, and Temperature.

VS-800
Vital Signs Monitor



★ 3-Year Manufacturer Warranty on Main Unit!

MINDRAY logo and USA Corporate Office contact information: 8650 154th Ave., NE Redmond, WA 98052. Tel: 425-881-0361 Fax: 425-881-0879 E-mail: usa@mindray.com Toll free: 888-816-8188

Canada Office contact information: #206-4250 Kingsway, Burnaby, BC V5H 4T7. Tel: 604-451-1199 Fax: 604-451-1189 E-mail: canada@mindray.com Toll free: 800-656-9088

Taking high cost out of quality healthcare






Parameters: NIBP, SpO₂ , Temp



The Barcode Scanner
Let you admit patient quickly and reduce the chance for transcription errors by simply scanning the barcode to get patient information at the point of care.



SmarTemp™ Thermometer
2 types of probes: oral/axillary probe (blue) and rectal probe (red). Probe can auto-warm up to an initial temperature to be sure of an accurate measurement. Neonatal patients are suitable for axillary measurement.



Powerful Data Storage and Management Functions
Allow you to easily keep and view patient data.
Total storage capacity: 100 patient IDs; 1200 sets of data for each parameter
All data can be transferred to a PC through Windows-based software for review and A4 paper printing.



VS-800 Vital Signs Monitor

- Suitable for adult, pediatric and neonatal patients
- Clinical proven SpO₂ technology
- SmarTemp™ thermometer
- Adjustable audible and visual alarms
- Central Station networking
- Optional barcode scanner for easily admitting the patients
- Removable and rechargeable battery, up to 10-hour working time
- Integrated handle and full mounting solutions

