

## Microcentrifuge, MiniStar



This personal centrifuge has a fixed speed of 6000 min<sup>-1</sup>/2000 xg and is designed for applications requiring relatively low g-forces, such as microfiltration, cell separation and quick spin downs from the walls of tubes. The tube rotor spins up to eight individual 1,5/2,0 ml plastic microcentrifuge tubes. It will also accept 0,5 ml and 0,4 ml tubes with the supplied adapters. The additional PCR rotor is designed to centrifuge four 8 x 0,2 ml PCR strip tubes. It will also accept up to 32 individual 0,2 ml tubes.

- Quick release rotors: No tools are required to exchange the rotors or attach them to the motor shaft
- Easy to operate: Close the lid of the centrifuge to start a run; press and release the tab on lid to stop

**Delivery information:** Supplied with 8x1,5/2,0 ml rotor, 4x(8x0,2 ml) strip tube rotor, 8 adapters for 0,4 ml tubes and 8 adapters for 0,5 ml tubes.

Model	MiniStar
Max. capacity	8x1,5/2,0 ml tubes or 4x(8x0,2 ml) tube strips
Max. speed (min <sup>-1</sup> )	6000
Max. RCF (xg)	2000
Noise level (dB(A))	56
Power consumption (W)	5
WxDxH (mm)	132x132x112
Weight (kg)	0,55
Acceleration/deceleration time (s)	5/1,5



Description	Pk	Cat. No.
Microcentrifuge MiniStar with rotors and adapters, EU-plug	1	521-2161
Microcentrifuge MiniStar with rotors and adapters, UK-plug	1	521-2162
Microcentrifuge MiniStar with rotors and adapters, CH-plug	1	521-2163



## Microcentrifuge, MiniStar silverline



This centrifuge was developed for applications that do not require high speed centrifugation, such as microfiltration and cell separation, and can also be used in a cold room. On closing the lid, the rotor quickly accelerates to maximum speed. Pressing the button on the front of the unit stops operation and opens the lid. Alternatively, the unit can be turned on and off with the switch at the back.

- Metal base with ABS body and transparent safety cover
- Quiet and compact design
- Starts and stops in seconds
- Easy to operate, for eight microtubes or four strip tubes

**Ordering information:** MiniStar silverline supplied with 8x1,5/2,0 ml rotor, 4x(8x0,2 ml) strip tube rotor, eight adapters for 0,5 ml and eight adapters for 0,4 ml tubes.

Model	MiniStar silverline
Max. capacity	8x1,5/2,0 ml tubes or 4x(8x0,2 ml) tube strips
Max. speed (min <sup>-1</sup> )	6000
Max. RCF (xg)	2000
Noise level (dB(A))	57
WxDxH (mm)	132x155x110
Weight (kg)	0,75

Description	Pk	Cat. No.
MiniStar silverline with rotors and adapters, EU-plug	1	521-2844
MiniStar silverline with rotors and adapters, UK-plug	1	521-2845
MiniStar silverline with rotors and adapters, CH-plug	1	521-2846

Description	Pk	Cat. No.
<b>Accessories</b>		
Rotor for 8x1,5/2,0 ml tubes	1	521-2847
Rotor for 4 strips of 8x0,2 ml tubes	1	521-2848