

Tissue-Tek® DRS™ 2000

Random Access & Multiple Slide Stainer



Tissue-Tek® DRS™ 2000 Multiple Slide Stainer

The Tissue-Tek® DRS™ 2000 Multiple Slide Stainer offers multiple staining protocols in a true random-access mode. Based on current workload demands, the system can handle up to 11 sets of 40 slides at a time. The operator selects from 20 methods stored in memory, presses START and walks away. It is a simple but intelligent instrument for unsurpassed productivity - slide after slide, shift after shift.

Specifications

Throughput	Up to 440 slides per hour
Dimensions	W x D x H 115 x 49 x 66 cm
Weight	90 kg
Number of reservoirs	
Start	1 to 3 stations
Reagent	27 reservoirs
Drying	1 station
Wash	0 to 4 stations
Capacity	
Reagent volume	650 mL
Staining protocols	11 runs continuously loading, and simultaneously running
Baskets	20 slide basket
Loading	Up to 40 slides per run
Temperature	
Drying station	30 ~ 65°C

Item code	Description
4931	Tissue-Tek® DRS™ 2000 with Dryer 220 V, 50/60 Hz
2008	Tissue-Tek® DRS™ Carbon Filter; 2/cs
4768	Tissue-Tek® 20-Slide Basket
4974	Tissue-Tek® Solution Reservoir 650 mL
4975	Tissue-Tek® Wash Reservoir 650 mL
4976	Tissue-Tek® Reservoir Lid
4977	Tissue-Tek® Basket Hook; 10/cs
4978	Tissue-Tek® Basket Adapter
4979	Tissue-Tek® Wash Station Plug
4981	Tissue-Tek® Drain Hose
4985	Tissue-Tek® Wash Reservoir Grommet
4986	Tissue-Tek® Lid 6 Large Solution Reservoir
4993	Tissue-Tek® Lid 7 Large Solution Reservoir
HIS40614	Tissue-Tek® Insert for slide basket (for 5 large slides maximum 60 mm wide)
6506	Tissue-Tek® Exhaust hose assembly 38 mm Ø
6507	Tissue-Tek® Exhaust hose assembly 75 mm Ø
SA-4931-UT	User Training Tissue-Tek® DRS™ 2000
SC-4931-G	Service Agreements Tissue-Tek® DRS™ 2000 Gold
SC-4931-P	Service Agreements Tissue-Tek® DRS™ 2000 Platinum
SC-4931-S	Service Agreements Tissue-Tek® DRS™ 2000 Silver

