

ELECTRODE PRODUCT GUIDE



Product Name	Differentiation Features	Location Used	Defibrillator
CPR-D-padz 	- Enables Real CPR Help® and See-Thru CPR® - One-piece design simplifies proper positioning on the chest - One size fits all, with ability to expand for larger chest	- EMS & Fire - Public Access	- M Series® - E Series® - R Series® - AED Plus® - AED Pro® - X Series®
CPR Stat-padz 	- Enables Real CPR Help and See-Thru CPR - Multi-function, two-piece design ensures all therapies can be delivered - Accommodates AP and AA pad placement	- Hospital - Med/Surg floors - ER & ICU - EMS & Fire	- M Series - E Series - R Series - AED Plus - AED Pro - X Series
OneStep Basic 	- Multi-function electrode - Automated self-test with R Series	Hospital	- M Series - R Series - X Series
OneStep Pacing 	- OneStep Basic features plus: – OneStep Pacing	Hospital	R Series
OneStep CPR (AP) 	- OneStep Basic features plus: – Real CPR Help – See-Thru CPR – Anterior–Posterior placement	Hospital	- M Series - R Series - X Series
OneStep CPR (AA) 	- OneStep Basic features plus: – Real CPR Help – See-Thru CPR – Anterior–Anterior placement	Hospital	- M Series - R Series - X Series
OneStep Complete 	All OneStep Basic, Pacing, and CPR capabilities	Hospital	R Series
OneStep Pediatric 	- Designed for patients 0–8 years of age and less than 55 pounds (25 kg) - Automatically adjusts joules setting to pediatric levels	Hospital	- M Series - R Series - X Series
OneStep Pediatric CPR 	- Designed for patients 0–8 years of age and less than 55 pounds (25 kg) - Automatically adjusts joules setting to pediatric levels - Real CPR Help - See-Thru CPR	Hospital	- M Series - R Series - X Series
Pedi-padz® II 	- Connects directly to AEDs & ALS devices - Pediatric intelligence function with ZOLL AEDs	- EMS & Fire - Public Access	- M Series - E Series - R Series - AED Plus - AED Pro - X Series
Pedi-padz Multi-function Liquid or Solid Gel 	- Multi-functional electrode for pediatric patients - Available in liquid or solid gel	- Hospital - EMS & Fire	- M Series - E Series - R Series - X Series
Pedi-padz Radiolucent Solid Gel 	- For use in cath and RP labs for procedures - Virtually invisible to x-ray	- Hospital – Cardiac cath lab – Electrophysiology labs	- M Series - E Series - R Series - X Series
Pro-padz® Multi-function Liquid or Solid Gel 	- Available in liquid or solid gel - Available with or without tray - Used in cath and EP labs for procedures	- Hospital – Cardiac catheterization lab – Electrophysiology labs	- M Series - E Series - R Series - X Series
Pro-padz Radiolucent Liquid or Solid Gel 	- For use in cath and EP labs for procedures - Virtually invisible to x-ray - Pre-Connect (solid only) - Speed pack	- Hospital – Cardiac catheterization lab – Electrophysiology labs	- M Series - E Series - R Series - X Series
Pro-padz Sterile 	- For use in operating rooms - Helps maintain sterile barrier - Packaged in two pouches - Available in 54" or 10-foot lead wires	- Hospital – Operating rooms	- M Series - E Series - R Series - X Series
Stat-padz® 	- Multi-function - Pre-Connect - Speed pack - Different pad shapes and colors	- EMS & Fire - Hospital	- M Series - E Series - R Series - X Series
Stat-padz II 	Connects directly to AEDs and ALS devices	- EMS & Fire - Hospital	- M Series - E Series - R Series - AED Plus - AED Pro - X Series

See reverse side for proper cable/connector requirements.

ELECTRODE GUIDE



CPR-D-padz®



CPR Stat-padz®



OneStep™ Complete



OneStep™ CPR AP



OneStep™ CPR AA



OneStep™ Pediatric CPR



AED Plus®



- Connects directly for Real CPR Help® and defibrillation

- Connects directly for Real CPR Help and defibrillation

- Not compatible

- Not compatible

- Not compatible

- Not compatible

AED Pro®



- Connects directly for Real CPR Help, See-Thru CPR®, and defibrillation

- Connects directly for Real CPR Help, See-Thru CPR, and defibrillation

- Not compatible

- Not compatible

- Not compatible

- Not compatible

X Series®



- Requires CPR Connector (PN# 8000-0370) for CPR feedback
- Connects directly for defibrillation

- Requires CPR Connector (PN# 8000-0370) for CPR feedback
- Connects directly for defibrillation

- Connects directly for all capabilities
- No OneStep Pacing Capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

E Series®



- Requires CPR Connector (PN# 8000-0370) for CPR feedback
- Connects directly for defibrillation

- Requires CPR Connector (PN# 8000-0370) for CPR feedback
- Connects directly for defibrillation

- Not compatible

- Not compatible

- Not compatible

- Not compatible

R Series®



- Requires R Series® CPR Connector (PN# 8009-0020) for CPR Feedback
- Connects directly for defibrillation

- Requires R Series CPR Connector (PN# 8009-0020) for CPR Feedback
- Connects directly for defibrillation

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

Shelf Life

5 yrs

2 yrs

2 yrs

2 yrs

2 yrs

2 yrs

Quantity and Part Numbers

Individual: 8900-0800-01

Individual: 8900-0402
8 per case: 8900-0400

Individual: 8900-0224-01
8 per case: 8900-0214-01

Individual: 8900-0223-01
8 per case: 8900-0213-01

Individual: 8900-0225-01
8 per case: 8900-0217-01

Individual: 8900-000219-01
8 per case: 8900-000220-01