

BenchTop Pipettor

Accelerate Liquid Transfers



- 96/384 semi-automated pipettor
- Highly precise, accurate and convenient
- Portable and easy to use
- Aspirates and dispenses liquid samples with the touch of a button



Sorenson™
BioScience, Inc.

sorbio.com

BenchTop Pipettor

Accelerate Liquid Transfers



Additional Features:

- Calibrates to different fluid viscosities
- Easy glide 3-slot axis positioning provides easy movement
- Adjustable aspiration and dispensing speeds
- Height adjustment settings over plate or reservoir



Engaging



Aspirating



Dispensing

BenchTop Pipettor Specifications

Sorenson catalog #73960

Dispensing Precision	<2% @ 25 μ l – 250 μ l	<10% @ 5 μ l – 20 μ l
Dispensing Accuracy	+/- 4% @ 25 μ l – 250 μ l	+/- 6% @ 5 μ l – 20 μ l
Resolution	1 μ l through menu, or 5 μ l using rotary knob	
Pipette Volume	1.0 μ l – 250 μ l	
Dimensions	8.5" x 18" footprint; 17" high (21.6 cm x 45.7 cm footprint; 43.2 cm high)	
Weight	21 lb (9.5 kg)	
Operating Temperature	10 – 37 °C	



Contact the Sorenson marketing department at 801-266-9334 or email sorbio@corning.com for additional information.

sorbio.com