

Be Ready at All Times

SE-3

3-Channel Electrocardiograph



SE-3

3-Channel Electrocardiograph




SE-3 (Color)

SE-3 (Monochrome)



SE-1515 Data Management System is able to export multi-formats reports: XML/SCP/PDF/DICOM/JPG/BMP.

Patient information can be automatically transferred from SE-1515 to SE-3, typo error can be avoided. You can upload reports to SE-1515 automatically to manage your data, which provide a convenient method to share, review, edit and diagnose ECG reports between different doctors.



User-friendly Design

- Foldable LCD screen, vivid waveform display
- Built-in rechargeable Li-ion battery, AC/DC power supply

Streamlined Workflow

- Support barcode scanner to simplify patient information input
- One-button operation
- Auto/Manual/Rhythm/Off work modes
- Automatic baseline adjustment
- Lead-off detection and alarm

Accurate Clinical Performances

- 12-lead simultaneous acquisition with optional 24-bit A/D converter
- Anti-noise technology with comprehensive filters and the high CMRR over 110dB (optional)
- Automatic ECG measurements and interpretation tested with CSE/AHA/MIT databases
- Complete digital filters, resisting baseline drift, AC and EMG interference
- Extended printing once arrhythmia is detected to further your rhythm analysis

Flexible Printing Solutions

- Built-in high-resolution thermal printer
- Grid printing on normal thermal paper
- 12-lead grid printing on normal A4 paper through external printer

Data Management

- Internal storage up to 500 ECGs and can be enlarged by USB flash disk
- Data transmission to PC via series port/LAN/USB
- Report formats: PDF and optional SCP/FDA-XML/DICOM data export
- PC-based ECG data management software (optional)




Foldable Screen



Portable Design
(Handle Optional)



One-button Operation





SE-3
3-Channel Electrocardiograph



Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.

CE 0123
ENG-ECG-SE-3
V1.2-20240630